

INTERPRETATION TABLE																									
HLA-A low resolution SSP typing																									
Amplification patterns of the A*0101 to A*8001 alleles																									
		Well ⁶																							
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.		120	135	205	190	160	165	175	165	75	85	80	185	175	200	90	180	180	200	160	155	240	85	75	360
PCR product		145	200	235		535	200	205	200			175		225	240	135	210			200	240	375	240	160	
			255									500													
Length of int.		800	1070	1070	800	800	800	1070	800	800	1070	1070	800	1070	1070	1070	1070	1070	1070	800	800	800	800	800	1070
pos. control ¹																									
5'-primer(s) ²		98	78	363	98	144	106	98	98	266	257	302	103	98	98	203	97	180	98	180	78	28	78	176	341
		5'-CTT ^{3'}	5'-TCT ^{3'}	5'-ATA ^{3'}	5'-CTA ^{3'}	5'-GCC ^{3'}	5'-CCA ^{3'}	5'-CTC ^{3'}	5'-CTA ^{3'}	5'-ACG ^{3'}	5'-Cgg ^{3'}	5'-ggA ^{3'}	5'-CCT ^{3'}	5'-CTT ^{3'}	5'-CAC ^{3'}	5'-gAA ^{3'}	5'-TCA ^{3'}	5'-TTT ^{3'}	5'-CAC ^{3'}	5'-TTT ^{3'}	5'-TCT ^{3'}	5'-TCg ^{3'}	5'-TCT ^{3'}	5'-gCA ^{3'}	5'-ggA ^{3'}
		123	106		413	317	176	368	102	266	261	385	423	423		363					391	261	527	261	
		5'-AgT ^{3'}	5'-CCA ^{3'}		5'-CCg ^{3'}	5'-gCT ^{3'}	5'-gCA ^{3'}	5'-gTT ^{3'}	5'-ACA ^{3'}	5'-ACg ^{3'}	5'-AAC ^{3'}	5'-ggC ^{3'}	5'-gCT ^{3'}	5'-gCT ^{3'}		5'-ATA ^{3'}					5'-ACg ^{3'}	5'-AAC ^{3'}	5'-TgC ^{3'}	5'-AAC ^{3'}	
			144				368		413												391				
			5'-gCA ^{3'}				5'-gTT ^{3'}		5'-CCg ^{3'}												5'-ACg ^{3'}				
			144																						
			5'-TCA ^{3'}																						
3'-primer(s) ³		203	240	527	256	265	256	259	259	302	299	341	257	282	257	299	238	317	256	299	265	97	265	292	418
		5'-TCT ^{3'}	5'-ggA ^{3'}	5'-CCA ^{3'}	5'-CTg ^{3'}	5'-CCC ^{3'}	5'-CCC ^{3'}	5'-gTT ^{3'}	5'-gTT ^{3'}	5'-ggC ^{3'}	5'-TCg ^{3'}	5'-CgT ^{3'}	5'-gCA ^{3'}	5'-gAC ^{3'}	5'-gCA ^{3'}	5'-CCA ^{3'}	5'-CCT ^{3'}	5'-ggA ^{3'}	5'-CCC ^{3'}	5'-CCA ^{3'}	5'-CCC ^{3'}	5'-ggT ^{3'}	5'-CCC ^{3'}	5'-gTT ^{3'}	5'-gTC ^{3'}
			292	527	559	570	292	502	538			521	559	559	299	411	265		259	299	282	355	282	292	
			5'-gTg ^{3'}	5'-CCT ^{3'}	5'-CCg ^{3'}	5'-CCg ^{3'}	5'-gTT ^{3'}	5'-CTT ^{3'}	5'-CCA ^{3'}			5'-ggg ^{3'}	5'-CTC ^{3'}	5'-CCg ^{3'}	5'-TCg ^{3'}	5'-TCA ^{3'}	5'-CCC ^{3'}		5'-gTT ^{3'}	5'-CCA ^{3'}	5'-gAC ^{3'}	5'-gAC ^{3'}	5'-gAC ^{3'}	5'-gTg ^{3'}	
				555			502	539					559							341	506		570	299	
				5'-CCA ^{3'}			5'-CTg ^{3'}	5'-TCT ^{3'}					5'-CgT ^{3'}							5'-CgT ^{3'}	5'-TgT ^{3'}		5'-CAC ^{3'}	5'-TCT ^{3'}	
			555			538																			
			5'-gCA ^{3'}			5'-CAG ^{3'}																			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ⁴	ser ⁵																								
*01010101-0104N, 0106, 0109-0111N, 0116N-0118N, 0120, 0122N-0127N, 0129	A1, null	1			4																				24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*0107	-	1			4																				
*0108, 0114, 0115N	1	1																							24
*0112, 0119	-	1		3																					24
*0113	-	1			4						11			13		15									24
*0121	-	1		3	4																				24
*0128	-	1			4						11					15									24
*0130	-	1																							
*02010101-0209, 0211-0216, 0218- 0222, 0224-0233, 0236-0238, 0240- 0245, 0247, 0249, 0251, 0253N, 0254, 0257-0261, 0263- 0269, 0271, 0272, 027401-0277, 0279- 0292, 0294N, 0296, 0297, 0299, 9201, 9202, 9204-9209, 9211-9213N, 9215, 9216, 9218-9221, 9223-9228, 9	A2, A203 , null		2																						
*0210, 0239, 0250, 0252, 0273, 0293, 0295, 9210, 9214, 9217	A210		2																		20				
*021701-021702	A2		2					w																	
*0234, 023501, 023503, 0256, 0262, 0270	-		2											13							20		22		
*023502	-													13							20		22		
*0246, 0270	-		2																		20		22		
*0248	-																				20		22		
*0255	-		2						8														21		23
*0278	-		2		4																20		22		
*0293	-		2											13							20		22		
*9222	-					6															20				
*9224	-		2							9															
*9229	-																				20		22		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

[illegible]

[illegible]

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
----------	--	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*680101-6804, 6806-6814, 6816-6819, 6821-6828, 6830-6841	A68								8												20				
*6805, 6815, 6820	-								8												20			23	
*6829	-								8										18		20	21			
*6901	A69								8													21			
*7401-7409, 7411, 7412N	A74																			19					
*7410	-																			19				24	
*8001	A80						6																	23	24
HLA-A allele ⁴	ser ⁵																								
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
B*1827							6																		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24