

[illegible]

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*12:07	1	2		4	5	6	7		9			
*12:08	1	2		4	5	6		8	9			
*12:09	1	2	3	4		6	7	8	9			
*12:10	1	2		4	5	6	7	8	9			
*12:11	1	2		4	5	6	7	8	9			
*12:12	1	2		4	5	6	7		9			
*12:13	1	2		4	5	6	w		9	10		
*12:14	1	2			5	6	7	8				
*12:15	1	2			5	6	7			10		
*12:16	1	2		4	5					10		
*12:17	1	2		4	5	6	7	8	9			
*12:18	1	2		4		6	7		9	10		
*12:19	1	2		4	5				9	10	11	
*12:20	1	2		4	5	6	7		9	10		
*12:21		2		4	5	6	7		9	10		
*12:22	1			4	5			8				
*12:23	1	2		4	5	6			9	10		
*12:24N	1	2		4	w	6	7	8	9			
*12:25	1	2		4	5	6	7	8	9			
*12:26	1	2		4	5	6	7		9	10		
*01:01:01-01:34												
*08:01:01-08:01:05, 08:07-08:08, 08:11, 08:16, 08:26, 08:39										10		
*08:02:01-08:02:03, 08:09, 08:21			3							10		
*08:03:02, 08:14-08:15, 08:23, 08:27, 08:29, 08:30:02, 08:33, 08:35-08:36, 08:38								8				
*08:04:01, 08:04:04-08:04:05, 14:15			3							10	11	
*08:04:02-08:04:03			3							10	w	
*08:05		2								10		
*08:06, 14:73										10	11	
*08:10								8			11	
*08:12						6		8				
*08:13			3									
*08:17	1									10		
*08:18, 08:40		2						8				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	135	215	165	105	165	250	215	195	165	195	250	170
PCR product(s)												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

*08:19, 08:34					5			8				
*08:22						6				10		
*08:24		2	3							10		
*08:25		2			5			8				
*08:28	1		3							10	11	
*08:30:01			3					8				
*08:31		2								10	11	12
*08:32				4				8				
*08:37	1							8				
*08:41		2								10		12
*11:67	1									10	11	12
*11:76, 13:34, 13:64, 14:41, 14:77												
*13:17	1	2	3								11	
*14:04, 14:50, 14:76, 14:79											11	
*14:11											11	12
*14:28						6						
*14:31		2									11	
*14:52		2	3								11	
*15:01:01:01-15:51, 16:01:01-16:05:02, 16:07-16:17												
DRB1 allele⁴												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

INTERPRETATION TABLE												
DRB1*12 SSP subtyping												
Amplification patterns of the DRB1*12:01 to 12:26 alleles												
Well ⁵												
13	14	15	16	17	18	19	20	21	22	23	24	
115	135	200	75	90	90	195	220	105	105	110	170	Length of spec.
185					135				220			PCR product(s)
255												
515	515	515	430	430	515	430	430	430	430	430	430	Length of int.
												pos. control ¹
13 (125)	149 (534)	16 (133)	25 (160)	-16 (40)	16 (133)	16 (133)	25 (161)	16 (133)	26 (165)	65 (280)	16 (133)	5'-primer ²
5'-gTg ³	5'-Cag ³	5'-gTT ³	5'-TgT ³	5'-CAA ³	5'-gTT ³	5'-gTT ³	5'-gCT ³	5'-gTT ³	5'-TTA ³	5'-AgC ³	5'-gTT ³	
37 (197)			26 (165)	65 (280)						152 (543)		
5'-gTT ³			5'-TTC ³	5'-AgC ³						5'-gAT ³		
61 (269)												
5'-CTA ³												
85 (341)	181 (630)	69 (293)	38 (199)	-3 (79)	32 (182)	67 (286)	85 (341)	38 (199)	47 (227)	85 (341)	59 (262)	3'-primer ³
5'-Cag ³	5'-CTT ³	5'-CTC ³	5'-Cag ³	5'-AgC ³	5'-TAC ³	5'-gAg ³	5'-Cag ³	5'-Cag ³	5'-ggT ³	5'-Cag ³	5'-CTg ³	
				85 (341)	47 (226)				86 (344)	179 (624)		
				5'-Cag ³	5'-gAg ³				5'-CAC ³	5'-ACA ³		
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
												DRB1 allele ⁴
								21				*12:01:01
								21				*12:01:02-12:01:03
								21				*12:02:01, 12:02:03-12:02:05
								21				*12:02:02
								21				*12:03:02
								21				*12:04
13								21				*12:05
	14							21				*12:06
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
115	135	200	75	90	90	195	220	105	105	110	170	Length of spec.
185					135				220			PCR product(s)

255												
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
		15						21				*12:07
			16					21				*12:08
								21				*12:09
				w				21				*12:10
					18			21				*12:11
						19		21				*12:12
							20	21				*12:13
13												*12:14
13								21				*12:15
								21	22			*12:16
								21		23		*12:17
								21			24	*12:18
								21				*12:19
13								21				*12:20
								21	22			*12:21
								21	22			*12:22
			16					21				*12:23
13								21				*12:24N
				17				21		23		*12:25
					18			21				*12:26
										23		*01:01:01-01:34
												*08:01:01-08:01:05, 08:07-08:08, 08:11, 08:16, 08:26, 08:39
												*08:02:01-08:02:03, 08:09, 08:21
												*08:03:02, 08:14-08:15, 08:23, 08:27, 08:29, 08:30:02, 08:33, 08:35-08:36, 08:38
												*08:04:01, 08:04:04-08:04:05, 14:15
												*08:04:02-08:04:03
												*08:05
												*08:06, 14:73
												*08:10
												*08:12
						19						*08:13
												*08:17
												*08:18, 08:40
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
115	135	200	75	90	90	195	220	105	105	110	170	Length of spec.
185					135				220			PCR product(s)
255												
13	14	15	16	17	18	19	20	21	22	23	24	Well No.

												*08:19, 08:34
												*08:22
												*08:24
												*08:25
												*08:28
												*08:30:01
												*08:31
								21	22			*08:32
												*08:37
												*08:41
												*11:67
			16									*11:76, 13:34, 13:64, 14:41, 14:77
												*13:17
												*14:04, 14:50, 14:76, 14:79
												*14:11
												*14:28
												*14:31
												*14:52
	w											*15:01:01:01-15:51, 16:01:01-16:05:02, 16:07-16:17
												DRB1 allele ⁴
13	14	15	16	17	18	19	20	21	22	23	24	Well No.