

INTERPRETATION TABLE								
DRB5 SSP subtyping								
Amplification patterns of DRB5*0101 to 0205 alleles								
	Well ⁶							
	1	2	3	4	5	6	7	8
Length of spec.	255	210	225	100	150	145	150	215
PCR product								
Length of int.	515	515	430	430	515	430	430	430
pos. control ¹								
5'-primer(s) ²	13	13	13	38	37	38	37	13
	5'-g TA 3'	5'-g TA 3'	5'-g TA 3'	5'-AC T 3'	5'-Ag g 3'	5'-AC g 3'	5'-Ag A 3'	5'-g TA 3'
	13	13				38		13
	5'-g TA 3'	5'-g TA 3'				5'-AC g 3'		5'-g TA 3'
3'-primer(s) ³	85	67	71	57	72	72	72	71
	5'-C AA 3'	5'-gAA 3'	5'-g CC 3'	5'-gC g 3'	5'-gC g 3'	5'-gC g 3'	5'-gC g 3'	5'-g Cg 3'
		72	74					71
		5'-gC g 3'	5'-CAg 3'					5'-g Cg 3'
			78					
			5'-CAC 3'					
Well No.	1	2	3	4	5	6	7	8
DRB5 allele ⁴								
*010101	1	2	3	4	5			
*010102	1	2	3	4				
*0102	1	2	3			6	7	
*0103	1	2				6	7	8
*0104	1	2	3	4				
*0105	1	2	3		5	6		
*0106				4				8
*0107	1	2	3	4	5			
*0108N	1	2	3			6	7	
*0109	1	2	3	4	5			
*0110N	1	2	3			6	7	
*0111	1			4				8
*0112, DRB1*0907 ⁵	1	2	3					
*0113	1	2	3		5			
*0202								8
*0203	1							8
*0204		2						8
*0205			3					
DRB5 allele ⁴								
Well No.	1	2	3	4	5	6	7	8

INTERPRETATION TABLE								
DRB5 SSP subtyping								
Amplification patterns of DRB5*0101 to 0205 alleles								
Well ⁶								
9	10	11	12	13	14	15	16	
175	130	200	185	195	210	240	140	Length of spec.
225								PCR product
430	430	430	515	430	430	430	430	Length of int.
								pos. control ¹
13	38	13	37	108	13	13	37	5'-primer(s) ²
5'-g TA ^{3'}	5'-AC T ^{3'}	5'-g TA ^{3'}	5'-Ag A ^{3'}	5'-Ag A ^{3'}	5'-g TA ^{3'}	5'-g TA ^{3'}	5'-Ag A ^{3'}	
58	67	67	85	160	70	78	70	3'-primer(s) ³
5'-C CT ^{3'}	5'-gAT ^{3'}	5'-gAT ^{3'}	5'-C Ag ^{3'}	5'-CAT ^{3'}	5'-gTT ^{3'}	5'-CAC ^{3'}	5'-CTg ^{3'}	
74						79		
5'-CAg ^{3'}						5'-T gC ^{3'}		
9	10	11	12	13	14	15	16	Well No.
								DRB5 allele ⁴
								*010101
								*010102
								*0102
								*0103
9								*0104
								*0105
		11						*0106
	10							*0107
				13				*0108N
					14			*0109
						15		*0110N
		11						*0111
						15		*0112, DRB1*0907 ⁵
9								*0113
		11	12					*0202
		11						*0203
			12					*0204
			12				16	*0205
								DRB5 allele ⁴
9	10	11	12	13	14	15	16	Well No.