

Well No.										1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	2	2	2	3	3	3	3	3	3	4	4	4	4	4	4		
HLA-B allele ^{2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	2	2	2	2	2	2	3	3	3	3	3	3	3	4	4	4	4	4	4	
*15:31	+	+			+				+																												+	+		+		+
*15:32	+	+	+																			+				+											+		+	+	+	
*15:33	+	+																				+		+			+										+	+	+	+	+	
*15:34	+	+	+																	+		+																+	+	+	+	
*15:35	+	+	+																			+					+										+	+	+	+	+	
*15:36	+	+				+																												+			+		+	+	+	
*15:37					+		+			+			+	+									+						+								+	+	+	+	+	
*15:38:01-15:38:02	+	+	+							+													+						+								+	+	+	+	+	
*15:39:01-15:39:02, 15:219	+	+																				+														+	+	+	+	+	+	
*15:40	+	+																											+								+	+	+	+	+	+
*15:42	+	+																												+												
*15:43	+	+	+																			+	+									+						+	+	+	+	
*15:44		+				+	+						+																		+							+			+	
*15:45	+	+	+								+				+								+															+	+	+	+	+
*15:46	+	+	+					+															+			+											+	+	+	+	+	
*15:47:01-15:47:02	+							+																						+						+		+	+	+	+	
*15:48	+	+																					+									+					+	+	+	+	+	
*15:49	+							+				+																		+							+	+	+	+	+	
*15:50	+	+	+																										+								+	+	+	+	+	+
*15:51					+		+						+										+										+				+	+	+	+	+	
*15:52					+		+		+				+																	+							+	+	+	+	+	
*15:53, 15:212	+		+					+															+			+										+	+	+	+	+	+	
*15:54	+		+					+															+													+	+	+	+	+	+	
*15:55	+	+			+				+		+																									+	+	+	+	+	+	
*15:56		+	+									+											+														+	+	+	+	+	
*15:57	+	+	+																				+														+	+	+	+	+	
*15:58	+	+	+																				+													+	+	+	+	+	+	
*15:60	+	+	+																				+														+	+	+	+	+	
*15:61	+							+															+													+	+	+	+	+	+	
*15:62	+						+	+													+															+	+	+	+	+	+	
*15:63	+	+	+																				+													+	+	+	+	+	+	
*15:64	+				+		+																+													+	+	+	+	+	+	
*15:65	+	+																					+													+	+	+	+	+	+	
*15:66		+	+										+										+													+	+	+	+	+	+	
*15:67								+										+		+	+	+	+	+											+	+	+	+	+	+		
*15:68			+							+													+													+	+	+	+	+	+	
*15:69	+							+																						+						+	+	+	+	+	+	
*15:70	+	+	+																				+													+	+	+	+	+	+	
*15:71	+	+	+																				+													+	+	+	+	+	+	
*15:72					+		+						+											+		+									+	+	+	+	+	+		
*15:73	+	+	+																				+								+		+			+	+	+	+	+	+	
*15:74								+												+			+													+	+	+	+	+	+	
*15:75	+	+	+																				+						+							+	+	+	+	+	+	
*15:76		+	+			+																	+													+	+	+	+	+	+	
*15:77	+	+	+	+			+																														+	+	+	+	+	
*15:78:01-15:78:03	+	+	+																			+		+												+	+	+	+	+	+	
*15:79N	+	+	+										+										+													+	+	+	+	+	+	
*15:80					+	+	+						+										+												+	+	+	+	+	+	+	
*15:81	+	+	+	+																			+													+	+	+	+	+	+	
*15:82	+	+	+																				+													+	+	+	+	+	+	
*15:83	+	+												+																+						+	+	+	+	+	+	
*15:84	+	+							W	+													+												+	+	+	+	+	+		
*15:85	+	+	+				+													+		+													+	+	+	+	+	+		
*15:86	+	+	+						+												+									+		+			+	+	+	+	+	+		
*15:87	+	+	+															+					+													+	+	+	+	+	+	
*15:88	+	+				+	+																													+	W			+	+	
*15:89	+	+				+	+																													+	+	+	+	+	+	
*15:90			+			+		+					+	+									+												+	+	+	+	+	+		
*15:91	+							+		+									+				+												+	+	+	+	+	+		
*15:92	+	+	+	+																		+														+	+	+	+	+	+	
Well No.										1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	2	2	2	2	3	3	3	3	3	3	3	4	4	4	4	4	4
HLA-B allele ^{2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	2	2	2	2	2	2	2	3	3	3	3	3	3	3	4	4	4	4	4	4
Well No.					</																																					

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

Well No.									1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	2	2	2	2	3	3	3	3	3	3	4	4	4	4	4	4	
HLA-B allele ^{2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7				
*40:158																																									
*40:174																						+						+													
*41:12, 54:02																																									
*42:09, C*07:68, C*07:140, C*07:151																													+												
*44:22, 44:105																																									
*44:43:01-44:43:02																	+					+																			
*44:144																	+																								
*45:09								+									+				+								+												
*46:01:01, 46:01:03-46:01:07, 46:04, 46:07N, 46:15N-46:16, 46:20, 46:22-46:23, 46:27-46:28		+	+																			+												W		+		+	+	+	+
*46:01:02, 46:02, 46:05		+	+																			+														+		+	+	+	+
*46:03		+																																+	W		+		+	+	+
*46:06		+							+							+								+											W		+	+	+		
*46:08		+																					+												W		+	+	+	+	+
*46:09		+	+																				+												W		+		+	+	+
*46:10		+	+																				+												W		+		+	+	+
*46:11		+							+										+			+								+					W		+	+	+	+	+
*46:12		+	+							+													+												W		+		+	+	+
*46:13:01-46:13:03, 46:26		+							+														+												W		+	+			+
*46:14		+	+																				+												W		+		+	+	+
*46:17		+	+		</																																				

Well No.											1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	2	2	2	3	3	3	3	3	3	3	3	3	4	4	4	4	4	4	4	4		
HLA-B allele ^{2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	
*51:64											+					+							+						+																			
*51:82																																																
*51:97								+															+																									
*51:126																							+					+																				
*52:01:01-52:01:10, 52:04-52:06:02, 52:09, 52:11, 52:13, 52:15, 52:17, 52:20, 52:22-52:23, 52:25, 52:27, 52:29, 78:05-78:06																																																
*52:01:11, 52:26																																																
*52:07																																																
*52:08																																																
*52:14								+																																								
*52:21			+																																													
*52:24																																																
*52:28																																																
*53:22																																																
*54:03, 55:08, 55:51, 56:01:01-56:02, 56:04, 56:08-56:09, 56:11, 56:13-56:17, 56:20:01-56:20:02, 56:24, 56:26-56:30, 56:33-56:34																																																
*54:06			+																																													
*54:14, 55:13, 55:23, 55:32								+																																								
*54:20, 55:09, 55:52																																																
*55:01:01-55:01:07, 55:01:09, 55:05, 55:11, 55:15, 55:17, 55:24, 55:28-55:29, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:53-55:55N																																																
*55:16																																																
*55:21			+																																													
*56:05:01-56:05:02, 78:01:01-78:01:02, 78:03																																																
*56:06, 56:25																																																
*56:21																																																
*56:22								+																																								
*56:31																																																
*56:35																																																
*57:01:01-57:01:02, 57:01:04-57:01:13, 57:03:01-57:03:02, 57:06-57:10, 57:14-57:15, 57:17-57:18, 57:20-57:23, 57:27, 57:31, 57:33, 57:35-57:38, 57:40, 57:44, 57:46-57:48, 57:50, 57:52-57:53, 57:55-57:56																																																
*57:02:01-57:02:02, 57:12, 57:19, 57:28N, 57:30, 57:42																																									</							

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]


[illegible]

[illegible]

[illegible]

[illegible]

[illegible]



[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

Lengt of spec.
PCR product(s)
Length of int. pos control ¹
Well No.
HLA-B allele ^{2,3}
*15:01:01:01, 15:01:03, 15:01:06-15:01:13, 15:01:15- 15:01:16, 15:01:18-15:01:28, 15:217, 15:231-15:232, 15:244
*15:01:01:02N
*15:01:02
*15:01:04
*15:01:14
*15:01:17
*15:02:01-15:02:05
*15:03:01-15:03:02, 15:210, 15:220, 15:242 ¹⁰
*15:03:03
*15:04
*15:05:01-15:05:02
*15:06
*15:07:01-15:07:02
*15:08
*15:09
*15:10:01
*15:10:02
*15:11:01-15:11:02, 15:11:05
*15:11:03
*15:11:04
*15:12
*15:13:01-15:13:02
*15:14
*15:15
*15:16:01-15:16:03
*15:17:01:01-15:17:02, 15:216
*15:18:01, 15:18:04
*15:18:02
*15:18:03
*15:19
*15:20
*15:21
*15:23
*15:24
*15:25:01-15:25:03, 15:240
*15:26N
*15:27:01-15:27:03
*15:28
*15:29
*15:30
Well No.
HLA-B allele ^{2,3}

Well No.
HLA-B allele ^{2,3}
*15:31
*15:32
*15:33
*15:34
*15:35
*15:36
*15:37
*15:38:01-15:38:02
*15:39:01-15:39:02, 15:219
*15:40
*15:42
*15:43
*15:44
*15:45
*15:46
*15:47:01-15:47:02
*15:48
*15:49
*15:50
*15:51
*15:52
*15:53, 15:212
*15:54
*15:55
*15:56
*15:57
*15:58
*15:60
*15:61
*15:62
*15:63
*15:64
*15:65
*15:66
*15:67
*15:68
*15:69
*15:70
*15:71
*15:72
*15:73
*15:74
*15:75
*15:76
*15:77
*15:78:01-15:78:03
*15:79N
*15:80
*15:81
*15:82
*15:83
*15:84
*15:85
*15:86
*15:87
*15:88
*15:89
*15:90
*15:91
*15:92
Well No.
HLA-B allele ^{2,3}
Well No.

HLA-B allele ^{2,3}
*15:93
*15:94N
*15:95
*15:96
*15:97
*15:98
*15:99
*15:101
*15:102
*15:103
*15:104
*15:105
*15:106
*15:107
*15:108
*15:109
*15:110
*15:111N
*15:112
*15:113
*15:114
*15:115
*15:116
*15:117
*15:118
*15:119
*15:120, 15:228 ⁵
*15:121
*15:122
*15:123
*15:124
*15:125
*15:126
*15:127
*15:128
*15:129
*15:131
*15:132
*15:133
*15:134
*15:135
*15:136
*15:137
*15:138
*15:139
*15:140
*15:141
*15:142
*15:143
*15:144
*15:145
*15:146
*15:147
*15:148
*15:149N
*15:150
*15:151
*15:152
*15:153
*15:154
Well No.
HLA-B allele ^{2,3}
Well No.
HLA-B allele ^{2,3}

*15:155
*15:156
*15:157
*15:158
*15:159
*15:160
*15:161
*15:162
*15:163
*15:164
*15:165
*15:166, 15:193⁶
*15:167
*15:168
*15:169, 15:225⁷
*15:170
*15:171
*15:172
*15:173
*15:174
*15:175
*15:176
*15:177
*15:178
*15:179
*15:180, 56:03⁴
*15:181N
*15:182N
*15:183
*15:184
*15:185
*15:186
*15:187
*15:188
*15:189
*15:190N
*15:191
*15:192
*15:194
*15:195
*15:196
*15:197
*15:198
*15:199
*15:200
*15:201
*15:202, 15:239
*15:203
*15:204
*15:205
*15:206
*15:207
*15:208, 15:230⁸
*15:209N
*15:211
*15:213
*15:214
*15:215
*15:218Q
Well No.
HLA-B allele^{2,3}
Well No.
HLA-B allele^{2,3}
*15:221
*15:222

*15:223
*15:224
*15:226N
*15:227
*15:229
*15:233
*15:234
*15:235
*15:236
*15:237
*15:238
*15:241
*15:243
*07:02:01-07:02:10, 07:02:12-07:02:17, 07:02:19-07:02:25, 07:02:27-07:03, 07:07-07:13, 07:15-07:16, 07:18:01-07:18:02, 07:21-07:22:02, 07:27, 07:30-07:31, 07:33, 07:35-07:39, 07:41-07:42, 07:44, 07:46-07:47, 07:49N-07:51, 07:54, 07:56:01-07:59, 07:61-07:64, 07:67N, 07:70-07:71, 07:73-07:77, 07:79, 07:81-07:85, 07:87-07:89, 07:91-07:95, 07:98-07:99, 07:101-07:104, 07:106-07:111N, 07:113-07:114, 07:116-07:119, 07:121-07:122, 07:124-07:130, 07:132-07:136, 07:139, 07:141-07:145, 07:148-07:154, 13:18, 18:21, 40:32, 40:98, 40:128, 48:08, C*04:108, C*15:08
*07:02:11
*07:05:01-07:05:02, 07:05:04-07:06, 07:32, 07:40, 07:53, 07:69, 07:78, 07:90, 07:97, 07:105, 07:112, 07:123, 07:137-07:138, 07:140, 40:15-40:16, 40:23, 40:161, 48:05
*07:05:03
*07:20, 07:24, 07:60
*07:23
*07:29, 08:20, 08:79
*07:34, 07:80, 13:16, 27:20, 27:24, 38:10, 38:32, 39:35, 40:95, 40:148, 40:183, 48:15, 55:01:08, 55:31, C*15:24
*07:45, 07:52, 39:08, 39:11
*07:65, 07:72, 07:86, 54:10, 55:37
*07:68:01-07:68:03
*07:100
*07:120
Well No.
HLA-B allele ^{2,3}
Well No.
HLA-B allele ^{2,3}

*07:131, 13:03, 18:58, 35:01:28, 35:57, 38:30, 40:63, 40:92, 41:17, 44:10, 44:14, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 45:14, 48:25, 56:07, 57:11, 57:16, 57:45, 57:51, 58:01:10-58:01:11, 58:11, 59:04, 82:01-82:03, C*03:81, C*15:15
*08:01:01-08:01:10, 08:01:12- 08:01:17, 08:04-08:05, 08:07- 08:19N, 08:22, 08:24, 08:27, 08:29-08:34, 08:36, 08:38- 08:39, 08:41-08:49, 08:51- 08:53, 08:55, 08:57-08:65, 08:67N-08:69, 08:73-08:78, 08:80-08:88, 40:43, 45:06, 54:12, 55:12, 55:49
*08:01:11, 08:56, 55:03, 55:25
*08:21, 08:25
*08:26, 08:50, 14:11, 14:28, 39:04, 39:63, 40:166, 44:90, 54:01:01-54:01:02, 54:04- 54:05N, 54:07-54:09, 54:11, 54:13, 54:15-54:19, 54:21- 54:22, 54:24-54:26, 55:02:01- 55:02:06, 55:04, 55:07, 55:14, 55:19-55:20, 55:22, 55:26- 55:27, 55:30, 55:34-55:35, 55:39, 55:41-55:43, 55:46- 55:48, 55:50, 55:56, 56:10, 56:12, 56:18-56:19N, 56:32
*08:28, 08:35, 08:37
*08:40
*08:70
*08:71
*13:01:01-13:01:06, 13:06- 13:07N, 13:12-13:13, 13:17, 13:21-13:23, 13:25-13:26, 13:28 13:29, 13:36, 13:43, 13:50- 13:52, 40:137, 44:08, 44:57, 44:60, 57:01:03, 57:04-57:05, 57:25-57:26, 57:32, 57:41, 57:43
*13:04, 49:04
*13:20
*13:31, 39:36
*13:39
*13:41
*13:46, 39:06:01-39:06:02, 39:34, 39:50, 39:57, 39:62, 39:64, 40:86, 40:159, 45:02, 49:05, 50:09, 59:03, 73:01- 73:02
*13:48, 40:05, 40:28
*14:05, 18:14
Well No.
HLA-B allele ^{2,3}
Well No.
HLA-B allele ^{2,3}

*14:06:01, 18:01:01:01- 18:01:05, 18:01:07-18:01:16, 18:03-18:07:02, 18:09, 18:15, 18:17N, 18:20, 18:23N-18:26, 18:28, 18:31-18:34, 18:37- 18:51, 18:53-18:55, 18:59- 18:62, 18:64-18:71, 39:43, 39:48, 40:113, 40:163
*14:13
*18:12, 18:63
*18:13, 18:30
*18:19
*18:29
*18:35, 40:77, 40:87:01- 40:87:02, 40:121, 44:150, 48:26
*18:36, 18:57, 39:32
*18:56, 44:17
*27:25, 44:146, 45:05, A*26:68, C*02:56, C*06:20, C*12:50
*27:32, 38:19, 39:03, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 40:35, 42:10, 44:54, 44:106
*27:69
*35:01:01:01-35:01:27, 35:01:29-35:03:05, 35:03:07- 35:04:03, 35:06-35:07, 35:09:01 35:09:03, 35:11:01-35:12:03, 35:17, 35:19-35:21, 35:23- 35:27, 35:30, 35:32:01- 35:32:02, 35:34, 35:36, 35:38- 35:42:02, 35:47-35:49, 35:52- 35:56, 35:59, 35:63-35:64, 35:68:01-35:68:02, 35:70, 35:74 35:75, 35:77-35:78, 35:81- 35:85, 35:90-35:92, 35:94- 35:95, 35:98, 35:101:01- 35:101:02, 35:103-35:104, 35:106-35:109, 35:111-35:114, 35:116, 35:118-35:130N, 35:132-35:134N, 35:137- 35:141, 35:143-35:155, 35:159-35:160, 35:162, 35:165N-35:171, 35:173N, 35:175, 35:177- 35:184, 35:189, 35:191, 35:193- 35:197, 35:200, 48:02:01- 48:02:02, 51:01:31, 51:02:03, 51:08-51:09:02, 51:15, 51:19, 51:40, 51:44N, 51:72-51:73, 51:87, 51:90-51:91, 51:94, 51:108, 51:115, 51:122, 53:01:01-53:01:06, 53:02-53:06, 53:08:01-53:10, 53:12-53:13, 53:15, 53:18-53:21, 53:25- 53:27, 58:01:01-58:01:09, 58:02, 58:04, 58:06, 58:09- 58:10N, 58:12-58:17N, 58:19, 58:21, 58:23-58:26, 58:28- 58:34, 78:04
Well No.

HLA-B allele ^{2,3}
Well No.
HLA-B allele ^{2,3}
*35:03:06, 35:15, 35:33, 35:35, 35:65Q, 35:88, 35:93, 35:96, 35:131, 35:157, 35:164, 35:188, 35:190, 51:10, 51:17-51:18, 51:31, 51:36, 51:81, 51:85, 51:93, 53:01:07, 53:07, 58:05, 58:27, 58:35, 81:03
*35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:199, 53:14, 58:18
*35:08:01-35:08:05, 35:18, 35:61, 35:99-35:100, 35:105, 35:142, 35:156, 35:158, 35:172, 35:174, 35:176, 35:187, 35:192, 51:05, 51:54, 53:11, 53:16, 53:24, 58:22
*35:10, 35:13, 35:16, 35:28, 52:03, 52:10:01-52:10:02, 52:12, 52:19, 53:17:01-53:17:02
*35:14:01-35:14:02, 35:62, 35:67, 35:102, 35:135, 35:163
*35:29:01-35:29:02, 35:76, 35:136, 35:198
*35:37, 51:37, 51:92
*35:43, 35:185
*35:44, 51:61
*35:45
*35:46, 58:36
*35:50
*35:66
*35:69
*35:71, 58:07
*35:79
*35:80
*35:86
*35:110
*35:115
*35:117
*35:161
*35:186, 51:20, 53:23
*37:01:01-37:01:05, 37:01:07-37:05, 37:09-37:21, 37:24-37:27, 37:29-37:33N, 40:185, 44:83, 44:134, 50:11
*37:07
*37:22, 49:03
*37:23
*39:33
*40:01:06
*40:12
*40:26
*40:51
*40:52, 40:59, 48:23
*40:71, 44:135
*40:79, 40:190, 41:14
*40:150
Well No.
HLA-B allele ^{2,3}

Well No.
HLA-B allele ^{2,3}
*40:158
*40:174
*41:12, 54:02
*42:09, C*07:68, C*07:140, C*07:151
*44:22, 44:105
*44:43:01-44:43:02
*44:144
*45:09
*46:01:01, 46:01:03-46:01:07, 46:04, 46:07N, 46:15N-46:16, 46:20, 46:22-46:23, 46:27-46:28
*46:01:02, 46:02, 46:05
*46:03
*46:06
*46:08
*46:09
*46:10
*46:11
*46:12
*46:13:01-46:13:03, 46:26
*46:14
*46:17
*46:18
*46:19
*46:21:01-46:21:02
*46:24
*46:25
*46:29
*46:30
*48:12, 48:14
*48:19
*51:01:01-51:01:18, 51:01:20- 51:01:30, 51:01:33-51:02:02, 51:02:04-51:04, 51:06:02- 51:07:02, 51:11N-51:12, 51:14, 51:22-51:24:04, 51:26-51:28, 51:30, 51:32-51:33, 51:35, 51:38-51:39, 51:41N, 51:43, 51:46, 51:48-51:53, 51:55- 51:58, 51:60, 51:65-51:71, 51:74-51:80, 51:83-51:84, 51:86, 51:88-51:89, 51:95- 51:96, 51:98N-51:100, 51:102- 51:107, 51:109-51:114, 51:116- 51:121, 51:123-51:125, 51:127- 51:131, 52:02, 52:16, 52:18, 58:08:01-58:08:02, 78:02:01- 78:02:02
*51:01:19
*51:01:32, 51:16, 51:21, 51:34, 51:101
*51:06:01
*51:13:01-51:13:02, 51:63
*51:29
*51:45
*51:59
*51:62
Well No.
HLA-B allele ^{2,3}

Well No.
HLA-B allele ^{2,3}
*51:64
*51:82
*51:97
*51:126
*52:01:01:01-52:01:10, 52:04-52:06:02, 52:09, 52:11, 52:13, 52:15, 52:17, 52:20, 52:22-52:23, 52:25, 52:27, 52:29, 78:05-78:06
*52:01:11, 52:26
*52:07
*52:08
*52:14
*52:21
*52:24
*52:28
*53:22
*54:03, 55:08, 55:51, 56:01:01-56:02, 56:04, 56:08-56:09, 56:11, 56:13-56:17, 56:20:01-56:20:02, 56:24, 56:26-56:30, 56:33-56:34
*54:06
*54:14, 55:13, 55:23, 55:32
*54:20, 55:09, 55:52
*55:01:01-55:01:07, 55:01:09, 55:05, 55:11, 55:15, 55:17, 55:24, 55:28-55:29, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:53-55:55N
*55:16
*55:21
*56:05:01-56:05:02, 78:01:01-78:01:02, 78:03
*56:06, 56:25
*56:21
*56:22
*56:31
*56:35
*57:01:01-57:01:02, 57:01:04-57:01:13, 57:03:01-57:03:02, 57:06-57:10, 57:14-57:15, 57:17-57:18, 57:20-57:23, 57:27, 57:31, 57:33, 57:35-57:38, 57:40, 57:44, 57:46-57:48, 57:50, 57:52-57:53, 57:55-57:56
*57:02:01-57:02:02, 57:12, 57:19, 57:28N, 57:30, 57:42
*57:13
*57:24
*57:29
*57:34
*57:39
*57:49, 57:57
*57:54
*78:07
*83:01
Well No.
HLA-B allele ^{2,3}
Well No.
HLA-B allele ^{2,3}

A*68:56
C*03:48
C*03:87, C*05:27, C*05:39
C*06:33
C*07:02:30
C*07:46
C*07:186
C*07:213
C*08:16:02
C*12:36
Well No.
HLA-B allele ^{2,3}