

INTERPRETATION TABLE						
KIR HLA Ligand						
	Well					
	1	2	3	4	5	6
Length of spec. PCR product	340	340	350	350	370	350
Length of int. pos. control ¹	800	800	800	1070	1070	1070
5'-primer(s) ²	1 st I	1 st I	1 st I	1 st I	1 st I	1 st I
	5'-CgA ³	5'-CgA ³	5'-CAg ³	5'-CAg ³	5'-gCA ³	5'-CAg ³
3'-primer ³	302	302	309	309	317	310
	5'-ggC ³	5'-ggT ³	5'-gTg ³	5'-ATC ³	5'-ggA ³	5'-ggT ³
Well No.	1	2	3	4	5	6
HLA allele						
C*01:02:01-01:13, 01:15-01:46, 01:48-01:58, 01:60, 01:62-01:75, 02:27-01-02:27:02, 03:02:01-03:03:14, 03:03:16-03:04:16, 03:04:18-03:06, 03:08-03:09, 03:11:01-03:11:02, 03:13-03:14, 03:16-03:28, 03:30-03:44, 03:46-03:114, 03:116:01-03:129, 03:131-03:133, 03:135-03:139, 03:141-03:162, 03:164-03:189N, 04:11, 04:29, 04:36, 04:55, 06:11, 07:01:01-07:02:40, 07:03-07:06, 07:08, 07:10-07:33N, 07:35-07:48, 07:50-07:75, 07:77-07:114, 07:116-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:246, 07:248-07:294, 07:296-07:307, 08:01:01-08:09, 08:11-08:63, 08:65-08:77, 12:02:01-12:03:25, 12:06-12:08, 12:10:01-12:20, 12:22-12:26, 12:28-12:32, 12:34-12:40, 12:42Q-12:53, 12:55-12:59, 12:61-12:71, 12:73-12:95, 14:02:01-14:03, 14:05-14:11, 14:13-14:48, 14:50-14:51, 15:07, 15:25, 15:43, 16:01:01-16:01:13, 16:04:01, 16:06-16:08, 16:10-16:11, 16:13-16:18, 16:20-16:24, 16:26-16:36, 16:38-16:45, 16:49-16:56	1					
C*01:14, 01:59, 02:02:01-02:02:03, 02:02:05-02:02:11, 02:02:13-02:11, 02:13-02:26:03, 02:28-02:40, 02:42-02:67, 03:07, 03:15, 03:45, 03:130, 03:140, 03:163, 04:01:01-04:01:28, 04:01:30-04:01:51, 04:03-04:10, 04:12-04:20, 04:23-04:28, 04:30-04:35, 04:37-04:54, 04:56-04:145, 05:01:01-05:01:26, 05:03-05:91N, 06:02:01-01-06:02:01:02, 06:02:03-06:02:11, 06:02:13-06:10, 06:12-06:51, 06:53:01-06:99, 07:07, 07:09, 07:49, 07:76, 07:210, 07:238, 07:247, 08:10, 12:04:01-12:05, 12:09, 12:21, 12:33, 12:41, 12:54, 12:60, 14:04, 14:12, 14:49, 15:02:01-15:06:03, 15:08-15:13, 15:15-15:19, 15:22-15:24, 15:26-15:42, 15:44-15:66, 16:02:01-16:02:09, 16:09, 16:12, 16:19, 16:25, 16:46-16:48, 16:57, 17:01:01-01-17:19, 18:01-18:06	2					
C*02:12, 03:03:15, 03:10, 03:29, 12:72, 15:21, 16:37	w					

Well No.	1	2	3	4	5	6
Length of spec. PCR product	340	340	350	350	370	350
Well No.	1	2	3	4	5	6
B*07:27, 15:43, 18:54, 27:03-27:07:03, 27:09-27:11, 27:13-27:17, 27:19-27:21, 27:23-27:25, 27:27-27:29, 27:31-27:32, 27:34-27:39, 27:41, 27:43, 27:45-27:48, 27:50-27:52, 27:54-27:56, 27:58-27:61, 27:63-27:64N, 27:66N-27:74, 27:76, 27:78-27:82, 27:84-27:88, 27:90:01-27:94N, 27:96-27:101, 27:103-27:105, 37:01:01, 37:01:03-37:04:02, 37:06-37:09, 37:12-37:13, 37:15-37:33N, 37:35-37:36, 38:17, 40:188, 47:01:01:01-47:01:02, 47:05-47:08, 53:03						6
B*07:36, 07:38, 07:81, 07:180, 08:03, 08:52, 08:78, 15:13:01-15:13:02, 15:16:01-15:17:02, 15:23-15:24:02, 15:67, 15:87, 15:95, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 18:67, 27:02:01-27:02:02, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 27:95, 27:102, 37:34, 38:01:01-38:01:07, 38:05-38:07, 38:09-38:14, 38:16, 38:19-38:22, 38:24-38:28, 38:30-38:34N, 38:36-38:42, 40:13, 40:19, 40:109, 40:117, 44:06, 44:18, 44:25, 44:50, 44:95, 48:18, 49:01:01-49:01:04, 49:03-49:26, 51:01:01-51:01:39, 51:01:41-51:24:04, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:152, 52:01:01:01-52:19, 52:21-52:31, 53:01:01-53:02, 53:04-53:08:02, 53:10, 53:14-53:30, 54:12, 56:21, 57:01:01-57:11, 57:13-57:64, 58:01:01-58:02, 58:04-58:16, 58:18-58:29, 58:31N-58:40, 58:42, 59:01:01:01-59:05	4					
B*07:149, 08:02, 13:01:01-13:04, 13:06-13:08Q, 13:10-13:12:01, 13:13-13:23, 13:25-13:38, 13:40-13:66, 15:36, 15:89, 15:115, 15:256, 18:09, 27:01, 37:10, 38:02:01-38:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 38:43, 40:47, 40:96, 40:110, 40:157, 40:201, 44:02:01:01-44:02:18, 44:02:20-44:05:02, 44:05:04, 44:07-44:08, 44:10, 44:12-44:17, 44:19N-44:24, 44:26-44:45, 44:47-44:49, 44:51-44:74, 44:76-44:89, 44:91-44:94, 44:96-44:128, 44:130, 44:132-44:151, 44:153-44:169, 44:171N-44:175, 47:04, 49:02, 51:54, 51:78:01-51:78:02, 52:20, 53:09, 53:11-53:13, 56:07			3			
A*01:95, 02:81, 02:87, 02:112, 02:124, 02:129, 02:136, 03:152, 23:01:01-23:58, 24:02:01:01-24:03:02, 24:05-24:11N, 24:13:01-24:15, 24:17-24:18, 24:20-24:27, 24:29-24:43, 24:45N-24:64, 24:66-24:88, 24:90N-24:99, 24:101-24:108, 24:110-24:128, 24:130-24:210, 24:212-24:232N, 25:01:01-25:22, 29:13, 31:07-31:08, 31:10, 32:01:01-32:57, 68:36	5					
HLA allele						
Well No.	1	2	3	4	5	6