

Olerup SSP[®] DRB1*03

Product number:	101.113-24/04 – including <i>Taq</i> pol. 101.113-24u/04u – without <i>Taq</i> pol.
Lot number:	36E
Expiry date:	2010-March-01
Number of tests:	24 tests – Product No. 101.113-24 4 tests – Product No. 101.113-04
Number of wells per test:	29
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 36E.

**CHANGES COMPARED TO THE PREVIOUS *OLERUP SSP[®]*
DRB1*03 Lot**

The DRB1*03 specificity and interpretation tables have been updated for the DRB alleles described since the previous *Olerup SSP[®]* DRB1*03 lot was made (**Lot No. X31**).

Five wells have been added to the DRB1*03 kit,
wells **25 to 29**.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
2	-	Removed	Primer moved to new vial 26.
18	-	Modified	Primer modified for the DRB1*030104 allele.
19	Exchanged	-	Primer exchanged for less primer dimer formation tendency.
25	New	New	New primer pair for the DRB1*0332 allele.
26	New	New	Primer pair from vial 2.
27	New	New	New primer pair for the DRB1*0334 allele.
28	New	New	New primer pair for the DRB1*0336 allele.
29	New	New	New primer pair for the DRB1*0339 allele.

Change in revision R01 compared to R00:

- Footnote has been to the Specificity Table for primer mix 28.

PRODUCT DESCRIPTION

DRB1*03 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DRB1*0301 to DRB1*0339 alleles.

PLATE LAYOUT

Each test consists of 29 PCR reactions in a 32 well cut PCR plate. Wells 30 to 32 are empty.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	empty	empty	empty

The 24 well cut PCR plate is marked with 'DRB1*03'.

Well No. 1 is marked with the Lot No. '36E'.

The PCR plates are heat-sealed with a PCR-compatible foil.

Please note: When removing each 32 well PCR plate, make sure that the remaining plates stay sealed. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

The interpretation of DRB1*03 PCR-SSP subtypings will be influenced by one DRB1*01 allele, one, DRB1*04 allele, one DRB1*08 allele, two DRB1*09 alleles, several DRB1*11 alleles, three DRB1*12 alleles, one DRB1*15 allele, one DRB1*16 allele, as well as by many DRB1*13 and DRB1*14 alleles, when present on the other haplotype and also by the DRB3*0107, DRB3*0108, DRB3*0109, DRB3*0204, DRB3*0206, DRB3*0219, DRB3*0220, DRB3*0222 and DRB3*0303 alleles.

UNIQUELY IDENTIFIED ALLELES

All the phenotypically different DRB1*03 alleles, i.e. **DRB1*0301 to DRB1*0339**, recognized by the HLA Nomenclature Committee in January 2008¹ will give rise to unique amplification patterns by the primers in the DRB1*03 subtyping kit.

The DRB1*03 subtyping kit cannot distinguish the DRB1*030101 and 030106 alleles, the DRB1*030201 and 030202 alleles or the DRB1*030501 and 030502 alleles.

¹DRB1 alleles listed on the IMGT/HLA web page 2008-January-11, release 2.20.0, www.ebi.ac.uk/imgt/hla.

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com

RESOLUTION IN HOMO- AND HETEROZYGOTES

The 39 DRB1*03 alleles can be combined in 780 homozygous and heterozygous combinations. 210 of these genotypes do not give rise to unique amplification patterns. The different sizes of the specific PCR products generated by the primers in vials 11, 14, 20 and 23 have not been considered in these calculations.

++++-+-+	+++-----	-+-----++	-----	0302, 0323 = 0319, 0338 = 0323, 0338
++++-+-+	+++-----	-+-----+	-----	0301, 0338 = 0302, 0313 = 0313, 0338
++-+-+--	++-----	-+-+----+	-----	0303, 0325 = 0303, 0328
+--+-+--	+--+-----	-+-----+	-----	0309, 0315 = 0314, 0323
+--+-+--	+--+-----	-+-----+	-----	0301, 0309 = 0305, 0323 = 0309, 0313 =
				0309, 0319 = 0309, 0323
+--+-+--	+-----+	-+-----	-----	0301, 0314 = 0305, 0315 = 0314, 0315
+-----	+--+-----	-+-+----	-----	0304, 0311 = 0304, 0324
+-----	+-----	-+-+----	-----	0301, 0304 = 0304, 0328
+-----	+-----	-+-+----	-----	0308, 0325 = 0308, 0328
+-----	+-----	-+-+----	-----	0307, 0325 = 0307, 0328
+-----	+-----	-+-+----	-----	0301, 0321 = 0307, 0330
+-----	++--+--	-+-----	-----	0301, 0317 = 0317, 0337 = 0335, 0337
+-----	++--+--	-+-----	-----	0306, 0311 = 0311, 0335
+-----	++--+--	-+-----	-+---	0312, 0326 = 0316, 0326
+-----	++--+--	-+-----	-----	0306, 0312 = 0306, 0316
+-----	++-----	-+-+----	-+---	0325, 0326 = 0326, 0328 = 0331, 0333
+-----	++-----	-+-+----	-----	0301, 0331 = 0306, 0325 = 0306, 0328 =
				0328, 0331
+-----	++-----	-+-----	-+---	0301, 0326 = 0306, 0333 = 0326, 0333
+-----	+--+-----	-+-----+	-----	0311, 0323 = 0323, 0324
+-----	+--+-----	-+-----	-----	0311, 0313 = 0313, 0324
+-----	+--+-----	-+-----+	-----	0312, 0323 = 0316, 0323
+-----	+--+-----	-+-----	-----	0312, 0313 = 0313, 0316
+-----	+--+-----	-+-+----	-----	0323, 0325 = 0323, 0328
+-----	+--+-----	-+-+----	-----	0313, 0325 = 0313, 0328
+-----	+--+-----	-+-----+	-----	0301, 0323 = 0313, 0319 = 0313, 0323 =
				0319, 0323 = 0323, 0323
+-----	+--+-----	-+-----	-----	0301, 0313 = 0313, 0313
+-----	+--+-----	-+-+----	-----	0310, 0325 = 0310, 0328
+-----	+--+-----	-+-----	-----	0311, 0312 = 0311, 0316
+-----	+--+-----	-+-----	-----	0311, 0337 = 0324, 0337
+-----	+--+-----	-+-----	-----	0311, 0315 = 0315, 0324
+-----	+--+-----	++-----	-----	0311, 0320 = 0320, 0324
+-----	+--+-----	++-----	-----	0311, 0322 = 0322, 0324
+-----	+--+-----	-+-+----	-----	0311, 0325 = 0311, 0328 = 0324, 0325 =
				0324, 0328
+-----	+--+-----	-+-+----	-----	0311, 0330 = 0324, 0330
+-----	+--+-----	-+-----+	-----	0311, 0318 = 0318, 0324
+-----	+--+-----	-+-----+	-----	0311, 0319 = 0319, 0324
+-----	+--+-----	-+-----	+-----	0311, 0332 = 0324, 0332
+-----	+--+-----	-+-----	-+---	0311, 0333 = 0324, 0333
+-----	+--+-----	-+-----	---+	0311, 0334 = 0324, 0334
+-----	+--+-----	-+-----	-----	0311, 0336 = 0324, 0336
+-----	+--+-----	-+-----	-----	0311, 0339 = 0324, 0339
+-----	+--+-----	-+-----	-----	0301, 0311 = 0301, 0324
+-----	+--+-----	-+-----	-----	0311, 0311 = 0311, 0324

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com

+--+--+--	+-----++-	-+-----	-----	0312, 0337 =	0316, 0337	
+--+--+--	+-----++	-+-----	-----	0312, 0315 =	0315, 0316	
+--+--+--	+-----+-	++-----	-----	0312, 0320 =	0316, 0320	
+--+--+--	+-----+-	-+-----	-----	0312, 0322 =	0316, 0322	
+--+--+--	+-----+-	-+-+----	-----	0312, 0325 =	0312, 0328 =	0316, 0328
+--+--+--	+-----+-	-+----+-	-----	0312, 0318 =	0316, 0318	
+--+--+--	+-----+-	-+-----	-----	0312, 0319 =	0316, 0319	
+--+--+--	+-----+-	-+-----	+-----	0312, 0332 =	0316, 0332	
+--+--+--	+-----+-	-+-----	-+-----	0312, 0333 =	0316, 0333	
+--+--+--	+-----+-	-+-----	--+---	0312, 0334 =	0316, 0334	
+--+--+--	+-----+-	-+-----	---+-	0312, 0336 =	0316, 0336	
+--+--+--	+-----+-	-+-----	-----+	0312, 0339 =	0316, 0339	
+--+--+--	+-----+-	-+-----	-----	0301, 0312 =	0301, 0316 =	0312, 0316
+--+--+--	+-----+-	-+-+----	-----	0325, 0337 =	0328, 0337	
+--+--+--	+-----+-	-+-----	-----	0301, 0337 =	0337, 0337	
+--+--+--	+-----++	-+-+----	-----	0315, 0325 =	0315, 0328	
+--+--+--	+-----++	-+-----	-----	0301, 0315 =	0315, 0315	
+--+--+--	+-----	++-+----	-----	0320, 0325 =	0320, 0328	
+--+--+--	+-----	++-----	-----	0301, 0320 =	0320, 0320	
+--+--+--	+-----	+++-----	-----	0322, 0325 =	0322, 0328	
+--+--+--	+-----	++-----	-----	0301, 0322 =	0322, 0322	
+--+--+--	+-----	-+-++-	-----	0318, 0325 =	0318, 0328	
+--+--+--	+-----	-+-++-	-----	0319, 0325 =	0319, 0328	
+--+--+--	+-----	-+-+----	+-----	0325, 0332 =	0328, 0332	
+--+--+--	+-----	-+-+----	-+-----	0325, 0333 =	0328, 0333	
+--+--+--	+-----	-+-+----	--+---	0325, 0334 =	0328, 0334	
+--+--+--	+-----	-+-+----	---+-	0325, 0336 =	0328, 0336	
+--+--+--	+-----	-+-+----	-----+	0325, 0339 =	0328, 0339	
+--+--+--	+-----	-+-+----	-----	0301, 0325 =	0301, 0328 =	0325, 0328 =
				0328, 0328		
+--+--+--	+-----	-+----+-	-----	0301, 0318 =	0318, 0318	
+--+--+--	+-----	-+----+-	-----	0301, 0319 =	0319, 0319	
+--+--+--	+-----	-+-----	+-----	0301, 0332 =	0332, 0332	
+--+--+--	+-----	-+-----	-+-----	0301, 0333 =	0333, 0333	
+--+--+--	+-----	-+-----	--+---	0301, 0334 =	0334, 0334	
+--+--+--	+-----	-+-----	---+-	0301, 0336 =	0336, 0336	
+--+--+--	+-----	-+-----	-----+	0301, 0339 =	0339, 0339	
+--+--+--	++-+----	-----	-+-----	0324, 0326 =	0326, 0335	
+--+--+--	++-+----	-----	-----	0306, 0324 =	0306, 0335	
+--+--+--	++-+----	-----	-+-----	0306, 0326 =	0326, 0326	
+--+--+--	+-----	-+-+----	-----	0304, 0304 =	0304, 0325	
-+++--++	-+-+----	-----+	-----	0302, 0324 =	0303, 0327	
-+++--++	-+-+----	-+-----+	-----	0302, 0307 =	0303, 0329	
-+++--++	+++-----	-+-----+	-----	0302, 0309 =	0309, 0338	
-+++--++	++-+----	-----+	-----	0327, 0338 =	0335, 0338	
-+++--++	++-+----	-----+	-----	0302, 0338 =	0338, 0338	
-+++--++	-+-+----	-+-----+	-----	0327, 0329 =	0329, 0335	
-+++--++	-+-+----	-----+	-----	0302, 0327 =	0302, 0335	
-+++--++	-+-+----	-----+	-----	0327, 0327 =	0327, 0335	
-+-+----+	-+-+----	-----+	-----	0303, 0324 =	0303, 0335	
--+--+-	-+-+----	-+-----+	-----	0305, 0309 =	0309, 0309	
--+--+-	+-----+	-+-----	-----	0305, 0314 =	0314, 0314	
-----	-+-+----+	-----	-----	0317, 0317 =	0317, 0335	

SPECIFICITY TABLE

DRB1*03 SSP subtyping

Specificities and sizes of the PCR products of the 29 primer mixes used for DRB1*03 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DRB1*03 alleles ³	Other amplified DRB alleles ⁴
1	220 bp	515 bp	030101-030106, 0304, 0306, 0308, 0310-0313, 0315, 0318-0320, 0322, 0323, 0325, 0326, 0328, 0330-0334, 0336, 0337, 0339	0906, 1159, 1327, 1371
2⁵	85 bp	430 bp	030201-0303, 0327, 0329, 0338	1315, 1319, 1326, 1353, 1357, 1402-140302, 1406, 1412, 1413, 1418-1420, 1424, 1427, 1429, 1440, 1441, 1447-1449, 1451, 1463, 1467
3⁶	190 bp	515 bp	030201-030202, 030501-030502, 0309, 0314, 0327, 0329, 0338	1109, 1120, 130201-130203, 130501-130502, 1326, 1329, 1331, 1334, 1336, 1339, 1341, 135002, 1356, 1363, 1365, 1367, 1373, 1374, 1402-140302, 1409, 1413, 1419, 1424, 1427, 1430, 1446-1448, 1451, 1463, 1467, 1608, DRB3*0108, DRB3*0206, DRB3*0220
4	190 bp	430 bp	030101-030106, 0303, 0306-0308, 0310-0313, 0315, 0318-0320, 0322-0324, 0326, 0328, 0332-0334, 0336, 0337, 0339	1116, 1140, 1158, 1159, 130101-130103, 1306, 1309, 1310, 1315, 1318, 1320, 1327, 1328, 1332, 1335, 1340, 1342, 1343, 1351,

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com

				1353, 1359, 1361, 1369, 1371, 1377- 1380, 1406, 1412, 1417, 1418, 1421, 1429, 1433, 1452, 1459, 1464
5⁵	75 bp	430 bp	0304	
6⁵	105 bp	430 bp	030101-030106, 0304-030502, 0308-0315, 0316 ^{weakly} , 0318- 0320, 0322, 0323, 0325, 0328, 0330, 0332-0334, 0336, 0337, 0339	1159, 1327, 1341, 1371
7⁵	135 bp	430 bp	0308	1159
8⁷	185 bp	430 bp	030201-0303, 0307, 0321, 0329, 0338	0116, 0422, 1107, 1525, DRB3*0109, DRB3*0204, DRB3*0219, DRB3*0222, DRB3*0303
9⁵	135 bp	430 bp	030101-030106, 0304-0306, 0309, 0311, 0313 ^{weakly} , 0314-0316, 0318- 0320, 0322, 0323, 0325, 0326, 0328, 0330-0334, 0336, 0337, 0339	1327, 1341, 1371, DRB3*0107
10⁵	145 bp	430 bp	030201-0303, 0306, 0317, 0326, 0327, 0331, 0335, 0338	0820, 1117, 1137, 1141, 1148, 1150, 1152, 1156, 130301- 130302, 130701- 1308, 1312, 1313, 1319, 1326, 1332, 1333, 1336-1338, 1340, 1347-1349, 1353, 1355, 1358, 1360, 1365, 1370, 1372, 1376, 1381, 140101-140301, 140501-1409, 1412- 1414, 1416, 1418- 1420, 1422-1427, 1429, 143201-

				143202, 1434, 1436-1441, 1443-1445, 1447-1449, 1451, 1454-1456, 1458-1460, 1462, 1463, 1467, 1469, 1470, 1474, 1475
11 ^{5,8}	105, 180 bp	515 bp	0309, 0313, 0323, 0338	1339, 1354, 1377-1379, 1448, 1464
12 ⁵	135 bp	430 bp	0310	
13	225 bp	430 bp	0311, 0317, 0324, 0327, 0335	1438
14 ^{5,9}	105, 135 bp	515 bp	0312, 0316	
15	155 bp	430 bp	0317, 0337	1139, 1143, 1150, 1459
16	195 bp	430 bp	0314, 0315	0905, 1159, 1327, 1341, 1371
17	255 bp	515 bp	0320	110601-110602, 1121, 1147, 1157, 1358, 1381, 1429
18 ⁵	120 bp	430 bp	030101-030106, 0304-030502, 0307-0310, 0312-0316, 0318-0323, 0325, 0328-0330, 0332-0334, 0336, 0337, 0339	1107, 1525
19 ⁶	175 bp	515 bp	0322	
20 ^{5,10}	75, 120 bp	430 bp	0304 ^{weakly} , 0325, 0328, 0331	
21 ⁵	70 bp	430 bp	0321, 0330	1110, 111201-111202 ^{weakly} , 1113, 1205, 1214, 1215, 1357, 1364, 1435, 1465, 1472
22 ⁵	95 bp	430 bp	0318	
23 ^{5,11}	65, 165 bp	430 bp	0309, 0319, 0323	
24 ^{5,7}	65 bp	430 bp	030201-0303, 0327, 0329, 0338	1315, 1326, 1353, 1402-140302, 1406, 1412, 1413, 1418, 1419, 1424, 1427, 1429, 1446-1448, 1451, 1452, 1463, 1467, DRB3*0220

25	180 bp	430 bp	0332	
26⁵	115 bp	515 bp	0326, 0333	1144, 1149, 1150
27	230 bp	515 bp	0334	
28¹²	160 bp	430 bp	0336	
29	180 bp	430 bp	0339	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DRB1*03SSP subtypings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective length of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the longer, 515 bp, internal positive control band in order to help in the correct orientation of the DRB1*03 subtyping.

In addition, wells number 3, 11, 14, 17, 19, 26 and 27 contain the primer pair giving rise to the longer, 515 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

³For several DRB alleles only partial second exon nucleotide sequences are available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in the first hyperpolymorphic region of the second exon of DRB alleles are conserved within allelic groups and that unknown sequences of codons 87 to 92 are identical with the DRB1*0101 consensus sequence.

⁴Due to the sharing of sequence motifs within the DR52 group of DRB1 alleles, several DRB1*11, several DRB1*12 as well as many DRB1*13 and DRB1*14 alleles are amplified by some of the DRB1*03 primer mixes. In addition, the DRB1*0116 allele, the DRB1*0422 allele, the DRB1*0820 allele, the DRB1*0905 and 0906 alleles, the DRB1*1525 allele and DRB1*1608 allele are amplified by a few of the DRB1*03 primer mixes. Finally, the DRB3*0107 allele is amplified by primer mix 9, the DRB3*0108 and DRB3*0206 alleles are amplified by primer mix 3, the DRB3*0109, DRB3*0204, DRB3*0219, DRB3*0222 and DRB3*0303 alleles by primer mix 8 and the DRB1*0220 allele by primer mixes 3 and 24.

⁵Specific PCR fragments shorter than 150 base pairs have a lower intensity than longer PCR bands.

⁶Primer mixes 3 and 19 have a tendency of giving rise to primer oligomer formation. Most pronounced for primer mix 19.

⁷The specific primer pairs in primer mixes 8 and 24 yield less product than the other DRB1*03 primer mixes.

⁸Primer mix 11: Specific PCR fragment of 105 bp in the DRB1*0309 and 0323 alleles. Specific PCR fragment of 180 bp in the DRB1*0313 and 0338 and the 1339, 1354, 1377-1379, 1448 and 1464 alleles.

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com

⁹Primer mix 14: Specific PCR fragment of 105 bp in the DRB1*0316 allele. Specific PCR fragment of 135 bp in the DRB1*0312 allele.

¹⁰Primer mix 20: Specific PCR fragment of 75 bp in the DRB1*0304^{weakly}, 0325 and 0331 alleles. Specific PCR fragment of 120 bp in the DRB1*0328 allele.

¹¹Primer mix 23: Specific PCR fragment of 65 bp in the DRB1*0309 and 0323 alleles. Specific PCR fragment of 165 bp in the DRB1*0319 allele.

¹²Primer mix 28 may give to nonspecific amplifications.

INTERPRETATION TABLE**DRB1*03 SSP subtyping****Amplification patterns of the DRB1*0301 to *0339 alleles**

	Well ⁵															
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	220	85	190	190	75	105	135	185	135	145	105	135	225	105	155	195
PCR product											180			135		
Length of int.	515	430	515	430	430	430	430	430	430	430	515	430	430	515	430	430
pos. control ¹																
5'-primer ²	26	13	37	37	26	26	26	26	26	13	13	26	13	26	13	26
	5'-g TA ³ 5'-g TC ³ 5'-Ag A ³ 5'-Ag A ³ 5'-g TA ³ 5'-g TA ³ 5'-g TA ³ 5'-g TT ³ 5'-g TA ³ 5'-g TC ³ 5'-g TC ³ 5'-g TA ³ 5'-g TC ³ 5'-g TA ³ 5'-g TC ³ 5'-g TA ³															
3'-primer(s) ³	86	28	86	86	37	47	58	74	57	47	34	57	74	48	51	77
	5'-C CA ³ 5'-CT C ³ 5'-C AC ³ 5'-C CA ³ 5'-C gg ³ 5'-g gA ³ 5'-C CT ³ 5'-C CC ³ 5'-C AT ³ 5'-g gT ³ 5'-C CC ³ 5'-C Ag ³ 5'-C CT ³ 5'-CCA ³ 5'-CCC ³ 5'-A gg ³															
											60			57		
											5'-A gg ³			5'-gCT ³		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB1 allele																
*030101-030106	1			4		6			9							
*030201-030202		2	3					8		10						
*0303		2		4				8		10						
*0304	1				5	6			9							
*030501-030502			3			6			9							
*0306	1			4					9	10						
*0307				4				8								
*0308	1			4		6	7									
*0309			3			6			9		11					
*0310	1			4		6						12				
*0311	1			4		6			9				13			
*0312	1			4		6								14		
*0313	1			4		6			w		11					
*0314			3			6			9							16
*0315	1			4		6			9							16
*0316						w			9					14		
*0317										10			13		15	
*0318	1			4		6			9							
*0319	1			4		6			9							
*0320	1			4		6			9							
*0321								8								
*0322	1			4		6			9							
*0323	1			4		6			9		11					
*0324				4									13			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

INTERPRETATION TABLE												
DRB1*03 SSP subtyping												
Amplification patterns of the DRB1*0301 to *0339 alleles												
Well ⁵												
17	18	19	20	21	22	23	24	25	26	27	28	29
255	120	175	75	70	95	65	65	180	115	230	160	180
			120			165						
515	430	515	430	430	430	430	430	430	515	515	430	430
13	47	32	26	37	26	26	28	27	13	13	33	27
5'-g TC ^{3'} 5'-g TT ^{3'} 5'-TC g ^{3'} 5'-g TA ^{3'} 5'-g TT ^{3'} 5'-g TA ^{3'} 5'-g TA ^{3'} 5'-g Ag ^{3'} 5'-AC g ^{3'} 5'-g TC ^{3'} 5'-g TC ^{3'} 5'-AAA ^{3'} 5'-C Cg ^{3'}												
85	73	77	37	47	45	34	37	73	38	76	73	73
5'-C Ag ^{3'} 5'-g gC ^{3'} 5'-A gT ^{3'} 5'-gTA ^{3'} 5'-g gA ^{3'} 5'-CCg ^{3'} 5'-C CC ^{3'} 5'-gTT ^{3'} 5'-g gC ^{3'} 5'-g Cg ^{3'} 5'-gTT ^{3'} 5'-g gC ^{3'} 5'-g gC ^{3'}												
			52			67						
			5'-CTT ^{3'}			5'-gAg ^{3'}						
17	18	19	20	21	22	23	24	25	26	27	28	29
	18											
							24					
							24					
	18		w									
	18											
	18											
	18											
	18					23						
	18											
	18											
	18											
	18											
	18											
	18				22							
	18					23						
17	18											
	18			21								
	18	19										
	18					23						
17	18	19	20	21	22	23	24	25	26	27	28	29

Length of spec.	220	85	190	190	75	105	135	185	135	145	105	135	225	105	155	195
PCR product											180			135		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*0325	1					6			9							
*0326	1			4					9	10						
*0327		2	3							10			13			
*0328	1			4		6			9							
*0329		2	3					8								
*0330	1					6			9							
*0331	1								9	10						
*0332	1			4		6			9							
*0333	1			4		6			9							
*0334	1			4		6			9							
*0335										10			13			
*0336	1			4		6			9							
*0337	1			4		6			9						15	
*0338		2	3					8		10	11					
*0339	1			4		6			9							
*0116, 0422								8								
*0820, 1117, 1137, 1141, 1148, 1152, 1156, 130301- 130302, 130701- 1308, 1312, 1313, 1333, 1337, 1338, 1347-1349, 1355, 1360, 1370, 1372, 1376, 140101- 140103, 140501- 140503, 140701- 1408, 1414, 1416, 1422-142302, 1425, 1426, 143201- 143202, 1434, 1436, 1437, 1439, 1443- 1445, 1454-1456, 1458, 1460, 1462, 1469, 1470, 1474, 1475										10						
*0905																16
*0906	1															
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

255	120	175	75	70	95	65	65	180	115	230	160	180	Length of spec.
			120			165							PCR product
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.
	18		20										*0325
									26				*0326
							24						*0327
	18		20										*0328
	18						24						*0329
	18			21									*0330
			20										*0331
	18							25					*0332
	18								26				*0333
	18									27			*0334
													*0335
	18										28		*0336
	18												*0337
							24						*0338
	18											29	*0339
													*0116, 0422
													*0820, 1117, 1137,
													1141, 1148, 1152,
													1156, 130301-
													130302, 130701-
													1308, 1312, 1313,
													1333, 1337, 1338,
													1347-1349, 1355,
													1360, 1370, 1372,
													1376, 140101-
													140103, 140501-
													140503, 140701-
													1408, 1414, 1416,
													1422-142302, 1425,
													1426, 143201-
													143202, 1434, 1436,
													1437, 1439, 1443-
													1445, 1454-1456,
													1458, 1460, 1462,
													1469, 1470, 1474,
													1475
													*0905
													*0906
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.

Length of spec.	220	85	190	190	75	105	135	185	135	145	105	135	225	105	155	195
PCR product											180			135		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*110601-110602, 1121, 1147, 1157																
*1107, 1525								8								
*1109, 1120, 130201- 130203, 130501- 130502, 1329, 1331, 1334, 135002, 1356, 1363, 1367, 1373, 1374, 1430, 1608			3													
*1110, 1113, 1205, 1214, 1215, 1364, 1435, 1465, 1472																
*111201-111202																
*1116, 1140, 1158, 130101-130103, 1306, 1309, 1310, 1318, 1320, 1328, 1335, 1342, 1343, 1351, 1359, 1361, 1369, 1380, 1417, 1421, 1433				4												
*1139, 1143															15	
*1144, 1149																
*1150										10					15	
*1159	1			4		6	7									16
*1315		2		4												
*1319, 1420, 1440, 1441, 1449		2								10						
*1326, 1402, 140301, 1413, 1419, 1424, 1427, 1447, 1451, 1463, 1467		2	3							10						
*1327, 1371	1			4		6			9							16
*1332, 1340				4						10						
*1336, 1365, 1409			3							10						
*1339			3								11					
*1341			3			6			9							16
*1353, 1406, 1412, 1418		2		4						10						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com

255	120	175	75	70	95	65	65	180	115	230	160	180	Length of spec.
			120			165							PCR product
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.
17													*110601-110602, 1121, 1147, 1157
	18												*1107, 1525
													*1109, 1120, 130201- 130203, 130501- 130502, 1329, 1331, 1334, 135002, 1356, 1363, 1367, 1373, 1374, 1430, 1608
				21									*1110, 1113, 1205, 1214, 1215, 1364, 1435, 1465, 1472
				w									*111201-111202
													*1116, 1140, 1158, 130101-130103, 1306, 1309, 1310, 1318, 1320, 1328, 1335, 1342, 1343, 1351, 1359, 1361, 1369, 1380, 1417, 1421, 1433
													*1139, 1143
								26					*1144, 1149
								26					*1150
													*1159
							24						*1315
													*1319, 1420, 1440, 1441, 1449
							24						*1326, 1402, 140301, 1413, 1419, 1424, 1427, 1447, 1451, 1463, 1467
													*1327, 1371
													*1332, 1340
													*1336, 1365, 1409
													*1339
													*1341
							24						*1353, 1406, 1412, 1418
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.

Length of spec.	220	85	190	190	75	105	135	185	135	145	105	135	225	105	155	195
PCR product											180			135		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*1354											11					
*1357		2														
*1358, 1381										10						
*1377-1379, 1464				4							11					
*140302		2	3													
*1429		2		4						10						
*1438										10			13			
*1446			3													
*1448		2	3							10	11					
*1452				4												
*1459				4						10					15	
DRB1 allele																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB3*0107									9							
DRB3*0108, 0206			3													
DRB3*0109, 0204, 0219, 0222, 0303							8									
DRB3*0220			3													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 430 base pairs, for most wells, or a band of 515 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the longer, 515 bp, internal positive control band in order to help in the correct orientation of the DRB1*03 subtyping.

In addition, wells number 3, 11, 14, 17, 19, 26 and 27 contain the primer pair giving rise to the longer, 515 bp, internal positive control band in order to allow kit identification.

²The codon, in the 2nd exon, matching the specificity-determining 3'-end of the primer is given. Codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given. Empty spaces indicate codon boundaries.

³The codon, in the 2nd exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given. Empty spaces indicate codon boundaries.

255	120	175	75	70	95	65	65	180	115	230	160	180	Length of spec.
			120			165							PCR product
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.
													*1354
				21									*1357
17													*1358, 1381
													*1377-1379, 1464
							24						*140302
17							24						*1429
													*1438
							24						*1446
							24						*1448
							24						*1452
													*1459
													DRB1 allele
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.
													DRB3*0107
													DRB3*0108, 0206
													DRB3*0109, 0204,
													0219, 0222, 0303
							24						DRB3*0220
17	18	19	20	21	22	23	24	25	26	27	28	29	Well No.

⁴Primer mix 11: Specific PCR fragment of 105 bp in the DRB1*0309 and 0323 alleles. Specific PCR fragment of 180 bp in the DRB1*0313 and 0338 and the 1339, 1354, 1377-1379, 1448 and 1464 alleles.

Primer mix 14: Specific PCR fragment of 105 bp in the DRB1*0316 allele. Specific PCR fragment of 135 bp in the DRB1*0312 allele.

¹⁰Primer mix 20: Specific PCR fragment of 75 bp in the DRB1*0304^{weakly}, 0325 and 0331 alleles. Specific PCR fragment of 120 bp in the DRB1*0328 allele.

Primer mix 23: Specific PCR fragment of 65 bp in the DRB1*0309 and 0323 alleles. Specific PCR fragment of 165 bp in the DRB1*0319 allele.

“w”, may be weakly amplified.

101.113-24/04 – including *Taq* polymerase
101.113-24u/04u – without *Taq* polymerase

Lot No.: **36E**

Lot-specific Information

www.olerup.com

CERTIFICATE OF ANALYSIS

Olerup SSP® DRB1*03 SSP

Product number: 101.113-24/04 – including *Taq* pol.
101.113-24u/04u – without *Taq* pol.
Lot number: 36E
Expiry date: 2010-March-01
Number of tests: 24 tests – Product No. 101.113-24
4 tests – Product No. 101.113-04
Number of wells per test: 29

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2007-296-01	13	2007-296-13	25	2008-441-25
2	2008-441-02	14	2007-296-14	26	2008-441-26
3	2007-296-03	15	2007-296-15	27	2008-441-27
4	2007-296-04	16	2007-296-16	28	2008-441-28
5	2007-296-05	17	2007-296-17	29	2008-441-29
6	2007-296-06	18	2008-441-18		
7	2007-296-07	19	2008-441-19		
8	2007-296-08	20	2007-296-20		
9	2007-296-09	21	2007-296-21		
10	2007-296-10	22	2007-296-22		
11	2007-296-11	23	2007-296-23		
12	2007-296-12	24	2007-296-24		

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 5, 7, 11, 12, 13, 14, 15, 16, 17, 19, 20, 21, 22, 23, 25, 26, 27, 28 and 29 were available. The specificities of the primers in primer solutions 5, 7, 11, 12, 13, 14, 15, 16, 17, 21 and 26 were tested by separately adding one additional 5'-primer, respectively one additional 3'-primer. In primer solutions 11 and 14 one of the 3'-primers were not possible to test. In primer solution 21 one of the 5'-primers was not possible to test. In primer solutions 19, 25, 28 and 29 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 20, 22, 23 and 27 it was only possible to test the 5'-primers, the 3'-primers were not possible to test.

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2008-April-01

Approved by:

Quality Control, Supervisor

Declaration of Conformity

Product name: Olerup SSP® DRB1*03
Product number: 101.113-24/04, 101.113-24u/04u
Lot number: 36E
Intended use: DRB1*03 high resolution histocompatibility testing

Manufacturer: Olerup SSP AB
Hasselstigen 1
SE-133 33 Saltsjöbaden, Sweden
Phone: +46-8-717 88 27
Fax: +46-8-717 88 18

We, Olerup SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2000 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at Olerup SSP AB, Hasselstigen 1, SE-133 33 Saltsjöbaden, Sweden.

The Authorized Representative located within the Community is: Olerup SSP AB.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Saltsjöbaden, Sweden
2008-April-01

Olle Olerup
Managing Director

101.113-24/04 – including *Taq* polymerase

101.113-24u/04u – without *Taq* polymerase

Lot No.: **36E**

Lot-specific Information

www.olerup.com

101.113-24/04 – including *Taq* polymerase

101.113-24u/04u – without *Taq* polymerase

Lot No.: **36E**

Lot-specific Information

www.olerup.com

101.113-24/04 – including *Taq* polymerase101.113-24u/04u – without *Taq* polymeraseLot No.: **36E**

Lot-specific Information

www.olerup.com**ADDRESSES:****Manufacturer:****Olerup SSP AB**, Hasselstigen 1, SE-133 33 Saltsjöbaden, Sweden.**Tel:** +46-8-717 88 27**Fax:** +46-8-717 88 18**E-mail:** info-ssp@olerup.com**Web page:** <http://www.olerup.com>**Distributed by:****Olerup GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.**Tel:** +43-1-710 15 00**Fax:** +43-1-710 15 00 10**E-mail:** support-at@olerup.com**Web page:** <http://www.olerup.com>**Olerup Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382**Tel:** 1-877-OLERUP1**Fax:** 610-344-7989**E-mail:** info.us@olerup.com**Web page:** <http://www.olerup.com>For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.