

INTERPRETATION TABLE				
HLA-B*82 SSP subtyping				
Amplification patterns of the B*82 alleles				
	Tube			
	1	2	3	4
Length of spec.	195	135	230	195
PCR product				
Length of int.	800	1070	1070	1070
pos. control¹				
5'-primer²	420	557	368	105
	5'-TTA ^{3'}	5'-ggA ^{3'}	5'-gTT ^{3'}	5'-gCT ^{3'}
3'-primer³	572	3rd I	557	259
	5'-gCg ^{3'}	5'-TAT ^{3'}	5'-ggC ^{3'}	5'-gTT ^{3'}
Tube No.	1	2	3	4
HLA-B allele				
*8201	+	+		+
*8202	+		+	+
*1506, 152701-152702, 1584, 9509			+	
*4410, 4415, 4418, 4501, 4505-4507, 5002	+			
HLA-B allele				
Tube No.	1	2	3	4
Cw*0408			+	
Tube No.	1	2	3	4