

INTERPRETATION TABLE

HLA-B*5701 SSP typing

	Well ⁴							
	1	2	3	4	5	6	7	8
Length of spec.	90	220	95	105	90	90	165	195
PCR product			180		250	210	215	
Length of int.	800	800	800	1070	800	1070	1070	1070
pos. control ¹								
5'-primer(s) ²	209	362	362	209	209	320	362	103
	5'-ggC ^{3'}	5'-ggT ^{3'}	5'-ggT ^{3'}	5'-ggC ^{3'}	5'-ggA ^{3'}	5'-CCC ^{3'}	5'-ggT ^{3'}	5'-CCT ^{3'}
					362	362		
					5'-ggT ^{3'}	5'-ggT ^{3'}		
3'-primer(s) ³	256	539	418	271	256	2 nd I	486	256
	5'-CCC ^{3'}	5'-TCA ^{3'}	5'-gTC ^{3'}	5'-CAC ^{3'}	5'-CCC ^{3'}	5'-TCg ^{3'}	5'-gCg ^{3'}	5'-CCC ^{3'}
			500		572	412	538	
			5'-ggA ^{3'}		5'-gCg ^{3'}	5'-gTT ^{3'}	5'-gTC ^{3'}	
Well No.	1	2	3	4	5	6	7	8
HLA-B allele								
*570101-570103	1	2						
*5702, 5712	1					6		
*570301-570302, 5708, 5717	1	2				6		
*5704	1		3					
*5705, 5711, 5713	1							
*5706	1	2	3					
*5707	1	2			5	6		
*5709	1					6	7	
*5710	1	2						8
*5714	1	2					7	
*5715	1	2		4				
*5716	1	2			5			
*4030, 4034		2				6		
*5514, 5814		2					7	
HLA-B allele								
Well No.	1	2	3	4	5	6	7	8