

INTERPRETATION TABLE		
DRB4*01030102N SSP typing		
	Well	
	1	2
Length of spec.	150	245
PCR product		
Length of int.	515	430
pos. control ¹		
5'-primer(s) ²	1 st ⁴	1 st ⁵
	5'-CAA ^{3'}	5'-ggg ^{3'}
3'-primer(s) ³	42	5 ⁶
	5'-TC A ^{3'}	5'-Tg C ^{3'}
Well No.	1	2
DRB4 allele		
*01030102N	1	
*01010101, 0102, 01030101, 010302, 010303, 010304, 0104, 0106, 0107, 0201N		2
*0105	?	
DRB4 allele		
Well No.	1	2