

INTERPRETATION TABLE		
DRB5*0108N SSP typing		
	Well	
	1	2
Length of spec.	195	260
PCR product		
Length of int.	515	430
pos. control¹		
5'-primer(s)²	108	13
	5'-Ag A ^{3'}	5'-g TA ^{3'}
3'-primer(s)³	160	87
	5'-CAT ^{3'}	5'-CTC ^{3'}
Well No.	1	2
DRB5 allele		
*0108N	1	2
*010101-0107, 0109-0113, 0202-0205		2
DRB5 allele		
Well No.	1	2