

[illegible]

[illegible]

*03:04:26, 03:18, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182				76															93					
*03:07, 03:45, 03:163			75	76						83				87					93	94				
*03:10, 03:29			75	76										87					93	94				
*03:13, 03:97, 03:132, 03:167, 03:177				76	77															94				
*03:15			75	76						83				87						94				
*03:16, 03:110			75																	94				
*03:58					77									87				92	93	94				
*03:67				76	w														93	w				
*03:81, 03:175				76	77											89			93	94				
*03:86					77									87				92		94				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ⁴																								
*03:94			75												87				92		94			
*03:99															87				92					
*03:130, 03:140			75	76						83											93	94		
*03:151, 03:165					77																93	94		
*04:01:01:01-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:51, 04:04:01-04:05, 04:07-04:09N, 04:12-04:13, 04:15:01-04:15:03, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:34-04:35, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:102, 04:104-04:106, 04:109-04:111, 04:113, 04:115N, 04:117-04:128, 04:131-04:139, 04:141-04:145						78				83			85					90			93			
*04:01:16, 04:23, 04:54, 04:72						78				83			85					90						
*04:01:23, 04:33, 04:116, 04:130						78				83			85								93			
*04:01:28, 04:140										83			85					90			93			
*04:03:01, 04:06, 04:80			75			78				83			85		87			90			93			
*04:10						78				83								90			93			
*04:11						78										88					93			
*04:14, 04:68						78				83			85					90			93		95	
*04:16			75							83			85		87			90			93			
*04:29						78							85			88		90			93			
*04:32, 04:77		74				78				83			85					90			93	94		
*04:36, 04:55						78										88		90			93			
*04:37						78				83			85		87			90		92	93			
*04:42										83			85	86				90			93		95	
*04:58						78				83			85					90	91		93			
*04:71	73					78				83			85					90			93			
*04:103										83			85		87			90			93			
*04:107			75			78				83			85		87						93			

Negative Control

*04:108						78					83		85				89	90							
*04:112						78					83		85				89	90			93				
*04:114						78					83		85			88		90			93				
*04:129						78	79				83		85					90			93				
*05:01:01:01-05:01:26, 05:03, 05:05-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:58:02-05:61, 05:63-05:67, 05:69-05:77, 05:80-05:84, 05:86-05:91N							79				83			87							93				
*05:04:01-05:04:02, 05:22:01-05:22:02							79				83			87											
*05:16, 05:85							79				83			87					92	93					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
HLA-C allele ⁴																									
*05:17, 05:25, 05:55, 05:68, 05:79							79				83			87				90			93				
*05:20, 05:43							79				83										93				
*05:23, 05:62							79				83			87				91			93				
*05:26						78	79				83			87							93				
*05:36							79				83			87			89				93				
*05:42							79				83		85	87			90				93				
*05:46							79				83		85	87							93				
*05:58:01			75				79				83			87							93				
*05:78						78	79				83										93				
*06:02:01:01-06:02:01:02, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:26, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98								80			83			86										95	
*06:02:12, 06:11, 06:52								80						86										95	
*06:02:16, 06:04, 06:40, 06:79N, 06:81								80			83													95	
*06:03:01			75					80			83			86	87						94				
*06:03:02, 06:76:01								80			83			86	87										
*06:05							79	80			83						90		92					95	
*06:06								80			83								92					95	
*06:08		74						80			83			86										95	
*06:09								80			83			86							93			95	
*06:14, 06:35, 06:72								80			83										93			95	
*06:17											83			86					92					95	
*06:31								80			83			86			90		92					95	
*06:32, 06:94-06:95, 06:99											83			86										95	
*06:67							79	80			83			86										95	

Negative Control

*06:76:02										83	84	85	86	87		90								
*06:77								80		83			86											
*06:89								80		83			86			89					95			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ⁴																								
*07:01:01:01-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:03-07:06, 07:08, 07:11-07:25, 07:27:01-07:27:02, 07:29-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153-07:155, 07:157-07:160, 07:162-07:172, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:266, 07:268-07:293, 07:297-07:307									81													95		
*07:02:08, 07:02:29, 07:26, 07:33N, 07:39-07:40, 07:92, 07:139, 07:152N, 07:156, 07:246, 07:267, 07:296									81															
*07:02:09						78			81				85				90						95	
*07:02:41, 07:223									81														w	
*07:07, 07:09, 07:49, 07:76, 07:247									81		83								92				95	
*07:10, 07:28, 07:41, 07:43, 07:184, 07:196									81												93		95	
*07:31:01-07:31:02									81					86			90						95	
*07:64									81						88									
*07:96:01-07:96:02			75						81													94		
*07:101, 07:148, 07:161		74							81	82													95	
*07:123, 07:173, 07:294									81								89						95	
*07:125						78			81				85	86									95	
*07:177									81					86			90							
*07:210, 07:238									81		83								92					
*07:242				76					81														95	
*08:01:01-08:01:09, 08:03:01-08:03:02, 08:06, 08:08, 08:14, 08:16:01-08:16:02, 08:20, 08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01											82							90			93			

Negative Control

*08:02:01-08:02:10, 08:04:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62-08:63, 08:66-08:71, 08:73-08:77										82										93				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ⁴																								
*08:05, 08:25										82			85								93			
*08:07, 08:47										82								91			93			
*08:09, 08:11										82								90						
*08:10							79			82	83				87			90			93			
*08:12										82										92	93			
*08:21										82			85					90			93			
*08:31		74								82											93			
*08:48										82														
*12:02:01-12:02:09, 12:10:01-12:10:02, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86												84	85											
*12:03:01:01-12:03:01:02, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:17-12:03:19, 12:03:21-12:03:22, 12:03:25, 12:06-12:07, 12:11-12:13, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82, 12:87-12:95												84	85	86										
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:26, 12:59													85	86										
*12:03:04												84		86										
*12:03:09								80				84		86										
*12:03:13, 12:03:24												84	w	86										
*12:03:23			75									84	85	86										
*12:04:01											83			86	87									
*12:04:02, 12:60											83	84	85	86	87									
*12:05											83	84		86	87									
*12:08												84	85				89							
*12:09											83				87					92				
*12:14:01													85					90	91					
*12:14:02												84	85					90	91					
*12:15								80				84	85	86			89							
*12:16												84	85									95		

Negative Control

*12:18:02, B*67:02 ⁷													85													
*12:21							79				83	84			87											
*12:24												84	85						92							
*12:28, 12:58												84	85	86				90								
*12:31												84	85	86								93				
*12:33							79				83			86	87											
*12:41											83	84	85		87											
*12:44													84	85								93				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
HLA-C allele ⁴																										
*12:49													84													
*12:54												83	84	85	86											
*12:55													84	85	86		88									
*12:72													84	85		87										
*12:81													84	85	86			89								
*14:02:01-14:02:07, 14:02:09-14:03, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:24:02, 14:26-14:27, 14:29-14:44, 14:46-14:48, 14:50-14:51																88				92						
*14:02:08													84			88				92						
*14:04												83				88				92						
*14:06, 14:08-14:09, 14:28, A*30:62 ⁸																88										
*14:10																88		90		92						
*14:12, 14:49												83								92						
*14:15														86		88		90								
*14:16														86		88										
*14:17																88			91	92						
*14:25		74														88				92		94				
*14:45																87	88			92						
*15:02:01-15:02:09, 15:02:11-15:02:15, 15:08, 15:10:02, 15:13, 15:15, 15:17-15:18, 15:26, 15:28, 15:31-15:35, 15:38-15:39, 15:41, 15:44-15:45, 15:47-15:53, 15:56-15:57, 15:60, 15:64			75					80				83			87		89				93					
*15:02:10			75					80				83			87		89				93	94				
*15:03, 15:16			75					80				83		85	87		89				93					
*15:04-15:06:03, 15:09, 15:19, 15:24, 15:27, 15:29-15:30, 15:40, 15:46, 15:54, 15:59, 15:66			75									83			87		89				93					
*15:07			75					80									89				93					
*15:10:01								80				83			87		89				93					
*15:11			75			78		80				83			87		89				93					
*15:12					77			80				83			87		89	90			93					
*15:21			75					80							87		89				93					
*15:22, 15:55			75									83			87		89									
*15:23												83					89				93					

Negative Control

*15:25			75											85					90	91		93		
*15:36						78						83							89			93		
*15:37, 15:58			75					80				83				87			89					
*15:42		74	75					80				83				87			89			93		
*15:43			75					80											89			93	94	
*15:61												83				87			89			93		
*15:62			75					80				83				87			89	90		93		
*15:63								80				83							89			93		
*15:65			75									83				87			89	90	91			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ⁴																								
*16:01:01-16:01:13, 16:06-16:08, 16:10-16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:01, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56														85					90					95
*16:02:01-16:02:09, 16:09, 16:12, 16:19, 16:25, 16:46-16:47											83		85		87			90						95
*16:04:01, 16:33, 16:42, 16:55														86				90		92				
*16:15:02												84	85					90						95
*16:18, 16:37													85		87			90						95
*16:20													85					89	90					95
*16:29		74												86				90		92				
*16:34			75			78							85					90						95
*16:35														86				90						
*16:40						79										88								
*16:43													85											95
*16:45													85					90			93			95
*16:48											83			86	87			90						
*16:50		74											85					90						95
*16:53						79										88				92				
*16:57											83		85					90						95
*17:01:01-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16, 17:18-17:19		74									83				87				91		93			
*17:01:04		74									83		85		87				91		93			
*17:05		74				79					83				87				91		93			
*17:07, 17:17											83				87				91		93			
*18:01-18:02, 18:04-18:06											83									92	93			95
*18:03		74									83									92	93	94		95
A*24:96, A*24:106, A*24:174, B*08:16, B*46:25																								95
B*14:03																				92				
B*35:222															87									
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

