

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01-01:02:17, 01:02:19-01:02:36, 01:06-01:07:01, 01:08, 01:11-01:13, 01:15-01:20, 01:23, 01:25-01:33, 01:37N-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:80-01:89N, 01:91-01:100, 01:103-01:107, 01:109N-01:111N	73																			92				
*01:02:18	73											84								92				
*01:03, 01:24, 01:34, 01:78	73																			92	93			
*01:04	73											84		86										
*01:05, 01:07:02, 01:22, 01:35-01:36, 01:79:01-01:79:02, B*39:76, B*54:18 ³	73																							
*01:09	73													86										
*01:10	73	74																		92				
*01:14, 01:59	73											83								92				
*01:21	73											84	85											
*01:43	73	74									82									92				
*01:49:01-01:50, 01:55	73																				93			
*01:60	73																		91	92				
*01:67, 14:58 ³	73															88				92				
*01:90, 01:101-01:102	73							80									89			92				
*02:02:01		74	75			78					83				w									
*02:02:02:01-02:02:03, 02:02:05-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38N, 02:40:01-02:40:02, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106		74	75			78					83				87									
*02:02:13, 02:02:26, 02:21, 02:26:02-02:26:03, 02:39, 02:45		74	75								83				87									
*02:02:29		74	75								83				87						94			
*02:02:30		74				78					83				87									
*02:04		74	w			78					83				87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06, 02:23, 02:36		74	75			78		80			83				87		89							
*02:12		74	75			78							85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01		74	w			78					83													
*02:14:02		74				78					83													
*02:22, 02:62		74	75			78					83				87						93			
*02:26:01		74									83				87									
*02:27:01-02:27:02		74	75			78																		
*02:47		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87			90						
*02:55:01		74	75			78					83		85		87									
*02:58			75			78					83				87									
*02:65		74	75			78					83													
*02:68		74	75			78					83				87		89							
*02:87		74								82														
*02:94		74	75			78	79				83				87									
*02:101		74	75			w					83				87									
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:02:14, 03:05, 03:17, 03:25, 03:27, 03:33, 03:35-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296			75	76																	94			

Negative Control

*03:02:10, 03:264				76																				
*03:02:13, 03:16, 03:110				75																		94		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*03:03:01-03:03:20, 03:03:22-03:03:29, 03:11:01-03:11:02, 03:18:02, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291				76	77																93	94		
*03:03:21, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195				76																	93			
*03:04:01:01-03:04:06, 03:04:08-03:04:24, 03:04:27-03:04:48, 03:06:01-03:06:02, 03:08-03:09, 03:14, 03:19, 03:23-03:24, 03:26, 03:28, 03:32, 03:34, 03:37:01-03:38:02, 03:40:01-03:40:04, 03:42, 03:44, 03:46-03:48, 03:51, 03:54, 03:57, 03:63-03:64:01, 03:65, 03:70, 03:72-03:74, 03:77-03:78, 03:82, 03:87:01-03:87:02, 03:90-03:93, 03:98, 03:100-03:101, 03:104-03:107, 03:109, 03:111, 03:114-03:115, 03:117-03:118, 03:123, 03:125, 03:128-03:129, 03:131, 03:134, 03:136, 03:138, 03:143, 03:147-03:149, 03:153, 03:157, 03:159, 03:162, 03:164, 03:173-03:174, 03:179-03:181, 03:183-03:184:02, 03:186:01-03:186:02, 03:191, 03:193, 03:208N, 03:210-03:213, 03:215, 03:218-03:219, 03:232-03:236, 03:239-03:240, 03:244Q, 03:246-03:247, 03:250, 03:252, 03:255-03:257, 03:259, 03:261, 03:263, 03:265N-03:266, 03:269-03:270, 03:277N, 03:279, 03:282-03:283, 03:287, 03:294				75	76																93	94		
*03:04:07, 03:41:01-03:41:02, 03:80, 03:113, 03:260, 03:286, 03:293				75																	93	94		
*03:04:25, 03:145, 03:172				75	76																	93		
*03:07, 03:45, 03:163, 03:297				75	76						83					87						93	94	
*03:10, 03:29, 03:209				75	76											87						93	94	
*03:13:01-03:13:02, 03:97, 03:132, 03:167, 03:177, 03:271				76	77																		94	
*03:15				75	76						83					87							94	
*03:58					77											87				92	93	94		
*03:67, 03:295				76	w																93	w		
*03:81, 03:175				76	77												89					93	94	
*03:86					77											87				92		94		
*03:94				75												87				92		94		
*03:99																87				92				
*03:130, 03:140				75	76						83										93	94		
*03:151, 03:165, 03:251, 03:267					77																93	94		
*03:155				75	76												88				93	94		
*03:199, 03:245				75	76												89					94		
*03:231					76		78										88				93			
*03:243				76	77						83										93	94		
*03:248				75	76																			

*03:262				76	77																?	94		
*03:268				76	77						83				87						93	94		
*03:278			w	76																	93	w		
*03:280				75	76						82										93	94		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*04:01:01:01-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:60, 04:01:62, 04:01:64, 04:01:66-04:01:68, 04:04:01-04:05, 04:07-04:09N, 04:12-04:13, 04:15:01-04:15:03, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:34-04:35, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:102, 04:104-04:106, 04:109-04:111, 04:113, 04:115N, 04:117-04:128, 04:131-04:139, 04:141-04:146, 04:148-04:152, 04:154-04:159, 04:161-04:168, 04:170N, 04:173N-04:177, 04:179-04:189, 04:191N, 04:193, 04:195, 04:197-04:206						78					83		85				90			93				
*04:01:16, 04:23, 04:54, 04:72						78					83		85				90							
*04:01:23, 04:33, 04:116, 04:130						78					83		85							93				
*04:01:28, 04:01:61, 04:01:63, 04:01:65, 04:140, 04:196											83		85				90			93				
*04:03:01-04:03:02, 04:06, 04:80, 04:147, 04:171			75			78					83		85		87		90			93				
*04:10, 04:153, 04:192, 04:194						78					83						90			93				
*04:11						78									88					93				
*04:14, 04:68						78					83		85				90			93		95		
*04:16			75								83		85		87		90			93				
*04:29						78							85		88		90			93				
*04:32, 04:77			74			78					83		85				90			93	94			
*04:36, 04:55						78									88		90			93				
*04:37						78					83		85		87		90		92	93				
*04:42:01-04:42:02											83		85	86			90			93		95		
*04:58						78					83		85				90	91		93				
*04:71	73					78					83		85				90			93				
*04:103											83		85		87		90			93				
*04:107			75			78					83		85		87					93				
*04:108						78		w			83		85				89	90						
*04:112						78					83		85				89	90		93				
*04:114						78					83		85		88		90			93				
*04:129						78	79				83		85				90			93				
*04:160			75			78					83		85		87		90	91		93				
*04:169						78					83						89			93				
*04:172						78							85		88					93				
*04:178						78		80			83		85				89	90		93				
*04:190						78					83		85		87		90			93				
*05:01:01:01-05:01:31, 05:03, 05:05-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:58:02-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106, 05:108-05:111, 05:114, 05:116-05:122						79					83				87					93				
*05:04:01-05:04:02, 05:22:01-05:22:02, 05:103						79					83				87									
*05:16, 05:85						79					83				87				92	93				
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115						79					83				87		90			93				
*05:20, 05:43, 05:113N						79					83									93				
*05:23, 05:62						79					83				87			91		93				
*05:26						78	79				83				87					93				
*05:36							79				83				87		89			93				
*05:42							79				83		85		87			90		93				

Negative Control

*05:46								79			83		85	87					93						
*05:58:01			75					79			83			87					93						
*05:78							78	79			83								93						
*05:105		74						79			83			87					93						
*05:107								79			83			87				92							
*05:112								79			83		85	87					93		95				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*06:02:01:01-06:02:01:03, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:01, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:164								80			83			86									95		
*06:02:12, 06:11, 06:52, 06:122, 06:124, 06:147								80						86										95	
*06:02:16, 06:02:43, 06:04:01, 06:40, 06:79N, 06:81, 06:101, 06:123, 06:127:01-06:127:02, 06:136, 06:153								80			83													95	
*06:02:38								80			83	84		86										95	
*06:02:39, 06:32, 06:94-06:95, 06:99											83			86										95	
*06:03:01, 06:132:01			75					80			83			86	87							94			
*06:03:02, 06:76:01								80			83			86	87										
*06:05							79	80			83						90		92					95	
*06:06								80			83								92					95	
*06:08		74						80			83			86										95	
*06:09								80			83			86							93			95	
*06:14, 06:35, 06:72, 06:143-06:144								80			83										93			95	
*06:17											83			86					92					95	
*06:31								80			83			86			90		92					95	
*06:67							79	80			83			86										95	
*06:76:02											83	84	85	86	87		90								
*06:77								80			83			86											
*06:89								80			83			86			89							95	
*06:118								80			83			86			90	91						95	
*06:132:02			75								83			86	87							94			
*07:01:01:01-07:01:40, 07:01:42-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:06, 07:08, 07:11-07:25, 07:27:01-07:27:02, 07:29:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:266, 07:268-07:293, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368-07:401, 07:404, 07:407-07:414, 07:417-07:446									81														95		
*07:01:41, 07:02:41, 07:223, 07:415-07:416									81															w	
*07:02:08, 07:02:29, 07:26, 07:33N, 07:39-07:40, 07:92, 07:139, 07:152N, 07:156, 07:246, 07:267, 07:296, 07:314:01, 07:317, 07:322, 07:335, 07:351, 07:405									81																
*07:02:09						78			81				85				90							95	

Negative Control

Negative Control

*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406										81	83								92			95		
*07:10, 07:28, 07:41, 07:43, 07:184, 07:196, 07:367										81										93		95		
*07:31:01-07:31:02										81				86				90				95		
*07:64, 07:402										81					88									
*07:96:01-07:96:02, 07:314:02			75							81											94			
*07:101, 07:148, 07:161		74								81	82											95		
*07:123, 07:173, 07:294										81							89					95		
*07:125, 07:356						78				81				85	86							95		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*07:154									81								90						95	
*07:177									81					86			90							
*07:210, 07:238									81		83								92					
*07:242				76					81														95	
*07:316, 07:338	73								81														95	
*08:01:01-08:01:17, 08:01:19, 08:03:01-08:03:03, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117, 08:119, 08:121N-08:122										82							90				93			
*08:01:18, 08:02:01:01-08:02:12, 08:04:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62-08:63, 08:66-08:71, 08:73-08:77, 08:90, 08:92-08:94, 08:100, 08:103, 08:107, 08:110-08:112, 08:114-08:116, 08:120, 08:123										82											93			
*08:05, 08:25										82			85								93			
*08:07, 08:104										82							91				93			
*08:09, 08:11, 08:83										82							90							
*08:10							79			82	83			87			90				93			
*08:12										82									92		93			
*08:21										82			85				90				93			
*08:31		74								82											93			
*08:47										82							90	91			93			
*08:48, 08:108										82														
*08:95										82					88		90				93			
*08:113								80		82							89							
*12:02:01-12:02:10, 12:02:12-12:02:13, 12:10:01-12:10:02, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N, 12:151, 12:153, 12:155Q, 12:161-12:162												84	85											
*12:02:11								80				84												

Negative Control

Negative Control

[illegible]

*14:10																		88		90		92				
*14:12, 14:49											83											92				
*14:15																		86		88		90				
*14:16																		86		88						
*14:17																				88		90	91	92		
*14:25			74																	88				92		94
*14:45																			87	88				92		
*14:54																				88					93	
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:20, 15:02:22-15:02:24, 15:06:01-15:06:03, 15:08-15:09, 15:10:02-15:10:03, 15:13, 15:15, 15:17-15:19, 15:26-15:28, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44-15:45, 15:47-15:53, 15:56-15:57, 15:60, 15:64, 15:67-15:68, 15:73-15:74, 15:78-15:84Q, 15:86-15:89, 15:91-15:101, 15:103, 15:106, 15:109, 15:112-15:114				75								80												93		
*15:02:10, 15:02:17				75								80								87		89			93	94
*15:02:21, 15:10:01, 15:61												80								87		89			93	
*15:03, 15:16				75								80								83		85	87		89	
*15:04:01-15:04:02, 15:24, 15:30, 15:40, 15:66, 15:76-15:77				75																83			87		89	
*15:05:01-15:05:08, 15:05:10, 15:29, 15:46, 15:54, 15:59, 15:69-15:70, 15:90, 15:105Q, 15:108, 15:110, 15:115N, 15:117				75									w												93	
*15:05:09													w							83			87		89	
*15:05:11, 15:116				75									w										87		89	
*15:07				75																				89		93
*15:11				75									80												89	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
*15:12						77		80				83				87		89	90				93			
*15:21, 15:71				75				80								87		89					93			
*15:22				75				w				83				87		89								
*15:23								w				83						89						93		
*15:25				75										85					90	91				93		
*15:31, 15:62				75				80				83				87		89	90					93		
*15:36						78		w				83						89						93		
*15:37, 15:55, 15:58, 15:72				75				80				83				87		89								
*15:42			74	75				80				83				87		89						93		
*15:43				75				80										89						93	94	
*15:63								80				83						89						93		
*15:65				75								83				87		89	90	91						
*15:75				75				80				83		85		87		89	90					93		95
*15:85				75				80				83						89						93		
*15:102				75								83				87		89					92			
*15:104		73		75				w				83				87		89						93		
*15:107							79	80				83				87		89						93		
*15:111			74	75				w				83				87		89						93		
*16:01:01:01-16:01:15, 16:01:18-16:01:19, 16:06-16:08, 16:10-16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56, 16:58-16:59, 16:62, 16:65, 16:67, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86															85										95	
*16:01:16								80							85				90						95	
*16:01:17, 16:15:02													84	85					90						95	
*16:02:01-16:02:12, 16:09, 16:12, 16:19, 16:25, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84													83		85		87								95	
*16:02:13													83	84	85		87								95	

ontrol

[illegible]