

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ^{1,2}	ser ³																								
*01:01:01-01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:55, 01:01:57-01:01:64, 01:01:67-01:04N, 01:06-01:07, 01:09:01-01:11N, 01:15N-01:18N, 01:20, 01:22N-01:27N, 01:29, 01:31N-01:33, 01:35-01:42, 01:44-01:50, 01:52:01N-01:78, 01:80-01:94, 01:96-01:98, 01:100-01:120, 01:122-01:125, 01:127-01:135, 01:137-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:175, 01:177, 01:180-01:191, 01:193, 01:195-01:199, 01:201-01:204, 01:206-01:207	A1, Null, -	1			4																				
*01:01:23, 01:01:48, 01:01:65-01:01:66, 01:08, 01:14, 01:30, 01:79, 01:99, 01:147Q, 01:149, 01:159, 01:167, 01:178N-01:179N, 01:205, 01:208Q	A1, Null, -	1																							
*01:01:56	-	1		4									w												
*01:12, 01:19, 01:200	A1, -	1	3																						
*01:13, 01:176	A1, -	1		4							11		13		15										
*01:21, 01:126	A1, -	1	3	4																					
*01:28	-	1		4							11				15										
*01:43	-	1		4						10															
*01:51	-	1		4				8		10															
*01:95	-	1		4														17							
*01:121	-	1		4																19					
*01:136, 01:192, 11:94, 11:112, 11:211 ⁴	-	1		4							11														
*01:145	-	1												w											
*01:194	-	1		4									13		15										
*02:01:01-01:02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:16, 02:18-02:22:02, 02:24:01-02:33, 02:36-02:45, 02:47, 02:49-02:54, 02:57-02:61, 02:63-02:64:02, 02:66-02:69, 02:71-02:77, 02:79:01-02:97:02, 02:99, 02:101:01-02:102, 02:104-02:107, 02:109, 02:111-02:128, 02:130-02:134, 02:136-02:168, 02:170-02:236, 02:238-02:242, 02:244-02:267, 02:269-02:299, 02:301N-02:302, 02:304-02:308, 02:310-02:337, 02:339-02:406, 02:408-02:448, 02:450-02:453, 02:455-02:526, 02:528-02:570, 02:572-02:579, 02:581, 02:583-02:616, 02:618Q-02:634	A2, Low A2, -, A203/A2, A210/A2, A2/blank, Null, -	2																							
*02:17:01-02:17:03	A2, -	2					w														21				
*02:34-02:35:01, 02:35:03, 02:56:01-02:56:02, 02:62, 02:103, 02:580	A2, -	2											13							20		22			
*02:35:02	A2												13							20		22			
*02:46, 02:70	A2	2																		20		22			
*02:48, 02:129, 02:571	A2, -																			20		22			
*02:55	A2/A28	2						8													21		23		
*02:65, 02:407, 02:449	A2/A19, -	2																19							
*02:78	-	w		4																20		22			
*02:108, 02:110, 02:268, 02:300, 02:303, 02:617	A2, -	2																			21				
*02:135	-	2						8					13											24	
*02:169	-	2		4																					
*02:237	-	2												14		16							23		
*02:243:01-02:243:02	-	2																18			21				
*02:309, 02:454	-	2										12		14										24	
*02:338	-	2	3																						
*02:527, 02:582	-	2						8																	
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*03:01:01-01:03:01:18, 03:01:20-03:01:22, 03:01:24-03:01:37, 03:01:39-03:07, 03:09-03:11N, 03:13-03:14, 03:16-03:17:02, 03:20-03:23:02, 03:25-03:29, 03:31, 03:33-03:35, 03:37-03:40, 03:42, 03:44-01:03:49, 03:51-03:56, 03:58, 03:60-03:62, 03:64-03:71, 03:73-03:74, 03:76-03:81, 03:83-03:87, 03:90-03:94, 03:96-03:106, 03:109-03:110, 03:112-03:123:01, 03:124-03:134, 03:136-03:141, 03:143-03:151, 03:153-03:170, 03:174-03:175, 03:179-03:185, 03:188-03:193, 03:195-03:197N, 03:200Q-03:202, 03:204, 03:206-03:210, 03:212-03:214, 03:216-03:218, 03:220-03:224, 03:226-03:230, 03:232-03:245, 03:247-03:248	A3, A3/blank, Null, -	3												13											
*03:01:19	-			3									12	13											24
*03:01:23, 03:08, 03:36N, 03:57, 03:59, 03:72, 03:89:01-03:89:02, 03:107-03:108, 03:111, 03:142, 03:172-03:173, 03:176, 03:178N, 03:198-03:199, 03:203, 03:205, 03:211	A3, Null, -			3																					
*03:01:38, 03:123:02, 03:171, 03:231	-			3										13	14										
*03:12	A3			3	4																				
*03:15, 03:19, 03:30	A3			3		5								13											
*03:18, 03:135	-	1		4										13											
*03:24, 03:50	A3, -			3				8						13											

[illegible]

*26:03:01, 26:06, 26:21, 26:111	A26, -								8			11	13										23		
*26:05	A26								8				13										23		
*26:07:01-26:07:02, 26:65	A26, -								8		10		13										24		
*26:09, 26:91	A26, -								8		10		12												
*26:16	A26						7	8		10			13												
*26:19	-			4				8		10			13	14											
*26:22	A26							8		10			13	14					21						
*26:29, 26:49	A26, -							8		10			w												
*26:30	A26							8					13								23	24			
*26:54	-							8		10		12	13				20								
*26:72	-			4				8		10			13												
*26:78	-							8				11									23				
*26:92	-							8				11	13									24			
*26:121, 26:123	-							8		10			13						21						
*26:122	-							8		10			13				20	21							
*29:01:01:01-29:06, 29:08N-29:12, 29:15-29:27, 29:29-29:38, 29:40-29:47, 29:50-29:66, 29:68-29:90	A29, Null, -													14											
*29:07, 29:49	A29, -				5	6								14											
*29:13, 29:39	-													14		17									
*29:14	-													14	16										
*29:28	-							8						14											
*29:48	-													14			18								
*29:67	-			4										14											
*30:01:01:30:04:02, 30:06-30:07, 30:09-30:20, 30:22-30:54, 30:56-30:88, 30:90-30:98, 30:102, B*07:260 ⁴	A30, Null, -														15										
*30:08	A30			4											15										
*30:55	-												13	15											
*30:89	-													15		17									
*30:99	-			3										15											
*31:01:02:01-31:02, 31:05-31:07, 31:09-31:20, 31:22-31:23, 31:25-31:28, 31:30-31:34, 31:36-31:61, 31:63-31:84, 31:86-31:108, 31:110-31:111	A31, Null, -														16										
*31:03	A31											12	14	16											
*31:04	A31											12	14	16								24			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*31:08, 33:53 ⁵	A31/A24, -					5										16									
*31:21	-															16	17								
*31:24	-												13			16									
*31:29	-					5	6									16									
*31:35	-														15	16	17								
*31:62	-															16					22				
*31:85	-					6										16									
*31:109	-															16		18							
*32:01:01:32:02, 32:06-32:14, 32:16-32:37, 32:39-32:41, 32:43-32:59, 32:61, 32:63, 32:65-32:70, 32:73-32:77, 32:81-32:86	A32, Null, -								9								17								
*32:03, 32:38, 32:60, 32:71, 32:78, 32:80	A32, -																17								
*32:04	A32/A3			3					9								17								
*32:05, 32:79	A32, -					5										16	17								
*32:15	A32							8	9								17	18							
*32:42	-								9				14				17								
*32:62	-								9								17		19						
*32:64	-			4					9								17								
*32:72	-					6			9								17								
*33:01:01:33:01:08, 33:03:01-33:12, 33:14-33:18:02, 33:20-33:21, 33:23, 33:25-33:37, 33:39-33:47, 33:49-33:50, 33:52, 33:54-33:60, 33:62-33:114	A33, Null, -																	18							
*33:13, 33:48	-									10							18								
*33:19	-				5		7										18								
*33:22	-																18			21					
*33:24	-																18					23			
*33:51	-							8									18								
*33:61	-												13				18								
*34:01:01-34:01:02, 34:05-34:06, 34:11-34:12, 34:14	A34, A10, -							8			11	12													
*34:02:01-34:02:03, 34:03, 34:07, 34:10N, 34:13, 34:15	A34, Null, -			3				8			11	12													
*34:02:04	-			3				8			11	12						20		22					
*34:04	A34			3				8			11	12		14											
*34:08	A34			3				8			11	12	13												
*34:09	-			3				8				12												24	
*36:01, 36:03, 36:05	A36, -	1																				22			
*36:02	A36	1		3											14							22	23		
*36:04	A36	1			4																	22			

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

*68:71	-							8				13						20							
*68:84	-							8		10								20							
*68:130	-							8				12						20							
*69:01, 69:03	A69, -							8												21					
*69:02	-							8			11									21					
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*74:01:01-74:01:02, 74:02:01:01-74:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:27	A74, Null, -																		19						
*74:01:03	-											12							19					24	
*74:04, 74:21	A74, -																		18	19					
*74:07	A74																	17	19						
*74:13	-													13		16			19						
*74:19	-				4														19						
*74:23	-			3															19						
*80:01:01:01-80:01:01:02, 80:03	A80, -						6																23		
*80:02	-						6					11													
B*07:81, B*08:52, B*18:67, B*38:41, B*51:185, B*53:05, B*53:16, B*53:33	B53, -								9																
B*18:27	-						6																		
C*07:404	-											11													
C*07:481	-																						?		
C*12:131	-				4																				
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24