

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01:01-01:02:17, 01:02:19-01:02:42, 01:05-01:08, 01:11-01:13, 01:15-01:01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79:01-01:89N, 01:91-01:100, 01:103-01:112, 01:114-01:115, 01:117N, 01:119-01:121Q, 01:123-01:130, 01:132-01:135, 01:137N-01:144	73																			92				
*01:02:18	73											84								92				
*01:03, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78	73																			92	93			
*01:04	73											84		86						92				
*01:09	73													86						92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118	73											83								92				
*01:21	73											84	85											
*01:43	73	74									82									92				
*01:60	73																		91	92				
*01:67, 01:122	73			76												88				92				
*01:90, 01:101-01:102, 01:113	73							80									89			92				
*01:116	73																89			92				
*01:131	73			76																92	93			
*01:136, 01:145N	73							80												92				
*02:02:01-02:02:03, 02:02:05-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:31-02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:114, 02:116-02:117, 02:119-02:130		74	75			78					83				87									
*02:02:13, 02:02:29		74	75								83				87						94			
*02:02:26, 02:21, 02:26:02-02:26:03, 02:39, 02:45		74	75								83				87									
*02:02:30		74				78					83				87									
*02:04		74	w			78					83				87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87		89							
*02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01		74	w			78					83													
*02:14:02		74				78					83													
*02:22, 02:62, 02:118		74	75			78					83				87						93			
*02:26:01		74									83				87									
*02:27:01-02:27:02, 02:65		74	75			78					83													
*02:40:01-02:40:02		74	75			78					83				87				91					
*02:47		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87			90						
*02:58			75			78					83				87									
*02:68		74	75			78					83				87		89							
*02:87		74									82													
*02:94		74	75			78	79				83				87									
*02:101		74	75			w					83				87									
*02:107		74									83													
*02:115		74	75			78					83		85					90						
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

[illegible]

*03:04:25, 03:145, 03:172, 03:342			75	76																	93			
*03:07:01-03:07:02, 03:45, 03:163, 03:297			75	76						83				87							93	94		
*03:10, 03:29, 03:209			75	76										87							93	94		
*03:13:01-03:13:02, 03:97, 03:132, 03:167, 03:177, 03:271				76	77																	94		
*03:15			75	76						83				87								94		
*03:58					77									87				92		93	94			
*03:67, 03:295				76	w																93	w		
*03:81, 03:175				76	77											89					93	94		
*03:86					77									87					92			94		
*03:94, 03:99:02			75											87					92			94		
*03:99:01														87					92					
*03:130, 03:140			75	76						83											93	94		
*03:151, 03:165, 03:251, 03:267					77																93	94		
*03:155			75	76												88					93	94		
*03:199, 03:245			75	76												89						94		
*03:231				76		78				83					88						93			
*03:243, 03:308				76	77					83											93	94		
*03:248			75	76						83														
*03:262				76	77																?	94		
*03:268				76	77					83				87							93	94		
*03:278			w	76																	93	w		
*03:280			75	76						82											93	94		
*03:302			75	76										87								94		
*03:317			75	76												89					93	94		
*03:324				76	77																93			
*04:01:01-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:60, 04:01:62, 04:01:64, 04:01:66-04:01:73, 04:01:75-04:01:83, 04:04:01-01-04:05, 04:07-04:09N, 04:12-04:13, 04:15:01-04:15:03, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:34-04:35, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:102, 04:104-04:106, 04:109-04:111, 04:113, 04:115N, 04:117-04:128, 04:131-04:139, 04:141-04:146, 04:148-04:152, 04:154-04:159, 04:161-04:168, 04:170N, 04:173N-04:177, 04:179-04:189, 04:191N, 04:193, 04:195, 04:197-04:209, 04:211, 04:216-04:219, 04:221-04:226, 04:228-04:229, 04:232-04:235, 04:237-04:248, 04:250-04:255N, 04:257-04:262, 04:264, 04:266-04:276							78			83		85					90				93			
*04:01:16, 04:23, 04:54, 04:72						78				83		85					90							
*04:01:23, 04:33, 04:116, 04:130, 04:231						78				83		85									93			
*04:01:28, 04:01:61, 04:01:63, 04:01:65, 04:140, 04:196, 04:213, 04:236N										83		85					90					93		
*04:01:74						78				83		85					w					93		
*04:03:01-04:03:03, 04:06, 04:80, 04:147, 04:171, 04:256			75			78				83		85		87			90					93		
*04:10, 04:153, 04:192, 04:194, 04:210, 04:214-04:215N, 04:227, 04:249, 04:278						78				83							90					93		
*04:11						78				83				88								93		
*04:14, 04:68						78				83		85					90					93	95	
*04:16			75							83		85		87			90					93		
*04:29, 04:114						78				83		85		88			90					93		
*04:32, 04:77, 04:265			74			78				83		85					90					93	94	
*04:36, 04:55						78				83				88			90					93		
*04:37						78				83		85		87			90		92			93		
*04:42:01-04:42:02, 04:220										83		85	86				90					93	95	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*04:58						78					83		85					90	91		93			
*04:71	73					78					83		85					90			93			

Negative Control

*04:103											83	85	87		90			93						
*04:107			75			78					83	85	87					93						
*04:108						78		w			83	85			89	90								
*04:112						78					83	85			89	90		93						
*04:129						78	79				83	85				90		93						
*04:160			75			78					83	85	87			90	91	93						
*04:169						78					83				89			93						
*04:172						78					83	85		88				93						
*04:178						78		80			83	85			89	90		93						
*04:190						78					83	85	87			90		93						
*04:212				76		78					83	85				90		93						
*04:230, 04:263						78					83	85				90		92	93					
*04:279N						78					83	w				90			93					
*05:01:01:01-05:01:35, 05:03, 05:05:01-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:58:02, 05:59-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106:01-05:106:02, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:142, 05:144-05:145, 05:148-05:150, 05:152-05:154N							79				83		87					93						
*05:04:01-05:04:02, 05:22:01-05:22:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147							79				83		87											
*05:16, 05:85							79				83		87				92	93						
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115							79				83		87		90			93						
*05:20, 05:43, 05:113N							79				83							93						
*05:23, 05:62, 05:146, 05:151							79				83		87			91		93						
*05:26						78	79				83		87					93						
*05:36							79				83		87		89			93						
*05:42							79				83	85	87		90			93						
*05:46							79				83	85	87					93						
*05:58:01, 05:58:03			75				79				83		87					93						
*05:78:01-05:78:02						78	79				83							93						
*05:105		74					79				83		87					93						
*05:107							79				83		87				92							
*05:112							79				83	85	87					93	95					
*05:129							79						87					93						
*05:132							79	80			83		87		89									
*05:134, 05:143							79				83		87			90	91	93						
*06:02:01:01-06:02:01:08, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44-06:02:50, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:179, 06:181-06:196, 06:198-06:202								80			83		86						95					
*06:02:12, 06:11, 06:52, 06:122, 06:124, 06:147								80					86							95				
*06:02:16, 06:02:43, 06:04:01-06:04:02, 06:40, 06:79N, 06:81, 06:101, 06:123, 06:127:01:01-06:127:02, 06:136, 06:153, 06:197, 06:203								80			83									95				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*06:02:38								80			83	84		86									95	
*06:02:39, 06:32, 06:94-06:95, 06:99, 06:180											83			86									95	

Negative Control

*06:03:01, 06:132:01			75					80			83			86	87					94				
*06:03:02, 06:76:01								80			83			86	87									
*06:05							79	80			83						90		92		95			
*06:06, 06:204								80			83								92		95			
*06:08			74					80			83			86							95			
*06:09								80			83			86						93	95			
*06:14, 06:35, 06:72, 06:143-06:144								80			83									93	95			
*06:17											83			86					92		95			
*06:31								80			83			86			90		92		95			
*06:67							79	80			83			86							95			
*06:76:02											83	84	85	86	87		90							
*06:77								80			83			86										
*06:89								80			83			86			89				95			
*06:118								80			83			86			90	91			95			
*06:132:02			75								83			86	87						94			
*07:01:01:01-07:01:40, 07:01:42-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:06, 07:08, 07:11-07:25, 07:27:01-07:27:02, 07:29:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:215, 07:217-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:266, 07:268-07:293, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:401, 07:404, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518-07:530, 07:532-07:535, 07:537-07:562, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599									81											95				
*07:01:41, 07:02:41, 07:223, 07:415-07:416											81										w			
*07:02:08, 07:02:29, 07:26:01-07:26:02, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405, 07:493, 07:517											81													
*07:02:09							78				81			85			90				95			
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598											81		83						92		95			
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568											81									93	95			
*07:31:01-07:31:02, 07:514											81			86			90				95			
*07:39-07:40, 07:335											81								92					
*07:53, 07:216											81		83								95			
*07:64, 07:402, 07:449											81					88								
*07:96:01-07:96:02, 07:314:02, 07:578			75								81										94			
*07:101, 07:148, 07:161			74								81	82									95			
*07:123, 07:173, 07:294											81						89				95			
*07:125, 07:356, 07:531							78				81		83	85	86						95			
*07:154											81						90				95			
*07:177											81			86			90		92					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*07:210, 07:238									81		83								92					
*07:242				76					81														95	
*07:316, 07:338, 07:579	73								81														95	

*07:563										81										92			95	
*07:583		74						78		81	82			85						90				
*08:01:01-08:01:17, 08:01:19-08:01:21, 08:03:01-08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153											82								90			93		
*08:01:18, 08:02:01-08:02:15, 08:04:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62-08:63, 08:66-08:71, 08:73-08:77, 08:90, 08:92-08:94, 08:100, 08:103, 08:107, 08:110-08:112, 08:114-08:116, 08:120, 08:123, 08:125-08:126, 08:132, 08:134, 08:138-08:140, 08:146, 08:150-08:151											82											93		
*08:05, 08:25											82			85								93		
*08:07, 08:104, 08:142											82								91			93		
*08:09, 08:11, 08:83											82								90					
*08:10								79			82	83			87				90			93		
*08:12											82									92		93		
*08:21, 08:137											82			85					90			93		
*08:31		74									82											93		
*08:47											82								90	91		93		
*08:48, 08:108, 08:152											82													
*08:95											82					88			90			93		
*08:113								80			82							89						
*12:02:01-12:02:10, 12:02:12-12:02:13, 12:10:01-12:10:02, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N, 12:151, 12:153, 12:155Q, 12:161-12:162, 12:164, 12:166, 12:168, 12:175, 12:177, 12:179, 12:183, 12:193, 12:196, 12:200, 12:204, 12:207, 12:212, 12:214												84	85											
*12:02:11								80				84												
*12:02:14-12:02:15, 12:18:02, B*67:02 ⁴													85											
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

[illegible]

[illegible]

*15:111					74	75					w					83					87					89					93						
*15:144							77					80															89	90					93				
*16:01:01:01-16:01:15, 16:01:18-16:01:22, 16:06-16:08, 16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98, 16:100, 16:105-16:106, 16:111-16:112																																					
*16:01:16																																					
*16:01:17, 16:15:02																																					
*16:02:01-16:02:12, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108																																					
*16:02:13																																					
*16:04:01:01-16:04:01:02, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82																																					
*16:10, 16:67																																					
*16:20																																					
*16:29																																					
*16:34																																					
*16:35																																					
*16:40, 16:110																																					
*16:43																																					
*16:45																																					
*16:48																																					
*16:50																																					
*16:53																																					
*16:64, 16:87																																					
*16:68																																					
*16:70																																					
*16:88																																					
*16:91																																					
*16:103																																					
*16:109																																					
*16:113																																					
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:38																																					
*17:01:04																																					
*17:05																																					
*17:07, 17:17, 17:23, 17:33																																					
*17:20																																					
*17:22																																					
*18:01-18:02, 18:04-18:08, 18:10																																					
*18:03																																					
*18:09																																					
A*24:64, A*24:96, A*24:106, B*08:16, B*08:123, B*08:145																																					
A*24:174																																					
B*07:239, B*14:03																																					
B*13:13:02, B*40:01:49																																					
B*15:277, B*35:222, B*42:22, B*54:38																																					
B*39:76, B*54:18																																					
B*44:02:44, B*56:01:09, B*58:74																																					
B*58:02:01																																					
B*73:01																																					
Well No.																																					

</

Negative Control