

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01:01-01:02:17, 01:02:19-01:02:55, 01:05-01:08, 01:11-01:13, 01:15-01:01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01:01:89N, 01:91-01:100, 01:103-01:112, 01:114-01:01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:146-01:01:155, 01:157-01:159, 01:161-01:167, 01:170-01:175	73																			92				
*01:02:18	73											84								92				
*01:03, 01:24, 01:34, 01:49-01:01:50, 01:55, 01:78, 01:176	73																			92	93			
*01:04	73											84		86						92				
*01:09	73													86						92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118	73										83									92				
*01:21	73											84	85											
*01:43	73	74								82										92				
*01:60	73																		91	92				
*01:67, 01:122	73		76													88				92				
*01:90, 01:101-01:102, 01:113	73							80									89			92				
*01:116	73																89			92				
*01:125, 14:58 ³	73															88				92				
*01:131	73		76																	92	93			
*01:136, 01:145N	73							80												92				
*01:156, B*07:239, B*14:03 ³																				92				
*01:160	73													86				90		92				
*01:168-01:169:01															87					92				
*01:169:02, 03:94 ³			75												87					92		94		
*01:177																88				92		94		
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:31-02:02:37, 02:02:39, 02:02:41-02:02:43, 02:02:45-02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:169N		74	75			78					83				87									
*02:02:13, 02:02:29		74	75								83				87							94		
*02:02:26, 02:21, 02:26-02:02:26:03, 02:39, 02:45		74	75								83				87									
*02:02:30, 02:02:40		74				78					83				87									
*02:02:38		74	75			78					83				87							94		
*02:02:44		74	75			78					83	84			87									
*02:04		74	w			78					83				87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87		89							
*02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01		74	w			78					83													
*02:14:02, 02:164		74				78					83													
*02:22, 02:62, 02:118		74	75			78					83				87						93			
*02:26:01, 02:134, 02:152, 02:167		74									83				87									
*02:27:01-02:27:02, 02:65, 02:131		74	75			78					83													
*02:40:01-02:40:02		74	75			78					83				87				91					
*02:47		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87			90						
*02:58			75			78					83				87									
*02:68		74	75			78					83				87		89							

Negative Control

*02:87				74								82															
*02:94				74	75				78	79			83				87										
*02:101				74	75				w				83				87										
*02:107				74									83														
*02:115				74	75				78				83		85					90							
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96			
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:02:14-03:02:19, 03:02:21-03:02:22, 03:05, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298:03:301, 03:315, 03:326, 03:329-03:330, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464				75	76																	94					
*03:02:10, 03:264, B*44:274 ³					76																						
*03:02:13, 03:02:20, 03:16, 03:110, 03:332				75																			94				
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:54, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460					76	77																93	94				
*03:03:21, 03:03:51, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195, 03:344:01:01-03:344:01:02, 03:454, 03:459					76																		93				
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451						77																	93	94			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96			

[illegible]

[illegible]

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*05:01:01:01-05:01:45, 05:03, 05:05:01-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:58:02, 05:58:04-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106:01-05:106:02, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:142, 05:144-05:145, 05:148-05:150, 05:152-05:155, 05:157-05:162, 05:164-05:190, 05:192-05:199, 05:201, 05:203-05:205, 05:208N-05:214							79				83				87						93				
*05:04:01-05:04:02, 05:22:01-05:22:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147, 05:191, 05:207							79				83				87										
*05:16, 05:85							79				83				87					92	93				
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115							79				83				87			90			93				
*05:20, 05:43, 05:113N, 05:156, 05:202Q							79				83										93				
*05:23, 05:62, 05:146, 05:151, 05:163							79				83				87				91		93				
*05:26						78	79				83				87						93				
*05:36							79				83				87		89				93				
*05:42							79				83		85		87			90			93				
*05:46							79				83		85		87						93				
*05:58:01, 05:58:03			75				79				83				87						93				
*05:78:01-05:78:02						78	79				83										93				
*05:105		74					79				83				87						93				
*05:107							79				83				87					92					
*05:112							79				83		85		87						93			95	
*05:129							79								87						93				
*05:132, 05:215							79	80			83				87		89								
*05:134, 05:143							79				83				87			90	91		93				
*05:200	73	74					79				83										93				
*05:206, 17:05 ³		74					79				83				87				91		93				
*06:02:01:01-06:02:01:17, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:179, 06:181-06:196, 06:198-06:202, 06:205-06:216, 06:218-06:223, 06:225-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N								80			83			86									95		
*06:02:12, 06:02:65, 06:11, 06:52, 06:122, 06:124, 06:147, 06:217, 06:248, 06:252								80						86										95	
*06:02:16, 06:02:43, 06:04:01-06:04:02, 06:40, 06:79N, 06:81, 06:101, 06:123, 06:127:01-01-06:127:02, 06:136, 06:153, 06:197, 06:203								80			83												95		
*06:02:38								80			83	84		86									95		
*06:02:39, 06:32, 06:94-06:95, 06:99, 06:180, 06:236, 06:244											83			86									95		
*06:03:01, 06:132:01			75					80			83			86	87							94			
*06:03:02, 06:76:01								80			83			86	87										
*06:05							79	80			83							90		92			95		
*06:06, 06:204								80			83									92			95		
*06:08		74						80			83			86									95		
*06:09:01-06:09:02								80			83			86							93		95		
*06:14, 06:35, 06:72, 06:143-06:144								80			83										93		95		
*06:17											83			86						92			95		
*06:31								80			83			86				90		92			95		
*06:67							79	80			83			86									95		

Negative Control

Negative Control

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]