

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01-01:02:17, 01:02:19-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72-01:02:80, 01:05-01:08, 01:11-01:13, 01:15-01:01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:47, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01-01:89N, 01:91-01:100, 01:103-01:112, 01:114-01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:214, 01:216-01:217:01:02. 01:219	73																			92				
*01:02:18	73											84								92				
*01:02:49, 01:136, 01:145:01N, 01:146:01-01:146:02	73							80												92				
*01:02:68, 01:02:71, 01:168-01:169:01	73														87					92				
*01:03:01-01:03:02, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:176	73																			92	93			
*01:04	73											84			86					92				
*01:09:01-01:09:02	73														86					92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118, 01:186	73										83									92				
*01:21	73											84	85											
*01:43	73	74								82										92				
*01:48	73		75																	92				
*01:60	73																			91	92			
*01:67, 01:122	73			76												88				92				
*01:90, 01:101-01:102, 01:113	73							80									89			92				
*01:116	73																89			92				
*01:125, 01:218, 14:58 ³	73															88				92				
*01:131	73			76																92	93			
*01:160	73													86				90		92				
*01:169:02, 03:94 ³	73		75												87					92			94	
*01:177	73															88				92			94	
*01:215	73												85							92				95
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:30-02:02:37, 02:02:39-02:02:43, 02:02:45-02:02:52, 02:02:54-02:04, 02:07-02:11, 02:15-02:16:02:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:110, 02:112-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186, 02:188, 02:190, 02:192N-02:196, 02:198-02:207		74	75				78				83				87									
*02:02:13, 02:02:29		74	75								83				87							94		
*02:02:26, 02:21, 02:26:01-02:26:03, 02:39, 02:45, 02:189, 02:191		74	75								83				87									
*02:02:38, 02:02:53		74	75				78				83				87							94		
*02:02:44		74	75				78				83	84			87									
*02:05:01-02:05:03, 02:17		74	75				78				83				86	87								
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75				78	80			83				87		89							
*02:12, 02:55:01-02:55:02		74	75				78				83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75				78				83				87			90						
*02:14:01-02:14:02, 02:27:01-02:27:02, 02:65, 02:131, 02:164, 02:187		74	75				78				83													
*02:22, 02:62, 02:118		74	75				78				83				87						93			
*02:40:01-02:40:02		74	75				78				83				87				91					
*02:47, 02:111		74	75				78	80			83				87									
*02:49		74	75				78				83		85		87			90						
*02:58			75				78				83				87									
*02:68		74	75				78				83				87		89							
*02:87		74								82														
*02:94		74	75				78	79			83				87									
*02:101		74	75				w				83				87									
*02:107		74									83													
*02:115		74	75				78				83		85					90						
*02:134, 02:152, 02:167		74									83				87									
*02:180		74	75				78				83				87			90	91					
*02:197			75				78				83				87						93			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*03:02:01-03:02:09, 03:02:11-03:02:26, 03:05:01:01-03:05:01:02, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505, 03:532, 03:544, 03:547, 03:549, 03:555, 03:561, 03:575			75	76																	94				
*03:02:10, 03:264			75	76																					
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:62, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119-01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495, 03:503, 03:506-03:507, 03:509N, 03:514-03:517, 03:519-03:520, 03:524, 03:527-03:529, 03:533, 03:537, 03:541-03:542, 03:545, 03:550-03:551, 03:556, 03:559-03:560N, 03:563, 03:567-03:568, 03:573-03:574			76	77																93	94				
*03:03:21, 03:03:51, 03:18:01, 03:39, 03:64:02, 03:144, 03:182, 03:195, 03:344:01-01-03:344:01:02, 03:504				76																	93				
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451, 03:534-03:535, 03:540					77																93	94			
*03:04:01-01-03:04:06, 03:04:08-03:04:14, 03:04:16-03:04:24, 03:04:27-03:04:94, 03:06:01-03:06:02, 03:08-03:09, 03:14, 03:19, 03:23N-03:24, 03:26:01-03:26:02, 03:28, 03:32:01-01-03:32:01:02, 03:34, 03:37:01-03:38:02, 03:40:01-03:40:04, 03:42:01-03:42:02, 03:44, 03:46-03:48, 03:51, 03:54, 03:57:01-03:57:02, 03:63-03:64:01, 03:65, 03:70, 03:72-03:74, 03:77-03:78, 03:80:02, 03:82, 03:87:01-03:87:02, 03:90-03:93, 03:98, 03:100-03:101, 03:104-03:107, 03:109, 03:111, 03:114-03:115:02, 03:117, 03:123, 03:125, 03:128-03:129, 03:131, 03:134, 03:136, 03:138, 03:143, 03:147-03:149, 03:153, 03:157, 03:159, 03:162, 03:164, 03:173-03:174, 03:179-03:181, 03:183-03:184:02, 03:186:01-03:186:02, 03:191, 03:193, 03:208N, 03:210-03:213, 03:215, 03:218-03:219, 03:232-03:236, 03:239-03:240, 03:244Q, 03:246-03:247, 03:250, 03:252, 03:255-03:257, 03:259, 03:261, 03:263:01-03:263:03, 03:265N-03:266, 03:269-03:270, 03:277N, 03:279, 03:281-03:283, 03:287:01-03:287:02, 03:294, 03:303, 03:305-03:306, 03:309-03:311, 03:313-03:314, 03:318N, 03:322-03:323N, 03:328, 03:331, 03:333-03:334, 03:337, 03:340, 03:343, 03:347-03:348, 03:353-03:355, 03:358-03:359, 03:361-03:362, 03:366N, 03:368-03:369, 03:376, 03:379, 03:381-03:382, 03:385, 03:387, 03:390, 03:392N-03:394, 03:396:01N-03:397, 03:399-03:401, 03:404:01-01-03:405, 03:408-03:410, 03:415, 03:417, 03:419-03:420, 03:423-03:424N, 03:426, 03:434, 03:437, 03:439, 03:441-03:443, 03:445N-03:446N, 03:448Q, 03:453, 03:455-03:456, 03:458, 03:462N-03:463N, 03:465-03:467, 03:469-03:470N, 03:478, 03:480, 03:483-03:488, 03:491:01-03:492, 03:496, 03:498-03:502, 03:508N, 03:510N-03:511N, 03:513, 03:518, 03:521-03:523, 03:525-03:526, 03:531N, 03:538-03:539, 03:543, 03:546, 03:548, 03:552-03:554, 03:558, 03:562, 03:564-03:565, 03:569-03:572, 03:576				75	76															93	94				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
Negative Control																									

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*03:04:07, 03:41:01-03:41:02, 03:80:01, 03:113:01-03:113:02, 03:260, 03:286, 03:293, 03:384, 03:435, 03:512			75																		93	94		
*03:04:15, 03:04:25-03:04:26, 03:118, 03:137, 03:145, 03:156, 03:166, 03:170, 03:172, 03:342, 03:454, 03:459			75	76																	93			
*03:07:01-01-03:07:02, 03:45, 03:163, 03:297			75	76						83					87						93	94		
*03:10, 03:29, 03:209			75	76											87						93	94		
*03:13:01-01-03:13:02, 03:97, 03:132, 03:167, 03:177, 03:271, 03:474				76	77																	94		
*03:15			75	76						83					87							94		
*03:16, 03:536, B*46:77 ³			75																			94		
*03:58	73				77										87					92	93	94		
*03:67, 03:295				76	w																93	w		
*03:81, 03:175				76	77												89				93	94		
*03:86	73				77										87					92		94		
*03:130, 03:140:01-01-03:140:01:02			75	76						83											93	94		
*03:155			75	76												88					93	94		
*03:199, 03:245			75	76													89					94		
*03:231				76		78				83						88					93			
*03:243, 03:308				76	77					83											93	94		
*03:248			75	76						83														
*03:262				76	77																?	94		
*03:268, 03:461				76	77					83					87						93	94		
*03:278, 03:566			75	76																	93	w		
*03:280			75	76						82											93	94		
*03:302			75	76											87							94		
*03:317, 03:388			75	76													89				93	94		
*03:324, 03:479				76	77																93			
*03:351				76	77											88					93	94		
*03:416	73			76	77																93	94		
*03:450				76	77										87						93	94		
*03:530				76						82											93			
*03:557				76												88					93	94		
*04:01:01-01-04:01:01:29, 04:01:01:31-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:60, 04:01:62, 04:01:64, 04:01:66-04:01:73, 04:01:75-04:01:141, 04:04:01:01-04:05, 04:07:01, 04:08-04:09N, 04:12-04:13:01:02, 04:15:01-04:15:03, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:34-04:35:02, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:99, 04:101-04:102, 04:104-04:106, 04:109-04:111, 04:113, 04:115N, 04:117-04:128, 04:131-04:139, 04:141-04:146, 04:148-04:152, 04:154-04:159, 04:161-04:166:01, 04:166:03-04:168, 04:170N, 04:173N-04:177, 04:179-04:189, 04:191N, 04:193, 04:195, 04:197-04:201:01, 04:202-04:209, 04:211, 04:216-04:219, 04:221-04:226, 04:228-04:229, 04:232-04:235, 04:237-04:248, 04:250-04:255N, 04:257-04:262, 04:264, 04:266-04:277, 04:280-04:282, 04:284-04:285, 04:287-04:293, 04:295-04:298, 04:300N-04:303, 04:305N-04:310, 04:312-04:334, 04:336, 04:338Q-04:340, 04:342-04:346, 04:348-04:350N, 04:352-04:356, 04:358-04:359, 04:361-04:362N, 04:364N-04:366, 04:369N-04:380, 04:382Q, 04:384-04:392, 04:394-04:396N, 04:398-04:399, 04:401, 04:403-04:405, 04:407-04:413, 04:415-04:421, 04:423-04:431, 04:433-04:435, 04:437-04:438, 04:440-04:441, 04:443-04:446, 04:448-04:459						78				83	85					90				93				
*04:01:16, 04:54:01-04:54:02, 04:72, 04:341, 04:367, 04:439						78				83			85				90							
*04:01:23, 04:33, 04:116, 04:130, 04:166:02, 04:231, 04:283						78				83			85								93			
*04:01:28, 04:01:61, 04:01:63, 04:01:65, 04:140, 04:196, 04:213, 04:236N, 04:304:01-04:304:02:02, 04:311										83			85				90				93			
*04:01:74						78				83			85				w				93			
*04:03:01-01-04:03:08, 04:06:01-04:06:03, 04:80, 04:147, 04:171, 04:190, 04:256, 04:286, 04:294, 04:337, 04:363, 04:381, 04:393, 04:400, 04:402			75			78				83			85		87		90				93			
*04:03:09			75							83					87						93	94		
*04:07:02										83											93	94		
*04:10, 04:153, 04:192, 04:194, 04:210, 04:214-04:215N, 04:227, 04:249, 04:278, 04:347						78				83							90				93			
*04:11						78				83						88					93			
*04:14, 04:68						78				83			85				90				93		95	
*04:16, 04:351:01-04:351:02			75							83			85		87		90				93			
*04:23						78		80		83			85				90							
*04:29, 04:114						78				83			85			88	90				93			
*04:32, 04:77, 04:265		74				78				83			85				90				93	94		
*04:36, 04:55						78				83						88	90				93			
*04:37	73					78				83			85		87		90		92		93			
*04:42:01-04:42:02, 04:220										83			85	86			90				93		95	
*04:58, 04:368						78				83			85				90	91			93			
*04:71, 04:436N	73					78				83			85				90				93			
*04:100						78		80		83			85				90				93			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*04:103:01-04:103:02											83		85		87			90			93			
*04:107			75			78					83		85		87						93			
*04:108						78		80			83		85				89	90						
*04:112						78					83		85				89	90			93			
*04:129						78	79				83		85					90			93			
*04:160:01-04:160:02			75			78					83		85		87			90	91		93			
*04:169						78					83						89				93			
*04:172						78					83		85			88					93			
*04:178, 04:422						78		80			83		85				89	90			93			
*04:201:02, 04:212					76	78					83		85					90			93			
*04:230, 04:263, 04:397, 04:432, 04:442						78					83		85					90		92	93			
*04:279N						78					83		w					90			93			
*04:299			75			78					83				87			90			93			
*04:335			75			78		80			83		85		87			90						
*04:357		74	75			78					83		85		87			90			93			
*04:360						78					83		85	86				90						
*04:383			75			78					83		85					90			93			
*04:406		74				78					83		85					90			93			
*04:414						78					83										93			
*04:447		74				78					83		85					90	91		93	94		
*05:01:01:01-05:01:50, 05:01:52-05:01:55, 05:01:57, 05:01:59-05:01:68, 05:03, 05:05:01-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:59-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106:01-05:106:02, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:142, 05:144-05:145, 05:148-05:150, 05:152-05:155, 05:157-05:162, 05:164-05:190, 05:192-05:199, 05:201, 05:203-05:205, 05:208N-05:214, 05:216:01:01-05:237, 05:239N-05:240, 05:242-05:244N, 05:246-05:251N, 05:253N-05:263N							79				83				87						93			
*05:01:51, 05:01:58, 05:04:01-05:04:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147, 05:191, 05:207							79				83				87									
*05:01:56							79	80			83				87						93			
*05:16							79				83				87					92	93			
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115							79				83				87			90			93			
*05:20, 05:43, 05:113N, 05:156, 05:202Q							79				83										93			
*05:22:01-05:22:02, 05:252							79	80			83				87									
*05:23, 05:62, 05:146, 05:151, 05:163							79				83				87				91		93			
*05:26						78	79				83				87						93			
*05:36							79				83				87		89				93			
*05:42							79				83		85	87			90				93			
*05:46:01:01-05:46:01:02							79				83		85	87							93			
*05:58:01-05:58:04			75				79				83			87							93			
*05:78:01-05:78:02						78	79				83										93			
*05:85	73						79				83				87					92	93			
*05:105		74					79				83				87						93			
*05:107	73						79				83				87					92				
*05:112							79				83		85	87							93		95	
*05:129							79							87							93			
*05:132, 05:215							79	80			83			87			89							
*05:134, 05:143							79				83			87				90	91		93			
*05:200	73	74					79				83										93			
*05:206, 17:05 ³		74					79				83			87					91		93			
*05:238							79				83			87			89							
*05:241	73						79	80			83				87					92				
*06:02:01:01-06:02:01:62, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73-06:02:96, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:69, 06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:179, 06:181-06:196, 06:198-06:202, 06:205-06:216, 06:218-06:223, 06:225-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N, 06:265-06:270, 06:272-06:279, 06:281N-06:303, 06:306-06:316N, 06:318-06:328, 06:330-06:334								80			83		86									95		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*06:02:12, 06:02:65, 06:11, 06:52, 06:122, 06:124, 06:147, 06:217, 06:248, 06:252, 06:280								80						86										95	
*06:02:16, 06:02:43, 06:04:01-06:04:02:02, 06:40, 06:81, 06:101, 06:123, 06:127:01-01-06:127:02, 06:136, 06:153:01-06:153:02, 06:197, 06:203, 06:264, 06:304, 06:317								80			83													95	
*06:02:38								80			83	84		86										95	
*06:02:39, 06:32, 06:94-06:95, 06:99, 06:180, 06:236, 06:244, 06:329:01-01-06:329:01:02											83			86										95	
*06:02:72											83		85	86				90						95	
*06:03:01, 06:132:01			75					80			83			86	87								94		
*06:03:02, 06:76:01								80			83			86	87										
*06:05							79	80			83							90		92				95	
*06:06:01-01-06:06:01:02, 06:79N, 06:271								80			83									92				95	
*06:08		74						80			83			86										95	
*06:09:01-06:09:02								80			83			86								93		95	
*06:14, 06:35, 06:72, 06:143-06:144								80			83											93		95	
*06:17											83			86						92				95	
*06:31								80			83			86				90		92				95	
*06:67							79	80			83			86										95	
*06:70:01-06:70:02, 06:305			75					80			83			86										95	
*06:76:02											83	84	85	86	87			90							
*06:77								80			83			86											
*06:89, 06:224								80			83			86			89							95	
*06:118								80			83			86				90	91					95	
*06:132:02			75								83			86	87							94			
*06:204	73							80			83									92				95	
*06:335											83		85	86										95	
*07:01:01-07:01:02:10, 07:01:04-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:01:94-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:06:04, 07:08, 07:11-07:19, 07:21-07:25, 07:27:01-07:27:02, 07:29:01-01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126, 07:128-07:132, 07:133:02-07:138, 07:140:01-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:215, 07:217-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:389, 07:391-07:401, 07:404, 07:405:02, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:530, 07:532-07:535, 07:537-07:561, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:713Q, 07:715:01-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844, 07:846:01-07:849N, 07:851-07:856N, 07:858-07:869, 07:871-07:890, 07:893, 07:895-07:931, 07:933-07:935, 07:937-07:939, 07:941-07:969, 07:971-07:974Q, 07:976-07:978									81														95		
*07:01:41, 07:02:41, 07:223, 07:415-07:416									81															w	
*07:01:74						78			81				85				89	90						95	
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:26:01-07:26:03, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405:01, 07:449, 07:493, 07:517, 07:714, 07:724, 07:821N, 07:845, 07:870, 07:891, 07:894, 07:936N, 07:970N									81																
*07:02:09						78			81				85					90						95	
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598									81	83										92				95	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568, 07:857, 07:892									81												93		95		
*07:20:01-07:20:02, 07:127:01-07:127:02, 07:263, 07:390, 07:562			75						81															95	
*07:31:01-07:31:02, 07:514:01-07:514:02									81				86				90							95	
*07:39-07:40, 07:335, 07:975									81										92						
*07:53, 07:216, 07:932									81		83													95	
*07:64, 07:402									81							88						94			
*07:96:01-07:96:02, 07:314:02, 07:578, 07:756			75						81													94			
*07:101, 07:148, 07:161			74						81	82														95	
*07:123:01-07:123:02, 07:173, 07:294, 07:626, 07:725, 07:775, 07:850									81								89							95	
*07:125:01, 07:356, 07:531						78			81	83			85	86										95	
*07:125:02						78			81	83			85	86				90						95	
*07:133:01, 07:242, 07:330:01, 07:748				76					81															95	
*07:154									81									90						95	
*07:177									81					86				90		92					
*07:210, 07:238									81	83										92					
*07:293, 07:519									81							88								95	
*07:316, 07:338, 07:579			73						81															95	
*07:563									81											92				95	
*07:583			74			78			81	82			85					90							
*07:723						78			81				85												
*07:940									81			84												95	
*08:01:01:01-08:01:17, 08:01:19-08:01:32, 08:03:01:01-08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22:01:01-08:22:01:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186-08:187, 08:189-08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203-08:204, 08:208N-08:210, 08:212, 08:217, 08:219-08:221, 08:228, 08:232										82								90			93				
*08:01:18, 08:02:01:01-08:02:34, 08:04:01:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19:02, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62:01-08:63, 08:66-08:71, 08:73, 08:75-08:77, 08:90, 08:92-08:94, 08:100, 08:103, 08:107, 08:110-08:112, 08:114-08:116, 08:120, 08:123, 08:125-08:126, 08:132, 08:134, 08:138-08:140, 08:146, 08:150-08:151, 08:156, 08:158-08:161N, 08:165-08:167, 08:169-08:172, 08:177, 08:179-08:182, 08:184, 08:191, 08:195, 08:198, 08:200-08:202, 08:206, 08:211, 08:213-08:216, 08:218, 08:222-08:227, 08:229-08:231, 08:233										82												93			
*08:05, 08:25										82			85										93		
*08:07, 08:104, 08:142										82									91			93			
*08:09, 08:11, 08:205										82								90							
*08:10							79			82	83				87			90					93		
*08:12:01:01-08:12:01:02										82										92			93		
*08:21, 08:137										82			85					90					93		
*08:31			74							82													93		
*08:47, 08:188										82								90	91			93			
*08:48, 08:152										82															
*08:74										82						88							93		
*08:83								80		82								90							
*08:95										82						88		90				93			
*08:105			75							82								90				93			
*08:108								80		82															
*08:113								80		82							89								
*08:163						78				82						88		90				93			
*08:183						78				82						88						93			
*08:185			74							82									91			93			
*08:207								80		82												93			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
Negative Control																									

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*12:02:01-12:02:10, 12:02:12-12:02:13, 12:02:16, 12:02:18, 12:02:20-12:02:27, 12:02:31, 12:02:33, 12:02:36, 12:02:38-12:02:39, 12:02:41-12:02:43, 12:03:47, 12:10:01-12:10:02, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N, 12:151, 12:153, 12:155Q, 12:161-12:162, 12:164, 12:166, 12:168, 12:175, 12:177, 12:179, 12:183, 12:193, 12:196, 12:200, 12:204, 12:207, 12:212, 12:214, 12:219N, 12:221, 12:226, 12:228, 12:231, 12:234, 12:236N, 12:240-12:241, 12:243, 12:247, 12:250-12:251, 12:255, 12:258, 12:261, 12:263, 12:268, 12:270N, 12:275, 12:279-12:281, 12:285, 12:294, 12:296, 12:303-12:304, 12:307-12:308, 12:315, 12:327N, 12:329N, 12:335-12:336, 12:338, 12:351Q												84	85												
*12:02:11								80				84													
*12:02:14-12:02:15, 12:02:19, 12:02:30, 12:02:32, 12:02:34, 12:02:37, 12:02:40, 12:18:02, 12:222, 12:224, B*67:02:01:01-67:02:01:02 ³													85												
*12:02:17, 12:239			75									84	85												
*12:02:28												84	w												
*12:02:29, 12:49, 12:123, B*18:01:42, B*57:01:24 ³												84													
*12:02:35								80				84	85												
*12:03:01:01-12:03:01:44, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:17-12:03:19, 12:03:21-12:03:22, 12:03:25, 12:03:27-12:03:29, 12:03:31-12:03:33, 12:03:35-12:03:46, 12:03:48-12:03:65, 12:03:67-12:03:77, 12:06-12:07, 12:11-12:13:01:02, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57:01-12:57:02, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82, 12:87-12:91, 12:93-12:95, 12:97, 12:99:01-12:102, 12:107-12:111, 12:115-12:116, 12:120-12:122, 12:125, 12:129, 12:131, 12:133, 12:138-12:139Q, 12:141, 12:143, 12:149-12:150, 12:152, 12:156-12:160, 12:163, 12:165, 12:167, 12:170-12:174, 12:180, 12:182, 12:184, 12:186-12:187, 12:189-12:192, 12:197, 12:199, 12:201-12:202, 12:205-12:206, 12:209-12:211, 12:213, 12:216, 12:218, 12:220, 12:223, 12:225, 12:227, 12:229-12:230, 12:232N, 12:235, 12:237-12:238, 12:242, 12:244-12:246, 12:248-12:249, 12:253-12:254, 12:256-12:257, 12:259-12:260:02, 12:262, 12:264-12:267, 12:271-12:274:01N, 12:276-12:278, 12:283-12:284, 12:286, 12:288-12:293, 12:295N, 12:297, 12:300, 12:302, 12:305-12:306, 12:310-12:314, 12:316-12:324N, 12:326, 12:328, 12:330N-12:332, 12:334, 12:337, 12:339-12:343N, 12:345N-12:348, 12:350												84	85	86											
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:03:30, 12:26, 12:59, 12:185													85	86											
*12:03:04, 12:98, 12:301												84	86												
*12:03:09								80				84	86												
*12:03:13, 12:03:24												84	w	86											
*12:03:23, 12:140:01-12:140:02			75									84	85	86											
*12:03:26, 12:03:66								80				84	85	86											
*12:03:34:01-12:03:34:02, B*15:510, B*35:475 ³													86												
*12:04:01											83		86	87											
*12:04:02:01-12:04:02:02, 12:60											83	84	85	86	87										
*12:05:01-12:05:02											83	84	86	87											
*12:08												84	85				89								
*12:09											83			87					92						
*12:14:01													85					90	91						
*12:14:02												84	85					90	91						
*12:15, 12:113								80				84	85	86			89								
*12:16:01, 12:217												84	85										95		
*12:16:02												84											95		
*12:21							79				83	84		87											
*12:24												84	85						92						
*12:28, 12:58, 12:325												84	85	86				90							
*12:31:01-12:31:02												84	85	86							93				
*12:33							79				83		86	87											
*12:41											83	84	85	87											
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*12:44, 12:344												84	85								93			
*12:54, 12:194											83	84	85	86										
*12:55												84	85	86		88								
*12:72												84	85		87									
*12:81, 12:282												84	85	86			89							
*12:92												84	85	86				91						
*12:119		74										84	85	86										
*12:127										82		84	85											
*12:132, 12:169, 12:252												84	85					90						
*12:135												84	85	86	87			90						
*12:144													85	86							93			
*12:146											83	84	85		87		90							
*12:147, 12:195:02												84	85	86									95	
*12:154, 12:178												84	85	86	87									
*12:176												84	85	86				90	91	92				
*12:181									81				85					90						
*12:188											83	84					89							
*12:195:01												84		86									95	
*12:195:03														86									95	
*12:198		74										84	85											
*12:203										82		84	85	86										
*12:208								80					85				89							
*12:215		73										84	85	86										
*12:233													85								93			
*12:269							78					84	85	86		88								
*12:287												84	85				89	90						
*12:298								80					85											
*12:309								80					85				89				93			
*12:333												84	85	86				90		92				
*12:349						78			81			84	85	86										
*14:02:01:01-14:02:07, 14:02:09-14:02:20, 14:02:22-14:02:40, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:21N, 14:23-14:24:02, 14:27, 14:29-14:34, 14:36-14:37, 14:39-14:40, 14:42-14:44, 14:46:01-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:60, 14:62, 14:65-14:69, 14:71-14:76, 14:78, 14:80-14:85, 14:88-14:89, 14:91, 14:93N-14:104, 14:106-14:107, 14:109, 14:111, 14:114-14:115, 14:117N, 14:119-14:120, 14:122-14:125, 14:127-14:138																	88			92				
*14:02:08, 14:02:21												84				88				92				
*14:03:01:01-14:03:04, 14:22, 14:35N, 14:38, 14:41, 14:61, 14:70, 14:79, 14:86, 14:112-14:113, 14:118, 14:126			75														88			92				
*14:03:05			75									84				88				92				
*14:04, 14:64											83					88				92				
*14:06, 14:08-14:09, 14:28:01-14:28:02, 14:63, 14:105Q, A*29:86, A*30:62, A*31:187, A*33:63, B*08:249, B*18:64, B*35:183, B*35:252, B*35:502, B*39:114, B*39:140, B*40:138, B*41:59 ³																88								
*14:10																88		90		92				
*14:12, 14:49, 14:108											83									92				
*14:15, 14:87													86		88			90						
*14:16													86		88									
*14:17															88			90	91	92				
*14:25		74													88					92		94		
*14:26								80							88					92				
*14:45		73													87	88				92				
*14:53			75													88								
*14:54			75													88					93			
*14:77											83					88								
*14:90			75								83		85			88							95	
*14:92																88	89							
*14:110, B*07:239, B*14:03 ³																				92				
*14:116			75													88		90		92				
*14:121									81							88				92				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*15:02:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:51, 15:05:01-15:05:07, 15:05:09-15:05:10, 15:05:13-15:05:17, 15:06:01-15:06:03, 15:08:01-15:09:02, 15:10:02-15:10:03, 15:13:01-15:13:02, 15:15, 15:17-15:19, 15:26-15:29, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44-01-15:54, 15:56-15:57, 15:59-15:61, 15:64, 15:67-15:70, 15:73-15:74, 15:78:01-15:84Q, 15:86-15:101, 15:103, 15:105Q-15:106, 15:108-15:110, 15:112-15:115N, 15:117, 15:119-15:124, 15:126-15:137, 15:139-15:143, 15:145N-15:148, 15:150-15:157, 15:159-15:161, 15:163-15:168, 15:170-15:177N, 15:179-15:180, 15:182-15:189N, 15:191-15:200, 15:202-15:207, 15:209-15:210, 15:212-15:217, 15:219-15:221, 15:223-15:224, 15:226-15:235Q, 15:237-15:243			75					80			83				87		89				93			
*15:02:10, 15:02:17			75					80			83				87		89				93	94		
*15:03, 15:16			75					80			83		85		87		89				93			
*15:04:01-15:04:05, 15:24, 15:30, 15:40, 15:66, 15:76-15:77, 15:211, 15:222, 15:236			75								83				87		89				93			
*15:05:11-15:05:12, 15:21, 15:71, 15:116, 15:244			75					80							87		89				93			
*15:05:18, 15:22, 15:37, 15:55, 15:58, 15:72, 15:125, 15:149, 15:225			75					80			83				87		89							
*15:07:01-15:07:01:02			75					80									89				93			
*15:10:01, 15:10:04, 15:208								80			83				87		89				93			
*15:11, 15:118			75			78		80			83				87		89				93			
*15:12					77			80			83				87		89	90			93			
*15:23:01-15:23:02, 15:85, 15:138, 15:158, 15:162, 15:181			75					80			83						89				93			
*15:25			75										85					90	91		93			
*15:31, 15:62, 15:169			75					80			83				87		89	90			93			
*15:36			75			78		80			83						89				93			
*15:42, 15:111		74	75					80			83				87		89				93			
*15:43			75					80									89				93	94		
*15:63								80			83						89				93			
*15:65			75								83				87		89	90	91					
*15:75			75					80			83		85		87		89	90			93		95	
*15:102	73		75								83				87		89			92				
*15:104	73		75					80			83				87		89				93			
*15:107						79		80			83				87		89				93			
*15:144					77			80									89	90			93			
*15:178			75			79		80			83				87		89				93			
*15:190			75								83				87		89							
*15:201		74	75					80			83				87		89		91		93			
*15:218			75	76				80			83				87		89				93			
*16:01:01-16:01:15, 16:01:18-16:01:41, 16:06-16:08, 16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98, 16:100, 16:105-16:106, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:180, 16:183, 16:185-16:187											83		85						90				95	
*16:01:16								80			83		85					90						95
*16:01:17, 16:15:02											83	84	85					90						95
*16:02:01-16:02:12, 16:02:14-16:02:20, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163, 16:166-16:167, 16:176, 16:179, 16:181, 16:184											83		85		87			90						95
*16:02:13											83	84	85		87			90						95
*16:04:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149, 16:182														86				90		92				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*16:10, 16:67													85					90						95
*16:20								80			83		85				89	90						95
*16:29		74												86				90		92				
*16:34			75				78				83		85					90						95
*16:35														86				90						
*16:40, 16:110							79									88								
*16:43											83		85											95
*16:45													85					90			93			95
*16:48											83			86	87			90						
*16:50		74									83		85					90						95
*16:53							79									88				92				
*16:64, 16:87											83		85				89	90						95
*16:68																		90		92				
*16:70			75								83		85		87		89	90						95
*16:88											83				87			90		92				
*16:91											83			86	87			90		92				
*16:103											83				87			90						
*16:109								80						86			89	90		92				
*16:113							79				83		85			88								95
*16:121			75				78				83		85		87			90						95
*16:139										82	83		85					90						95
*16:140											83				87									
*16:150							78							86		88		90		92				
*16:169								80					85				89	90						95
*16:170											83		85			88		90						95
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:60		74									83				87				91		93			
*17:01:04		74									83		85		87				91		93			
*17:07, 17:17, 17:23, 17:33											83				87				91		93			
*17:20		74									83								91		93			95
*17:22		74																	91		93			
*18:01:01:01-18:02:05, 18:04-18:08, 18:10-18:13, 18:15											83									92	93			95
*18:03		74									83									92	93	94		95
*18:09											83				87					92	93			
*18:14											83									92				95
A*01:110, A*02:636, A*03:210, A*11:135, A*24:137, A*32:96, A*68:160, B*07:234, B*15:01:72, B*15:125:02, B*15:315, B*18:126, B*39:138, B*44:02:44, B*44:85:01-44:85:02, B*46:01:33, B*50:56, B*51:226, B*55:30, B*58:74, B*58:107			75																					
A*24:64, A*24:96, A*24:106, B*08:16, B*08:123, B*08:145, B*44:527, B*50:76																								95
A*24:174, B*57:129											83													95
B*13:13:02, B*40:01:49-40:01:50																						94		
B*15:96, B*35:537, B*44:274				76																				
B*15:277, B*15:430, B*35:222, B*35:525, B*42:22, B*47:01:06, B*54:38															87									
B*15:539, B*39:76, B*54:18		73																						
B*46:84														86	87									
B*58:02:01:01-58:02:01:03								80																
B*58:02:02								80												92				
B*73:01:01:01-73:02																					w			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control