

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
HLA-B allele ⁴	ser. ⁵																							
*07:02:01-07:02:19, 07:05:01-07:06, 07:10, 07:12-07:15, 07:18:01- 07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34- 07:35, 07:39, 07:41-07:46, 07:49N, 07:52, 07:55- 07:59, 07:61-07:64, 07:66- 07:68, 07:70-07:72, 07:74- 07:76, 07:79-07:80, 07:82- 07:83, 07:86-07:99	B7, B7/ B42, Null, -	25																						
*07:03, 07:08, 07:16, 07:32, 07:37	B703, B7	25																						
*07:04, 07:25	B7	25													38									
*07:07, 07:77	B7, -	25																						
*07:09, 07:11, 07:17	B7	25											36											
*07:19, 07:33, 07:40, 07:53	B7/ 42																							
*07:20, 07:24	B7	25		27																				
*07:27	-	25										35								44				
*07:28	B7	25					30																	
*07:36, 07:81	B7, -	25																						
*07:38	-	25																						
*07:47	B7	25																						
*07:48, 07:51	B7																							
*07:50	-	25										35												
*07:54	B7	25																						
*07:60	-			27																				
*07:65	-	25																	42					
*07:69, 07:85	-	25																						47
*07:73	-	25															40							
*07:78	-	25										34												
*07:84	-	25																						
*08:01:01-08:01:10, 08:05, 08:08N, 08:10- 08:11, 08:15, 08:18- 08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33- 08:34, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61	B8, Null, -		26						32															
*08:02-08:03, 08:52	B8, -		26						32															
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*08:04	B8		26						32			35												
*08:06, 08:13, 08:20, 08:40	B8, –		26																					
*08:07, 08:14	B8		26																					
*08:09	B8		26						32															
*08:12:01-08:12:03, 08:16	B8		26						32															
*08:17, 08:54	–		26						32			35												
*08:21	B8		26	27					32															
*08:25	B8		26	27																				
*08:26, 08:50, 42:07 ⁶	B8, –								32															
*08:28, 08:37	B8, –		26									36												
*08:32	B8		26						32										42					
*08:35	–		26						32			36												
*08:36	–		26						32															
*08:38	–		26						32															
*08:49	–		26															41			44			
*08:55	–		26																			45		
*08:60	–		26						32															
*13:01:01-13:01:03, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29	B13, –			27	28						34							41						
*13:02:01-13:02:07, 13:08Q, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37	B13, –			27	28						34													
*13:03	B49/B15			27	28																			
*13:04, 13:10	B49/B15, B13			27	28							36												
*13:06	B13			27	28													41						
*13:07N	Null			27	28						34							41						
*13:09	–				28						34													
*13:11	B13			27	28						34													
*13:12-13:13	B13			27	28													41						
*13:14	B13			27							34													
*13:16	B13			27	28						34													
*13:21, 13:35	–			27	28																			
*13:26	–			27	28							36						41						
*13:31	–			27	28						34													
*13:36	–			27	28						34							41						
*14:01, 14:07N, 14:12, 14:14	B64, Null, –				28	29																		47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18	B65, B14, –				28		30												47					
*14:02:03	B65				28														47					
*14:02:05	B65				28		30												47					
*14:05-14:06:02, 14:13	B14, –						30												47					
*14:08, 39:25N, 39:30, 39:32-39:34, 39:47, 39:50 ⁷	B39, B16, Null, -						29												47					
*14:10	–						29									41			47					
*15:01:01:01- 15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:26N-15:27:03, 15:32, 15:34-15:35, 15:38:01- 15:38:02, 15:50, 15:56, 15:60, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96- 15:97, 15:102, 15:104- 15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128- 15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163- 15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:188N	B62, Null, –							31				36					43							
*15:01:02	B62							31		33		36						43						
*15:01:04	B62											36						43						
*15:02:01-15:02:04, 15:88, 15:112, 15:121, 15:139, 15:170	B75, –							31				36					41	43						
*15:03:01-15:03:02, 15:47, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156. 15:158. 15:173	B72, B70, –								32			35	36					43						
*15:03:03	B72							31	32			35	36					43						
*15:04	B62							31				36						43						
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*15:57	B62							31					36							43				
*15:62	-								32			35	36					41		43				
*15:64	B71								32			35	36							43				
*15:66	B62							31					36		38					43				
*15:68	-			27								35	36											
*15:69	B72								32			35	36							43		45		
*15:73	B62							31												43				
*15:74	B72								32				36							43				
*15:76	-							31					36							43				
*15:77	B62						30	31										41		43				
*15:80	B70								32				36		38			41		43				
*15:83	B62							31	32											43				
*15:85, 15:154, 15:165	B62, -							31					36					41		43				
*15:86	-							31												43				
*15:89	B77			w				31					36					41		43				
*15:91, 15:131	-								32			35	36	37						43				
*15:101	-							31					36							43				
*15:106	-				28			31	32				36					41		43				
*15:116	B62							31					36							43				
*15:124	B71								32				36		38					43				
*15:137	-							31												43				
*15:138	-	25						31					36							43				
*15:142	-		26					31					36							43				
*15:143	-							31	32				36							43				
*15:144	-							31					36					41		43				
*15:161	-								32				36	37	38					43				
*15:180	-		26						32				36							43				
*15:186	-								32				36		38					43		45		
*15:189	-						30	31					36		38					43				
*18:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:25, 18:27-18:28, 18:30-18:36, 18:38-18:40, 18:42-18:43, 18:45-18:47	B18, Null, -												36			39			42					
*18:04	B18												36			39			42					
*18:09	B18												36			39								
*18:12	B18												36			39			42					
*18:19	B18												36			39			42					
*18:22	B18												36			39		41	42					
*18:26, 18:41	B18, -															39			42					
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

*18:29, 35:32, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:99 ⁸	B35, Null, -											36					42							
*18:37	-										35	36												
*18:44	-				28							36			39		42							
*18:48	-											36												
*18:49	-											36			39									
*27:01	B27					29										40								
*27:02, 27:53, 27:57, 27:65N	B27, Null, -					29										40								
*27:03-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N	B27, Null, -					29										40								
*27:07, 27:11, 27:14-27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63	B27, -															40								
*27:08, 27:26, 27:40, 27:42, 27:44	B2708/ B27/ B7, -					29										40								
*27:12	-					29										40								
*27:16, 27:39, 27:59N	B27, Null, -					29										40								
*27:18	-					29					35					40								
*27:19	B27															40	41							
*27:23	-					29										40								
*27:24	B27															40								
*27:25	B27															40								
*27:29	B27					29					35					40								
*27:30	B27															40	41							
*27:33	-															40								
*27:41	-					29						36				40								
*27:52	-					29																		
*27:62	-															40								
*35:01:01-35:01:09, 35:01:11-35:01:19, 35:07-35:08:04, 35:11, 35:15, 35:20:01-35:20:02, 35:23, 35:29, 35:40N-35:42:01, 35:48, 35:50, 35:52, 35:54, 35:57, 35:61, 35:77-35:78, 35:82, 35:90-35:94	B35, Null, -											36					41	42						
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

*35:01:10	B35									33		36					41	42						
*35:02:01-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06, 35:12:01- 35:12:02, 35:33-35:34, 35:36, 35:39, 35:55, 35:59, 35:65Q, 35:70, 35:74, 35:83-35:85, 35:95, 35:98	B35, -																41	42						
*35:03:03, 35:56	B35																41							
*35:04:02	B35									33							41	42						
*35:05:01-35:05:02, 35:51, 35:58, 35:72, 35:89, 35:97	B35, -			27								36						42						
*35:09:01-35:09:02, 35:18, 35:31, 35:75, 35:88	B35, -																	42						
*35:10, 35:28, 35:69, 35:80	B35, -											36					41	42						
*35:13	B35																41	42						
*35:14:01-35:14:02, 35:62	B35											36					41	42						
*35:16	B35			27								36					41	42						
*35:17, 35:30	B35			27								36					41	42						
*35:19	B35											36					41							
*35:21, 35:24:01-35:24:02	B35/B 78											36					41	42						
*35:22	B35			27														42						
*35:25, 35:49	B35											36					41							
*35:26	B35					30						36					41							
*35:27	B35											36					41							
*35:35	B35											36					41	42			45			
*35:38	B35																41	42						
*35:42:02	-											36					41	42						
*35:43, 35:67, 35:79, 35:86	B35, -											36						42						
*35:44	-																	42						
*35:45, 35:71	B35											36	37				41	42						
*35:46	B35				28							36					41							
*35:47	B35				28							36					41							
*35:60	B35																	42						
*35:63	-				28							36					41							
*35:66	B35											36						42						
*35:76	B35/22											36					41	42						
*35:81, 35:96	-																41	42						
*35:87	-								32									42						
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21	B37, Null, -													37				41			44				
*37:02	B37/ B27					29								37							44				
*37:04:01-37:04:02	B37, -													37				41			44				
*37:05	-													37				41			44				
*37:07	-	25												37							44				
*37:08	B37													37				41			44				
*37:09	B37													37							44				
*37:10, 37:22	B37, -																	41			44				
*37:11	-																	41	42		44				
*37:12	B37								32					37							44				
*37:14	-																	41			44				
*38:01:01-38:02:02, 38:08- 38:09, 38:11-38:16, 38:18- 38:19, 38:21	B38, -					29																	45	46	47
*38:03	B38					29					35											45	46		
*38:04	B38																					45		47	
*38:05	B38						30															45	46	47	
*38:06	B38																					45			
*38:07	B38					29																45			
*38:10	B38					29																	46	47	
*38:17	-					29							37							44		45	46	47	
*38:20	-					29												41				45	46	47	
*38:22	-					29										40							46	47	
*39:01:01:01- 39:01:01:02L, 39:01:03- 39:01:08, 39:01:10- 39:01:11, 39:03, 39:05:01- 39:06:02, 39:09, 39:11, 39:14-39:15, 39:24, 39:26- 39:29, 39:31, 39:35, 39:37- 39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51- 39:57	B3901, B39, Null, -					29																45		47	
*39:01:09, 39:12, 39:45	B3901, B39, -																					45		47	
*39:02:01, 39:08	B3902/ B39					29					35											45			
*39:02:02, 39:13:01- 39:13:02, 39:23, 39:49	B3902, B39, -										35											45			
*39:04	B39						30															45		47	
*39:07	B16					?						36										45		47	
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	

[illegible]

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*40:25, 40:69	B40, B60				28																			
*40:26	B21									33														
*40:28	-									33								41						
*40:30, 40:34	B60				28																			
*40:31, 40:45, 40:80	B60, -				28																			
*40:38, 40:52	B60				28								36											
*40:46	B60											35					40							
*40:47	B40				28																			
*40:48	B60				28						34													
*40:58	B40				28																			
*40:59	B40				28								36											
*40:60	-				28								36											
*40:68	-																							
*40:71	B61			27																				
*40:74	B60																							
*40:76	B60				28																			
*40:83	-																							
*40:93	-											35					40							
*40:95	-																							
*40:96	-																							
*40:98	-	25																						
*41:01, 41:06-41:07, 41:12	B41, -				28				32															
*41:02:01-41:02:02, 41:11	B41, -				28				32															
*41:03:01-41:03:02	B41				28				32															
*41:04, 41:10	B41, -				28																			
*41:05	B41				28				32															
*41:08	B41	25			28				32							38								
*41:09, 45:02-45:03 ⁹	-, B45				28				32															
*42:01:01-42:02, 42:06, 42:10, 42:12	B42, -								32															
*42:04	B42								32															
*42:05:01-42:05:02	B42	25							32							38								
*42:08	B42								32															
*42:09	B42												36											
*42:11	-								32			35												
*42:13	-								32															
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

[illegible]

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*44:83	-																	41						
*44:90	-				28													41						
*44:95	-				28													41						
*44:97	-				28													41						
*45:01, 45:05, 45:07, 45:11	B45, -				28				32			35												
*45:04	B21				28				32															
*45:06	B45				28				32			35												
*45:08	-				28				32															
*45:09	-				28				32				36	37										
*45:10	-				28				32		34													
*46:01:01-46:01:04, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:22	B46, Null, -								31				36											
*46:02	B46								31				36											
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21	B46, -								31				36											
*46:06	-								31				36											
*46:11	B22								31															
*46:12, 46:20	-			27					31				36											
*46:17	-								31				36	37										
*46:18	-								31															
*47:01:01-47:01:01:02, 47:06-47:07	B47, -				28																			
*47:02	B47				28																			
*47:03	B47				28																			
*47:04-47:05	B17, -				28	29																		
*48:01:01-48:01:02, 48:07, 48:09, 48:11, 48:19-48:20, 48:22	B48, -											35												
*48:02:01-48:02:02	B48/B70											35	36					41						
*48:03:01-48:03:02	B48											35												
*48:04	B48											35												
*48:05, 48:08	B48, -	25										35												
*48:06	-																							
*48:10, 48:12, 48:17	B48, -											35												
*48:13	-											35												
*48:14	B48											35	36											
*48:15	B48	25										35												
*48:16	B48											35												
*48:18	B48											35												
*48:21	-					29						35												
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*49:01:01, 49:06, 49:08-49:10	B49, -				28																			
*49:01:02	B49				28																			
*49:02	B49				28																			
*49:03	B49																							
*49:04-49:05	B49				28								36											
*49:07	-				28						34													
*50:01:01-50:01:02, 50:04-50:05, 50:08	B50, -				28																			
*50:02	B45				28							35												
*50:06	-				28																			
*50:07	-																							
*50:09	-				28																			
*51:01:01, 51:01:03, 51:01:05-51:01:08, 51:01:10-51:02:01, 51:02:03, 51:03, 51:07:01, 51:09:01, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:24:03, 51:26-51:33, 51:38-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:88, 51:90-51:91, 51:94-51:95	B51, B5102, B5103, Null, -										33													
*51:01:02, 51:01:09, 51:02:02, 51:05, 51:07:02, 51:09:02, 51:89	B51, B5102, -										33													
*51:01:04, 51:02:04, 51:13:01, 51:24:01-51:24:02, 51:34-51:35, 51:92	B51, B5102, -																							
*51:04, 51:46, 51:56	B51, -										33							41						
*51:06, 51:59	B51, -										33													
*51:08, 51:20, 51:44N	B51, Null								32		33													
*51:10, 51:81	B51, -																							
*51:15	-																							
*51:36	B51								32															
*51:37, 51:63	B51, -										33		36											
*51:42	-																	41						
*51:45	-												36											
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

[illegible]

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*55:02:01-55:02:05, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:37, 55:39	B55, B54, –																							
*55:04, 55:08, 55:13, 55:27, 56:15, 56:19N, 56:22 ¹⁰	B55, B56, Null, –																							
*55:09, 55:24	B22, B55										34													
*55:14	-												36					41						
*55:17, 55:28	B55, –																							
*55:18	-											35												
*55:20	B22								32															
*55:21	B55																							
*55:22	B22										34													
*55:23, 55:32, 56:18 ¹¹	B55, –												36											
*55:34	–																							
*55:35	–																							
*55:40	–																							
*56:01:01-56:01:04, 56:08, 56:14, 56:16, 56:20, 56:24, 56:26-56:29	B56, Null, –																							
*56:02, 56:04	B56																							
*56:03	B22/15												36											
*56:05:01	B56										33													
*56:05:02	B56										33													
*56:06	B78										33								42					
*56:07	B56																							
*56:09	B56												36					41						
*56:10	B55																							
*56:11-56:12	B56, B22																	41						
*56:13	B56								32															
*56:17	B56																							
*56:21	–										33													
*56:23	–																							
*56:25	–																							
*57:01:01-57:01:04, 57:01:06-57:01:07, 57:06, 57:08, 57:10, 57:15-57:16, 57:18-57:20, 57:22-57:23, 57:25-57:27, 57:29-57:30	B57, –												36					41						
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
*57:01:05, 57:11, 57:13, 57:21, 57:31	B57, -												36											
*57:02:01-57:03:02, 57:05, 57:17, 57:28N	B57, Null, -																	41						
*57:04	B57						30											41						
*57:07	B57																	41						
*57:09	B57								32									41						
*57:12	B57																	41						
*57:14	-												36					41						
*57:24	-												36					41						
*58:01:01-58:01:02, 58:01:04, 58:04-58:05, 58:10N-58:15, 58:19, 58:21-58:24	B58, Null, -												36					41						
*58:01:03, 58:02, 58:06, 58:16, 58:25-58:26	B58, -												36											
*58:07	B58												36	37										
*58:08	B17/ B5									33														
*58:09	B17												36					41						
*58:17N	Null												36					41						
*58:18	-			27									36											
*58:20	-																				45			
*58:27	-																							
*58:28	-																	41						
*59:01, 59:05	B59, -																							
*59:02-59:03	B59, -																							
*59:04	-																							
*67:01:01	B67					29																45		
*67:01:02-67:02	B67, -																					45		
*73:01	B73																							
*78:01, 78:02:02-78:03, 78:07	B78, -									33								42						
*78:02:01, 78:04	B78, B35									33								42						
*78:05	-									33								42						
*78:06	-									33														
*81:01	B81																							
*81:02	B81																							
*81:03-81:04N	Null, -																							
*82:01-82:03	B82, -								32			35							43					
*83:01	-																	41						
C*01:30				28																				
C*02:06																								
C*02:23		25																						
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
C*03:05, C*03:25, C*03:27				27																				
C*03:12, 03:19																								
C*06:20																								
C*07:46									32															
C*15:02:04																								
C*15:25																								
C*16:10							30																	
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72					

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
																									ser. ⁵
48				52																			71	72	B7, B7/ B42, Null, -
48																							71	72	B703, B7
48			51	52																			71	72	B7
48				52																			71		B7, -
48				52																			71	72	B7
48				52																			71		B7/ 42
48				52																			71	72	B7
																						70		72	-
48				52																			71	72	B7
				52																		70		72	B7, -
																						70		72	-
48				w																			71	72	B7
48				52																			71	72	B7
48																							71	72	-
48				52									62										71	72	B7
48				52																			71		-
				52				56							64				w				71	72	-
48																							71	72	-
48				52																			71	72	-
48				52							59												71		-
48				52							59												71		-
48				51																			71	72	B8, Null, -
				51																		70		72	B8, -
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72		
48			51																				71	72	B8	
48																							71	72	B8, –	
48			51																				71	72	B8	
48			51			54																	71		B8	
48			51																				71		B8	
48			51											62									71	72	–	
48			51																				71	72	B8	
48																							71	72	B8	
48			51																				71	72	B8, –	
48			51																				71	72	B8, –	
48			51																	68			71	72	B8	
48			51																				71	72	–	
			51																				71	72	–	
48			51											62									71	72	–	
48			51			53																	71	72	–	
48																							71	72	–	
48			51			53																	71	72	–	
														62									70		B13, –	
														62									70		B13, –	
														62			65					70			B49/B 15	
														62								70			B49/B 15, B13	
														62		64						70			B13	
																						70			Null	
														62									71			-
														62				67				70			B13	
														62								70			B13	
														62								70			B13	
											59			62								70			B13	
														62								70			–	
														62								70			–	
											w			62							69	70			–	
									57													70			–	
48																							71		B64, Null, –	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72		
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72		

[illegible]

[illegible]

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
							w						62								69	71		B62	
													62									71		-	
																						71		B71	
													62								69	71		B62	
													62								69	71		-	
													62									71		B72	
													62				65	67		69	71			B62	
													62									71		B72	
			52						w												69	71		-	
													62								69	71		B62	
																						71		B70	
						54							62									71		B62	
													62								69	71		B62, -	
													62		64							71		-	
									57													70		B77	
													62									71		-	
			52										62								69	71		-	
													61	62								71		-	
	49												62									71		B62	
	49																					71		B71	
										w			62								69	71		-	
													62									71		-	
													62								69	71		-	
													61								69	71		-	
									57												69	71		-	
																						71		-	
																					69	71		-	
									57												69	71		-	
																				68		71		B18, Null, -	
																64				68		71		B18	
																				68		70		B18	
													62									71		B18	
																				68	69	71		B18	
							55													68		71		B18	
																				68		71		B18, -	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

																			68			71			B35
																			68			71			B35, –
																			68			71			B35
																			68			71			B35
																			68			71			B35, –
																			68			71			B35, –
													62									71			B35, –
													62									71			B35
																			68	69		71			B35
													62									71			B35
																			68			71			B35
													61									71			B35
								55												68		71			B35/B 78
																				68		71			B35
																						71			B35
																						71			B35
																			68			71			B35
																			68			71			B35
								53												68		71			B35
																64				68		71			–
																				68	69	71			B35, –
																				68	69	71			-
								53												68		71			B35
												57										71			B35
													61									71			B35
																64				68		71			B35
															61	62						71			-
																			68			71	72		B35
								52												w		71			B35/22
																				68		71			–
																				68		71	72		–
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48																				71	B39				
48												64			67					71	B39				
48																		69		71	B39				
48																		68		71	B39				
48																		68		71	B39				
48												62								71	B39				
48																			69	71	-				
48																				71	B39				
48																				71	B39				
48																				71	-				
48	49											61	62							71	B60, -				
48												61	62							71	B60				
		50										61	62							71	B61, B41, -				
		50										61	62							71	B61				
		50										61	62							71	B4005/ B50				
		50										61								71	B61				
?		50										61	62							71	B61, -				
?	49											61	62							71	B60, Null, -				
	49												62							71	B60/B 48				
		50							58		60		62						70		B61, -				
?												61	62							71	B60, B61				
?												61	62							71	72	B40, -			
48												61	62							71	72	B61			
?	49												62							71		B40/B 15			
48	49											61	62							71		B40			
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
?	49												61										71		B40, B60
?		50											61	62									71		B21
?		50					55						61	62									71		-
?													61	62									71	72	B60
?									57				61	62									71	72	B60, -
?	49												61	62									71		B60
?	49													62									71		B60
	49													62							70				B40
?	49												61	62				67					71		B60
?													61	62		64							71		B40
48													61	62									71	72	B40
48	49												61	62									71		-
?		50											61										71		-
?		50											61	62			65	67					71		B61
?	49	50											61	62									71		B60
	49													62											B60
?													61	62									71		-
?														62									71		-
?		50								59			61	62									71		-
		50												62							70				-
?		50											61	62									71	72	-
48			51			54							61	62									71		B41, -
48			51										61	62									71	72	B41, -
48			51										61	62									71		B41
48			51										61	62									71	72	B41, -
48			51			54							?	62									71		B41
48			51										61	62									71		B41
48						54							61	62									71		-, B45
48			51	52																			71	72	B42, -
48			51	52		54																	71		B42
48			51	52																			71	72	B42
48				52																			71		B42
48			51	52																			71	72	B42
48			51											62									71	72	-
48			51	52	53																		71	72	-
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

					53								62 63					70			B44, Null, –				
					53								62					70			B44, Null, –				
					53								62					70			B44, –				
					53 54					60			63				68	70			B44				
					53								62 63					70			B44				
					53							61 62 63							71		B45				
					53							62			65			70			B44				
					53								63					70			B44				
					53							62						70			B44				
					54							62			65			70			B12				
					53							62 63						70			B47/B 44, –				
					53							62						70			B44				
					54				58			62			65			70			B44				
					53 54							62 63						70			B44				
					53				58	60		62						70			B44				
49					53							62						70			B44				
												62						70			B44, –				
					53							62						70			B44				
					53							62 63						70			B44, –				
					53							62 63						70			B44				
					53						61 62								71		B12, –				
					53						62							70			B44, –				
					53						62 63							70			B44				
					53						62							70			B44				
					53						62							70			–				
					53						62							70			–				
					53						62							70			–				
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
					53									62	63							70			-
				52	53			56					61	62	63								71		-
					53					58		60		62	63							70			-
				52	53			56						62	63							70			-
48						54							61	62			65						71		B45, -
48						54							61	62			65						71		B21
48				52		54							61	62			65						71		B45
48						54							61	62		64							71		-
48						54							61	62									71		-
48						54							61	62									71		-
				52			55														69		71		B46, Null, -
				w			55														69		71		B46
				52			55																71		B46, -
	49			52			55																71		-
				52			55										65						71		B22
				52			55														69		71		-
				52			55														69		71		-
				52			55									64							71		-
								56						62								70			B47, -
								56					61	62									71		B47
								56					w	62								70	w		B47
													62									70			B17, -
48									57					62									71	72	B48, -
48														62									71		B48/B 70
48	49													62									71		B48
48									57					62									71		B48
48														62									71	72	B48, -
48									57														71	72	-
48														62									71	72	B48, -
48														62									71		-
48														62									71	72	B48
48									57					62									71	72	B48
48									57														71	72	B48
									57					62								70		72	B48
48														62									71		-
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
										58	59			62			65					70			B49, –
										58				62			65					70			B49
											59			62			65					70			B49
											59			62			65					70			B49
										58	59			62								70			B49
										58	59			62								70			–
48											59		61	62			65						71		B50, –
48											59		61	62			65						71		B45
48											59			62			65						71		–
48											59		61	62			65						71		–
48											59		61	62									71		–
						54					59	60								68		70			B51, B5102, B5103, Null, –
											59	60								68		70			B51, B5102, –
						54					59	60								68		70			B51, B5102, –
						54	55					60								68		70			B51, –
						54						60								68		70			B51, –
						54					59	60								68		70			B51, Null
						54						60								68		70			B51, –
												60					65			68		70			–
						54					59	60								68		70			B51
						54					59	60								68		70			B51, –
						53	54					60								68		70			–
						54						60								68		70			–
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
										59										68		70			–
						54				w	60										68	69	70		–
						54					60						65		67	68		70			–
						54					60										68	70			–
						54					59	60									68	70			–
											59	60									68	70			–
											59	60		62								70			B52, –
						54					59	60		62								70			B52
											59	60		62								70			B52, –
											59	60		62								70			–
						54					59	60		62								70			–
											59			62								70			–
											60										68		70		B53, –
						55					60										68		70		B53
																					68		70		B53/B 35
											60										68		70		B53, –
						54	55				60										68		70		B53/B 51
											60										68		70		–
											60										68		70	72	–
											60		62									70			–
					53						60										68		70		–
				52				56						63	64								71		B54, Null, –
				52				56						63		65							71		B22
				52				56						63									71		B54, –
				52				56						63							69		71		–
				52				56						63									71		–
				52				56		58				63	64								70		–
				52				56			59			63	64								71		–
				52				56			59				64								71		B55, –

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
				52				56								64							71		B55, B54, –
				52				56															71		B55, B56, Null, –
				52				56			59												71		B22, B55
				52				56									66						71		-
				52				56			59												71		B55, –
													62		64								71		-
			51	52		54		56															71		B22
				52				56			w					64					69		71		B55
				52				56															71		B22
				52				56															71		B55, –
				52				56					62		64								71		–
				52				56	w							64							71		–
				52				56			59		62		64								71		–
				52				56									65						71		B56, Null, –
				52				56									65		67				71		B56
				52				56													69		71		B22/15
				52				56			59				64								71		B56
				52				56			59												71		B56
				52				56			59				64				w				71		B78
				52				56									65					70			B56
				52				56															71		B56
				52				56								64		67					71		B55
				52				56															71		B56, B22
				52		54		56									65						71		B56
																	65						71		B56
				52				56		58	59											70			–
																64							71		–
				52				56		59							65						71		–
												60						66				70			B57, –
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
												60						66				70			B57, –
												60						66				70			B57, Null, –
												60						66				70			B57
					53							60						66				70			B57
					53							60						66				70			B57
																		66					71		B57
							55					60						66				70			–
					53							60						66				70			–
												60							67			70			B58, Null, –
												60							67			70			B58, –
												60							67			70			B58
						54					59	60							67			70			B17/ B5
							55					60							67			70			B17
																			67			70			Null
												60							67			70			–
												60							67			70			–
												60							67			70			–
												60							67			70			–
										58						64						70			B59, –
										58												70			B59, –
										58							65					70			–
48				52																			71		B67
48				52																			71		B67, -
48				52																	69		71		B73
											59					64				68			71		B78, –
											59									68			71		B78, B35
											59		62										71		-
											59		62										71		–
				52				56	57														71	72	B81
48				52					57														71	72	B81
				52					57														71	72	Null, -
				52				56									65						71		B82, –
				52	53			56						63									71		-
																			w						
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
																	65								
																				69					
																64									
					53																				
48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	

[illegible]

Well No.

Well No.
HLA-B allele ⁴
*07:02:01-07:02:19, 07:05:01-07:06, 07:10, 07:12-07:15, 07:18:01- 07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34- 07:35, 07:39, 07:41-07:46, 07:49N, 07:52, 07:55- 07:59, 07:61-07:64, 07:66- 07:68, 07:70-07:72, 07:74- 07:76, 07:79-07:80, 07:82- 07:83, 07:86-07:99
*07:03, 07:08, 07:16, 07:32, 07:37
*07:04, 07:25
*07:07, 07:77
*07:09, 07:11, 07:17
*07:19, 07:33, 07:40, 07:53
*07:20, 07:24
*07:27
*07:28
*07:36, 07:81
*07:38
*07:47
*07:48, 07:51
*07:50
*07:54
*07:60
*07:65
*07:69, 07:85
*07:73
*07:78
*07:84
*08:01:01-08:01:10, 08:05, 08:08N, 08:10- 08:11, 08:15, 08:18- 08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33- 08:34, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61
*08:02-08:03, 08:52
Well No.

Well No.
*08:04
*08:06, 08:13, 08:20, 08:40
*08:07, 08:14
*08:09
*08:12:01-08:12:03, 08:16
*08:17, 08:54
*08:21
*08:25
*08:26, 08:50, 42:07 ⁶
*08:28, 08:37
*08:32
*08:35
*08:36
*08:38
*08:49
*08:55
*08:60
*13:01:01-13:01:03, 13:17, 13:20, 13:22:01- 13:23, 13:25, 13:28-13:29
*13:02:01-13:02:07, 13:08Q, 13:15, 13:18- 13:19, 13:27, 13:30, 13:32 13:34, 13:37
*13:03
*13:04, 13:10
*13:06
*13:07N
*13:09
*13:11
*13:12-13:13
*13:14
*13:16
*13:21, 13:35
*13:26
*13:31
*13:36
*14:01, 14:07N, 14:12, 14:14
Well No.
Well No.

*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18
*14:02:03
*14:02:05
*14:05-14:06:02, 14:13
*14:08, 39:25N, 39:30, 39:32-39:34, 39:47, 39:50 ⁷
*14:10
*15:01:01:01- 15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:26N-15:27:03, 15:32, 15:34-15:35, 15:38:01- 15:38:02, 15:50, 15:56, 15:60, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96- 15:97, 15:102, 15:104- 15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128- 15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163- 15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:188N
*15:01:02
*15:01:04
*15:02:01-15:02:04, 15:88, 15:112, 15:121, 15:139, 15:170
*15:03:01-15:03:02, 15:47, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156. 15:158. 15:173
*15:03:03
*15:04
Well No.
Well No.

*15:05:01-15:06, 15:33, 15:39-15:40, 15:65, 15:78:01-15:78:02, 15:84, 15:107, 15:136, 15:141, 15:155, 15:179, 15:185
*15:07, 15:126
*15:08, 15:11:01- 15:11:03, 15:15, 15:148, 15:191
*15:09
*15:10:01, 15:37, 15:90, 15:99, 15:133
*15:10:02
*15:11:04
*15:12, 15:14, 15:19
*15:13
*15:16, 15:67, 15:95
*15:17:01:01-15:17:02, 15:162, 15:168, 15:177
*15:18:01-15:18:02, 15:18:04, 15:51-15:52, 15:72, 15:93, 15:108, 15:114, 15:119, 15:134, 15:153. 15:176
*15:18:03
*15:20, 15:25:01-15:25:02
*15:21, 15:44
*15:23, 15:115
*15:24, 15:87, 15:157
*15:28, 15:120, 15:149N
*15:29
*15:30, 15:58, 15:63
*15:31
*15:36
*15:42
*15:43
*15:45
*15:46
*15:48, 15:150, 15:188
*15:53
*15:54
*15:55
Well No.

Well No.
*15:57
*15:62
*15:64
*15:66
*15:68
*15:69
*15:73
*15:74
*15:76
*15:77
*15:80
*15:83
*15:85, 15:154, 15:165
*15:86
*15:89
*15:91, 15:131
*15:101
*15:106
*15:116
*15:124
*15:137
*15:138
*15:142
*15:143
*15:144
*15:161
*15:180
*15:186
*15:189
*18:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:25, 18:27-18:28, 18:30-18:36, 18:38-18:40, 18:42-18:43, 18:45-18:47
*18:04
*18:09
*18:12
*18:19
*18:22
*18:26, 18:41
Well No.
Well No.

*18:29, 35:32, 35:37, 35:53N, 35:64, 35:68:01- 35:68:02, 35:99 ⁸
*18:37
*18:44
*18:48
*18:49
*27:01
*27:02, 27:53, 27:57, 27:65N
*27:03-27:06, 27:09- 27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35 27:37, 27:45-27:46, 27:48- 27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N
*27:07, 27:11, 27:14- 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63
*27:08, 27:26, 27:40, 27:42, 27:44
*27:12
*27:16, 27:39, 27:59N
*27:18
*27:19
*27:23
*27:24
*27:25
*27:29
*27:30
*27:33
*27:41
*27:52
*27:62
*35:01:01-35:01:09, 35:01:11-35:01:19, 35:07- 35:08:04, 35:11, 35:15, 35:20:01-35:20:02, 35:23, 35:29, 35:40N-35:42:01, 35:48, 35:50, 35:52, 35:54, 35:57, 35:61, 35:77 35:78, 35:82, 35:90-35:94
Well No.
Well No.

*35:01:10
*35:02:01-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06, 35:12:01- 35:12:02, 35:33-35:34, 35:36, 35:39, 35:55, 35:59, 35:65Q, 35:70, 35:74, 35:83-35:85, 35:95, 35:98
*35:03:03, 35:56
*35:04:02
*35:05:01-35:05:02, 35:51, 35:58, 35:72, 35:89, 35:97
*35:09:01-35:09:02, 35:18, 35:31, 35:75, 35:88
*35:10, 35:28, 35:69, 35:80
*35:13
*35:14:01-35:14:02, 35:62
*35:16
*35:17, 35:30
*35:19
*35:21, 35:24:01-35:24:02
*35:22
*35:25, 35:49
*35:26
*35:27
*35:35
*35:38
*35:42:02
*35:43, 35:67, 35:79, 35:86
*35:44
*35:45, 35:71
*35:46
*35:47
*35:60
*35:63
*35:66
*35:76
*35:81, 35:96
*35:87
Well No.
Well No.

*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21
*37:02
*37:04:01-37:04:02
*37:05
*37:07
*37:08
*37:09
*37:10, 37:22
*37:11
*37:12
*37:14
*38:01:01-38:02:02, 38:08 38:09, 38:11-38:16, 38:18- 38:19, 38:21
*38:03
*38:04
*38:05
*38:06
*38:07
*38:10
*38:17
*38:20
*38:22
*39:01:01:01- 39:01:01:02L, 39:01:03- 39:01:08, 39:01:10- 39:01:11, 39:03, 39:05:01- 39:06:02, 39:09, 39:11, 39:14-39:15, 39:24, 39:26- 39:29, 39:31, 39:35, 39:37 39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51 39:57
*39:01:09, 39:12, 39:45
*39:02:01, 39:08
*39:02:02, 39:13:01- 39:13:02, 39:23, 39:49
*39:04
*39:07
Well No.
Well No.

*39:10, 39:16, 39:20
*39:17
*39:18
*39:19:01
*39:19:02
*39:22
*39:36
*39:39
*39:42
*39:43
*40:01:01-40:01:06, 40:01:08-40:01:10, 40:07, 40:51, 40:66, 40:73, 40:79
*40:01:07, 40:72:01- 40:72:02, 40:77
*40:02:01-40:02:10, 40:04, 40:06:01:01- 40:06:02, 40:39-40:40, 40:70, 40:89, 40:94, 40:97
*40:03, 40:20
*40:05
*40:08
*40:09, 40:11, 40:18, 40:24, 40:27, 40:29, 40:35, 40:37, 40:44, 40:50, 40:56-40:57, 40:64, 40:75, 40:78, 40:82, 40:85 40:86, 40:90-40:91, 40:99
*40:10, 40:22N, 40:33, 40:36, 40:42-40:43, 40:49, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87-40:88, 40:92
*40:12
*40:13, 40:19
*40:14:01-40:14:03, 40:53
*40:15, 40:32
*40:16
*40:21
*40:23
Well No.

Well No.
*40:25, 40:69
*40:26
*40:28
*40:30, 40:34
*40:31, 40:45, 40:80
*40:38, 40:52
*40:46
*40:47
*40:48
*40:58
*40:59
*40:60
*40:68
*40:71
*40:74
*40:76
*40:83
*40:93
*40:95
*40:96
*40:98
*41:01, 41:06-41:07, 41:12
*41:02:01-41:02:02, 41:11
*41:03:01-41:03:02
*41:04, 41:10
*41:05
*41:08
*41:09, 45:02-45:03 ⁹
*42:01:01-42:02, 42:06, 42:10, 42:12
*42:04
*42:05:01-42:05:02
*42:08
*42:09
*42:11
*42:13
Well No.
Well No.

*44:02:01-44:02:13, 44:11, 44:19N, 44:21- 44:24, 44:27, 44:33-44:34, 44:48, 44:52N, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93
*44:03:01-44:03:06, 44:03:08-44:05:02, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35-44:36, 44:38- 44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76-44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98
*44:03:07, 44:07, 44:81
*44:06
*44:08
*44:09
*44:10
*44:12
*44:14
*44:15
*44:16, 44:91
*44:17, 44:43
*44:18
*44:20
*44:25, 44:50
*44:31
*44:37, 44:64:01-44:64:02
*44:40
*44:41, 44:53, 44:99
*44:44
*44:46, 44:75
*44:47, 44:65
*44:49
*44:54
*44:57
*44:60
*44:62
Well No.

Well No.
*44:83
*44:90
*44:95
*44:97
*45:01, 45:05, 45:07, 45:11
*45:04
*45:06
*45:08
*45:09
*45:10
*46:01:01-46:01:04, 46:04 46:05, 46:07N, 46:09- 46:10, 46:14-46:16, 46:22
*46:02
*46:03, 46:08, 46:13:01- 46:13:03, 46:19, 46:21
*46:06
*46:11
*46:12, 46:20
*46:17
*46:18
*47:01:01:01-47:01:01:02, 47:06-47:07
*47:02
*47:03
*47:04-47:05
*48:01:01-48:01:02, 48:07, 48:09, 48:11, 48:19 48:20, 48:22
*48:02:01-48:02:02
*48:03:01-48:03:02
*48:04
*48:05, 48:08
*48:06
*48:10, 48:12, 48:17
*48:13
*48:14
*48:15
*48:16
*48:18
*48:21
Well No.

Well No.
*49:01:01, 49:06, 49:08-49:10
*49:01:02
*49:02
*49:03
*49:04-49:05
*49:07
*50:01:01-50:01:02, 50:04-50:05, 50:08
*50:02
*50:06
*50:07
*50:09
*51:01:01, 51:01:03, 51:01:05-51:01:08, 51:01:10-51:02:01, 51:02:03, 51:03, 51:07:01, 51:09:01, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:24:03, 51:26-51:33, 51:38-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:88, 51:90-51:91, 51:94-51:95
*51:01:02, 51:01:09, 51:02:02, 51:05, 51:07:02, 51:09:02, 51:89
*51:01:04, 51:02:04, 51:13:01, 51:24:01-51:24:02, 51:34-51:35, 51:92
*51:04, 51:46, 51:56
*51:06, 51:59
*51:08, 51:20, 51:44N
*51:10, 51:81
*51:15
*51:36
*51:37, 51:63
*51:42
*51:45
Well No.

Well No.
*51:54, 51:78
*51:61
*51:62
*51:64
*51:68
*51:93
*52:01:01, 52:01:03, 52:04-52:05, 52:07-52:08, 52:10-52:13, 52:15, 52:17- 52:18
*52:01:02, 52:01:04, 52:02-52:03, 52:06:01- 52:06:02, 52:09
*52:01:05, 52:14
*52:16
*52:19
*52:20
*53:01:01-53:01:05, 53:05, 53:08:01-53:08:02, 53:10, 53:16, 53:18, 53:20 53:21, 53:23
*53:02
*53:03, 53:09, 53:11- 53:13
*53:04, 53:07, 53:19
*53:06
*53:14
*53:15
*53:17
*53:22
*54:01-54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21- 54:22
*54:03
*54:04, 54:11, 54:15
*54:06
*54:09, 54:14
*54:12
*54:20
*55:01:01-55:01:06, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38

Well No.
Well No.
*55:02:01-55:02:05, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:37, 55:39
*55:04, 55:08, 55:13, 55:27, 56:15, 56:19N, 56:22 ¹⁰
*55:09, 55:24
*55:14
*55:17, 55:28
*55:18
*55:20
*55:21
*55:22
*55:23, 55:32, 56:18 ¹¹
*55:34
*55:35
*55:40
*56:01:01-56:01:04, 56:08, 56:14, 56:16, 56:20, 56:24, 56:26-56:29
*56:02, 56:04
*56:03
*56:05:01
*56:05:02
*56:06
*56:07
*56:09
*56:10
*56:11-56:12
*56:13
*56:17
*56:21
*56:23
*56:25
*57:01:01-57:01:04, 57:01:06-57:01:07, 57:06, 57:08, 57:10, 57:15-57:16, 57:18-57:20, 57:22-57:23, 57:25-57:27, 57:29-57:30
Well No.

Well No.
*57:01:05, 57:11, 57:13, 57:21, 57:31
*57:02:01-57:03:02, 57:05, 57:17, 57:28N
*57:04
*57:07
*57:09
*57:12
*57:14
*57:24
*58:01:01-58:01:02, 58:01:04, 58:04-58:05, 58:10N-58:15, 58:19, 58:21-58:24
*58:01:03, 58:02, 58:06, 58:16, 58:25-58:26
*58:07
*58:08
*58:09
*58:17N
*58:18
*58:20
*58:27
*58:28
*59:01, 59:05
*59:02-59:03
*59:04
*67:01:01
*67:01:02-67:02
*73:01
*78:01, 78:02:02-78:03, 78:07
*78:02:01, 78:04
*78:05
*78:06
*81:01
*81:02
*81:03-81:04N
*82:01-82:03
*83:01
C*01:30
C*02:06
C*02:23
Well No.

Well No.
C*03:05, C*03:25, C*03:27
C*03:12, 03:19
C*06:20
C*07:46
C*15:02:04
C*15:25
C*16:10
Well No.