

A-B-DR 2019 8H7 B

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
HLA-B allele ^{1,2}	ser. ³																																																	
*07:02:01-01-07:02:31, 07:02:33-07:02:38, 07:02:40-07:02:66, 07:05:01-01-07:06:03, 07:10, 07:13, 07:15, 07:21-07:23, 07:26, 07:30-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:49N, 07:51-07:52, 07:55-07:59, 07:61-07:64, 07:66-07:67N, 07:70-07:72, 07:74-07:76, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:136:02, 07:138, 07:140-07:145, 07:147-07:148, 07:150-07:161N, 07:166-07:173, 07:175-07:178, 07:181N, 07:183-07:185, 07:187-07:191, 07:194, 07:200-07:201N, 07:203-07:205, 07:208, 07:210-07:213, 07:215-07:217, 07:221, 07:224, 07:229-07:231N, 07:233-07:234, 07:238, 07:241-07:247, 07:249-07:254, 07:256, 07:258-07:261, 07:263-07:270, 07:272N, 07:274-07:278, 07:280-07:286, 07:288-07:300, 07:302-07:305, 07:307-07:335	B7, B42, Null, -	25																							48				52																		71	72		
*07:02:32	-	25				29																			48				52																		71	72		
*07:02:39, 07:07, 07:77, 07:193, 07:196, 07:198, 07:214, 07:239	B7, -	25																						48				52																			71			
*07:03, 07:08:01-07:08:02, 07:16, 07:32, 07:37:01-07:37:02, 07:207, 07:228:01-07:228:02, 07:232	B703, B7, -	25																						48																								71	72	
*07:04:01-07:04:02, 07:25, 07:146	B7, -	25													38									48			51	52																			71	72		
*07:09:01-07:09:02, 07:11, 07:17, 07:287	B7, -	25										36												48			52																				71	72		
*07:12, 07:14, 07:18:01-07:18:02, 07:163, 07:199, 07:225, 07:227, 07:248	B7, -	25																						48			52		55																		71	72		
*07:19, 07:33:01-07:33:03, 07:40, 07:53, 07:122, 07:139:01-07:139:02, 07:165	B7, 42, -																							48			52																				71			
*07:20, 07:24, 07:131	B7, -	25		27																				48			52																				71	72		
*07:27	-	25										35		37						44					48			52																		70	72			
*07:28	B7	25				30																		48			52																				71	72		
*07:29, 07:186	B7, -	25																						48			52												65							71	72			

*13:02:01-13:02:05, 13:02:07-13:02:24, 13:08, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40, 13:42, 13:44-13:45, 13:47, 13:49N, 13:54-13:56:02N, 13:64-13:66, 13:68-13:70, 13:74-13:75, 13:81-13:82, 13:84-13:85, 13:88-13:89, 13:91, 13:94, 13:96, 13:99-13:101, 13:103N-13:105, 13:111-13:119, 13:121, 13:123Q	B13, Null, -			27	28				34											62								70
*13:02:06	B13			27	28				34						50						62							70
*13:03, 13:120	B49, B15, -			27	28																62			65				70
*13:04, 13:10	B49, B15, B13			27	28					36											62							70
*13:06	B13			27	28							41									62		64					70
*13:07N	Null			27	28				34			41																70
*13:09, 13:67	-				28				34												62							70
*13:11	B13			27	28				34												62				67			70
*13:12:01-13:13:02, 13:21, 13:51, 13:79, 13:86	B13, -			27	28							41									62							70
*13:14	B13			27					34												62							70
*13:16	B13			27	28				34											59		62						70
*13:26:01	-			27	28					36			41		50						62							70
*13:26:02	-			27	28					36			41								62							70
*13:31	-			27	28				34											w	62					69	70	
*13:35, 13:59, 13:71, 13:108, 44:135, 44:158, 44:184 ⁴	-			27	28																62							70
*13:36	-			27	28				34			41							57									70
*13:39	-				28				34			41									62							71
*13:41	-			27	28				34												62					69	70	
*13:46	-			27	28			32								54					62							70
*13:48	-			27	28																							

*15:420	-					31					38		41	43							57									71
*15:434, 15:436	-								36		38			43																71
*15:449	-					31			36					43	45								62							71
*15:462	-					31	33				38			43										66				70		
*18:01:01:01-18:03:01, 18:04:02-18:08, 18:10-18:11, 18:13, 18:15, 18:17N-18:18:01:02, 18:20-18:21, 18:23N-18:25, 18:27-18:28, 18:30-18:36, 18:38-18:40, 18:42-18:43, 18:45-18:47, 18:49-18:51, 18:53, 18:55, 18:57:01-18:60, 18:62-18:65, 18:68, 18:70-18:71, 18:73-18:91, 18:93-18:101, 18:103-18:104, 18:106, 18:109-18:121, 18:123-18:124, 18:126-18:131:01:02, 18:133-18:135, 18:137-18:147, 18:151-18:158	B18, Null, -								36		39		42												68			71		
*18:03:02, 18:66	-								36		39																		71	
*18:04:01	B18								36		39		42										64			68		71		
*18:09	B18								36		39															68	70			
*18:12:01-18:12:02	B18, -								36		39		42										62					71		
*18:14, 18:44:01-18:44:02	B18, -			28					36		39		42													68		71		
*18:19	B18								36		39		42													68	69	71		
*18:22, 18:105	B18, -								36		39	41	42			50		55								68		71		
*18:26, 18:41, 18:61, 18:122, 18:148	B18, -										39		42													68		71		
*18:29, 18:72:01-18:72:03, 18:92, 18:102, 35:32:01-35:32:02, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:99, 35:118-35:119, 35:174, 35:369 ^d	B35, Null, -								36				42													68		71		
*18:37:01-18:37:02	-							35	36		39																	71		
*18:48	-								36		39					50							61					71		
*18:52	-								36		39		42			50										68		71		
*18:54	-								36	37	39			44												68	70			
*18:56, 18:108, 18:132	-								36	37	39		42													68		71		
*18:67, 18:136	-								36		39											60				68	70			
*18:69, 18:149-18:150	-								36		39	41	42				55									68		71		
*18:107	-			28					36														61					71		
*18:125	-								35	36		39												64				71		
*27:01	B27				29							40					52		w				62				70			

A-B-DR 2019 8H7 B

*46:67	-					31				36			41						52		55											69	71																
*46:71	-									36									52		55										69	71																	
*47:01:01-02-47:01:02, 47:06-47:10	B47, -				28													50			56										70																		
*47:01:03	-				28																56										70																		
*47:02	B47				28													50			56				61	62					71																		
*47:03	B47				28													50			56				w	62					70	w																	
*47:04	-				28	29																									70																		
*47:05	B17				28	29												50													70																		
*48:01:01-01-48:01:07, 48:07, 48:09, 48:19-48:20, 48:22, 48:27, 48:30, 48:32-48:38, 48:42-48:47	B48, -									35							48						57									71	72																
*48:02:01-48:02:03, 48:25	B48, B70, -									35	36			41			48															71																	
*48:03:01-48:03:02, 48:40	B48, -									35							48	49														71																	
*48:04:01-48:04:02, 48:24	B48, -									35							48						57									71																	
*48:05, 48:08	B48, -	25								35							48															71	72																
*48:06	-																48						57									71	72																
*48:10	B48									35							48															71	72																
*48:11	B48									35							48				53		57									71	72																
*48:12	B48				28					35							48															71	72																
*48:13	-									35							48															71																	
*48:14	B48				28					35	36						48															71	72																
*48:15	B48	25								35							48						57									71	72																
*48:16, 48:39	B48, -									35							48						57									71	72																
*48:17	-									35							48	49				55										71	72																
*48:18	B48									35													57	58								70	72																
*48:21, 48:26	-				29					35							48																71																
*48:23	-									35	36						48	49															71																
*48:28	-									35				39			48						57									71	72																
*48:29	-				28					35							48						57									71	72																
*48:31	-																48						57									71	72																
*48:41	-																48						57									71	72																
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72

*55:09, 55:24	B22, B55					34								52		56		59								71
*55:14	-						36			41				52		56					66					71
*55:17, 55:28, 55:58	B55, -													52		56		59								71
*55:18	-						35													62	64					71
*55:20	B22				32									51	52	54	56									71
*55:21	B55													52		56		w			64			69	71	
*55:22	B22					34								52		56										71
*55:23, 55:32, 56:18, 56:31-56:32, 56:50 ⁴	B55, -						36							52		56										71
*55:34	-													52		56				62	64					71
*55:40	-													52		56		59		62	64					71
*55:51	-		27											52		56										71
*55:56	-				32									52		54	56					64				71
*55:68	-						36							52		56		59								71
*55:78	-																	59			64				71	
*55:90	-													52		56		59			64				70	
*56:01:01-56:01:04, 56:01:06-56:01:12, 56:01:14, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:36, 56:38N, 56:40-56:41, 56:47-56:49, 56:52-56:53, 56:55:01:02, 56:58-56:59, 56:61-56:62	B56, Null, -													52		56						65				71
*56:02:01-56:02:02, 56:04:01-56:04:02, 56:51, 56:56-56:57	B56, -													52		56						65	67			71
*56:03	B22, 15						36							52		56								69	71	
*56:05:01	B56					33								52		56		59			64				71	
*56:05:02	B56					33								52		56		59			64				71	
*56:06	B78					33				42				52		56		59			64		w		71	
*56:07, 56:46	B56, -													52		56					65			70		
*56:09	B56						36			41				52		56									71	
*56:10	B55													52		56					64		67		71	
*56:11-56:12	B56, B22									41				52		56									71	
*56:13	B56				32									52		54	56					65			71	
*56:17, 56:33	B56, -																				65			71		
*56:20:02	-											48		52		56					65				71	
*56:21	-					33								52		56		58	59						70	

