

*11:07:02, 11:372	-			4					11								19								
*11:10, 11:183, 11:191, 11:269	A11, -			4			8		11																
*11:11, 34:19	Undefined, -			4					11	12															
*11:17, 11:40, 11:223	Undefined, -			4					10	11															
*11:113, 11:162	-			4					11				15												
*11:116, 11:140, 11:199:01-11:199:02, 11:222, 66:23	-			4					11			13													
*11:130	-			3	4				11			14							23						
*11:139	-			4	5		7		11																
*11:155	-			4					11										22						
*11:199:03	-			4					11		13					20			22						
*11:226	-	1		4					11										22						
*11:306, 11:440	-			4					11				16												
*11:325, 24:555 ^d	-			4	5				11																
*23:01:01:01-23:10:02, 23:12-23:13, 23:15-23:56, 23:58-23:63, 23:65, 23:67-23:68, 23:70-23:98, 23:100-23:103N, 23:106N-23:112, 23:114-23:127, 23:129-23:136, 23:138-23:140, 23:142N-23:150	A23, Null, -					5	6																		
*23:11N, 31:29	Null, -					5	6						16												
*23:14:01-23:14:02, 23:104-23:105, 23:141, 24:24, 24:71, 24:315, 24:392, 24:527, 24:569 ^d	A23, A9, A24, -					5	6	7																	
*23:57	-					5													21						
*23:64	-					5	6								17										
*23:66, 23:99, 23:128, 23:137, 24:14:01:01-24:15, 24:51-24:53, 24:57, 24:64, 24:94, 24:114, 24:138, 24:188, 24:222N, 24:228, 24:291, 24:296, 24:304, 24:316, 24:324, 24:412, 24:481, 24:515, 24:574, 24:610, B*51:421, B*53:72, C*04:01:03 ^d	A24, Null, Cw4, -					5																			
*24:02:01:01-24:11N, 24:13:01-24:13:03, 24:17:01:01-24:17:02, 24:20:01:01-24:23, 24:25-24:43, 24:45N-24:50, 24:54-24:56:02, 24:58-24:63, 24:66-24:70, 24:72-24:81, 24:83N-24:91, 24:93, 24:95-24:103, 24:105-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:203, 24:205-24:210, 24:212, 24:214-24:221, 24:223-24:227, 24:229-24:242, 24:244-24:290, 24:292-24:295, 24:297-24:303N, 24:305-24:314, 24:317-24:323N, 24:325-24:391, 24:393-24:411, 24:413-24:457, 24:459:01-24:470, 24:472-24:480, 24:482-24:514N, 24:516-24:526, 24:528-24:550, 24:552-24:554, 24:556N, 24:558-24:568N, 24:570-24:573, 24:575-24:585:02, 24:587-24:608N, 24:611-24:645, 24:647-24:654	A24, Low A24, A2403, A24, A9, Null, A24, A9, -					5	7																		
*24:18, 24:204, 24:213	A24, A3, -					5	7												23						
*24:19, 24:44	A9					5	7		11																
*24:82	-					5	7							18					21						
*24:92, 24:458, 24:551, 24:557	-			3		5																			
*24:104, 24:243	-					5	7							17											
*24:586N	Null					5	7							16											
*24:609	-					5	7												24						
*24:646	-					5	7							16	18										
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ^{1,2}	ser ³																								
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*25:01:01:01-25:02, 25:04, 25:07-25:12N, 25:14-25:16, 25:18-25:29, 25:31-25:33, 25:35, 25:37-25:42N, 25:44-25:45, 25:47-25:49N, 25:51, 25:53, 25:55-25:58, 25:60-25:64, 25:66-25:70, 25:72-25:76:02, 25:79-25:80, 25:84-25:87, 25:90-25:94	A25, Null, -								8	9				13											
*25:03	A25								8	9				13				17							
*25:05	A25								8	9		12		13							20				
*25:06	A25								8	9		12													
*25:13	-								8	9				13											24
*25:17, 25:34, 25:54, 25:71	-								8	9				13											
*25:30	-								8	9				13				17			21				
*25:36, 25:50, 25:52, 25:65	-								8	9															
*25:43	-				4				8	9															
*25:46	-													13											24
*25:59	-			1	4							10									20				
*25:78, 26:57, 66:15 ^d	-								8					13											
*25:81	-							5						13											
*25:82	-							5						13											24
*25:83	-											9		13											24
*25:88	-								8	9				13							20				

*68:05, 68:15, 68:20, 68:109, 68:136, 68:255	A68, -								8										20		23				
*68:08:02, 68:105:01:01-68:105:01:02, 68:157, 68:276Q	A68, -								8										20	21					
*68:11N, 68:41, 68:230	Null, -								8										20		22				
*68:13:01, 68:66:01-68:66:02, 68:134	A68, -				4				8										20						
*68:19	A68								8				14						20		23				
*68:25, 68:262	A68, -								8									19	20						
*68:29	A68								8								18	20	21						
*68:45, 68:117, 68:322	-						7	8										20							
*68:71	-							8				13						20							
*68:84, 68:227	-							8	10									20							
*68:103:01-68:103:02	-				3			8										20							
*68:130:01	-							8				12						20							
*68:130:02	-							8					14					20							
*68:178	-				4			8										20							
*68:249	-					5		8										20							
*68:250	-							8									18	20	21						
*68:271	-							8				11						20							
*68:294	-							8				12						20		24					
*68:330	-						7	8										20							
*69:01:01:01-69:01:06, 69:03-69:06, 69:08-69:09	A69, -							8											21						
*69:02	-							8				11							21						
*69:07	-							8					13						21						
*74:01:01:01-74:01:02, 74:01:04-74:03:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:29, 74:31-74:32N, 74:34, 74:36, 74:38-74:42, 74:45-74:47N, 74:49-74:53	A74, Undefined, Null, -																	19							
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ^{1,2}	ser ³																								
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*74:01:03	-												12							19					24
*74:04, 74:21, 74:43	A74, -																			18	19				
*74:07, 74:33, 74:48	A74, -																			17	19				
*74:13	-													13		16				19					
*74:19	-				4															19					
*74:23	-				3															19					
*74:30	-												11							19					
*74:35	-																			19		21			
*74:37	-													13						19					
*74:44	-															16				19					
*80:01:01:01-80:01:02, 80:03, 80:05-80:10	A80, Null, -						6																		23
*80:02	-						6						11												
*80:04	-						6												17						23
B*07:81, B*08:52, B*18:67, B*38:41, 38:67, B*49:67, B*51:01:110, 51:185, B*53:05:01-53:05:02, 53:16, 53:33, B*57:01:38, 57:60	B53, -									9															
B*07:439, B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569, 07:990, 07:994	-																								22
B*07:507, C*06:187	-																								23
B*15:01:83, B*35:01:60, B*58:01:24	-				w																				
B*18:27	-						6																		
C*12:131	-				4																				
Allele	Sero.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24