

INTERPRETATION TABLE																								
HLA-B low resolution SSP																								
Amplification patterns of the B*07:02 to B*83:01 alleles																								
		Well																						
		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	
Length of spec.		285	215	140	130	185	190	290	165	165	180	105	195	105	115	80	150	135	60	210	170	110	395	
PCR product(s)				235	265	235			220	190		195		395	150				245					
									330	390				435					400					
Length of int. pos. control ¹		800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070	1070	1070	1070	800	1070	1070	800	800	1070	1070	
5'-primer(s) ²		355	97	209	103	103	103	45	45	45	420	206	142	45	161	167	355	206	45	142	409	246	44	
		5'-TCA 3'	5'-TCg 3'	5'-ggC 3'	5'-CCg 3'	5'-CCT 3'	5'-CCg 3'	5'-ggA 3'	5'-ggA 3'	5'-ggA 3'	5'-TTA 3'	5'-AgA 3'	5'-TCT 3'	5'-ggA 3'	5'-Cgg 3'	5'-gCT 3'	5'-TCA 3'	5'-gAC 3'	5'-ggA 3'	5'-TCT 3'	5'-ggC 3'	5'-gAA 3'	5'-ggC 3'	
		363		363	103	363	418		357	435		420	419	540					368	368	420			
		5'-AgC 3'		5'-AgC 3'	5'-CCg 3'	5'-AAT 3'	5'-Agg 3'		5'-Tgg 3'	5'-AAA 3'		5'-TTA 3'	5'-gTC 3'	5'-gAC 3'					5'-gTT 3'	5'-gTC 3'	5'-TTA 3'			
				363	361				412										557					
				5'-AgC 3'	5'-AgT 3'				5'-ATA 3'										5'-ggA 3'					
3'-primer(s) ³		603	272	309	193	246	246	165	206	266	559	272	301	272	234	204	463	302	259	311	544	317	272	
		5'-gTg 3'	5'-Tgg 3'	5'-gTg 3'	5'-CgT 3'	5'-TAT 3'	5'-TAT 3'	5'-Tgg 3'	5'-CCT 3'	5'-TCC 3'	5'-CTC 3'	5'-Tgg 3'	5'-gTC 3'	5'-TgC 3'	5'-TCT 3'	5'-TCT 3'	5'-gCT 3'	5'-ggC 3'	5'-gTT 3'	5'-ggg 3'	5'-ggT 3'	5'-ggA 3'	5'-TgC 3'	
		605	272	559	193	246	572		538	559		572	301	309	272				259	311				
		5'-gCT 3'	5'-TgA 3'	5'-CAG 3'	5'-CgT 3'	5'-TAT 3'	5'-gCC 3'		5'-gTC 3'	5'-CAG 3'		5'-gCg 3'	5'-gTC 3'	5'-ATC 3'	5'-TgA 3'				5'-gTT 3'	5'-ggg 3'				
					583	559				583			570	605	272				266	538				
					5'-gTg 3'	5'-CTC 3'				5'-gTg 3'			5'-CCg 3'	5'-gCT 3'	5'-Tgg 3'				5'-TCC 3'	5'-gTC 3'				
													572						272					
													5'-gCg 3'						5'-Tgg 3'					
																			292					
																			5'-gTA 3'					
																			572					
																			5'-gCg 3'					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	

[illegible]

*08:04	B8, –	22	23					29			32												
*08:07, 08:14, 08:28, 08:37	B8, –	22	23																				
*08:09, 08:12:01-08:12:03, 08:16	B8		23					29															
*08:13, 08:20, 08:40	B8, –	22	23																				
*08:17, 08:54	–	22	23					29			32												
*08:21	B8	22	23	24				29															
*08:25	B8	22	23	24																			
*08:26, 08:50, 08:62, 42:07 ⁵	B8, –	22						29															
*08:32	B8	22	23					29								38							
*08:36	–	22	23					29															
*08:38	–	22	23					29															
*08:49	–	22	23											37			40						
*08:55	–	22	23															41					
*08:60	–	22	23					29															
*08:65	–	22						29															43
*13:01:01-13:01:05, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29	B13, –			24	25						31				37								
*13:02:01-13:02:09, 13:08Q, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40	B13, –			24	25						31												
*13:03	B49, B15			24	25																		
*13:04, 13:10, 13:21, 13:35	B49, B15, B13, –			24	25																		
*13:06	B13			24	25										37								
*13:07N	Null			24	25					31					37								
*13:09	–				25					31													
*13:12-13:13, 13:26	B13, –			24	25										37								
*13:14	B13			24						31													
*13:16	B13			24	25					31													
*13:31	–			24	25					31													
*13:36	–			24	25					31					37								
*13:39	–				25					31					37								
*14:01:01-14:01:02, 14:07N, 14:12, 14:14, 14:19	B64, –, Null				25	26																	43
*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18	B65, B14, –				25		27																43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40			

*14:02:03, 14:02:06	B65, –				25																		43
*14:02:05	B65				25		27																43
*14:05-14:06:02, 14:13	B14, –						27																43
*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50 ⁶	B39, B16, Null, –						26																43
*14:10	–						26									37							43
*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:190N, 15:192-15:193, 15:201,	B62, Null, –						28																
*15:01:02	B62						28		30								39						
*15:01:04	B62																39						
*15:01:17, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:02, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199	B62, –						28											39					
*15:01:19, 15:28, 15:120, 15:149N	B62, Null, –						28										39						
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

*15:02:01-15:02:05, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194	B75, –							28									37		39					
*15:03:01-15:03:02, 15:47, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173	B72, B70, blank, –							29		32									39					
*15:03:03	B72							28	29		32								39					
*15:04, 15:137	B62, –							28											39					
*15:07:01-15:07:02, 15:45, 15:126, 15:205, 15:207	B62, –			24				28											39					
*15:08, 15:11:01- 15:11:03, 15:11:05, 15:15, 15:148, 15:191, 15:209N	B75, B62, Null, –							28											39					
*15:09	B70								29	30			34						39					
*15:10:01, 15:18:01- 15:18:04, 15:37, 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133- 15:134, 15:153, 15:176, 15:197- 15:198, 15:200	B71, B70, –								29				34						39					
*15:10:02	B71							28	29				34						39					
*15:11:04	B75																		39					
*15:12, 15:14, 15:19	B76							28				33							39					
*15:13	B77							28					34			37			39					
*15:16, 15:67, 15:95	B63									30			34						39					
*15:17:01:01- 15:17:02, 15:162, 15:168, 15:177, 15:196, 15:208	B63, –							28		30			34						39					
*15:20, 15:25:01- 15:25:03, 15:204	B62, –							28								37			39					
*15:21, 15:44	B75							28					34			37			39					
*15:23, 15:115	B70, B5, –								29				34						39					
*15:24, 15:87, 15:157	B62, B15, –							28					34						39					
*15:29	B70								29										39					
*15:31	B75							28											39					
*15:36	–			w				28								37			39					
*15:42	B15							28											39					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
*15:43	B15							28											39				
*15:46	–				25			28	29										39				
*15:53	–				25				29										39				
*15:54	B72								29			32							39				
*15:55	B15			24				28											39				
*15:57	B62							28											39				
*15:62	–								29			32					37		39				
*15:64	B71								29			32							39				
*15:66	B62							28						34					39				
*15:68	–			24								32											
*15:69	B72								29			32							39		41		
*15:73	B62							28											39				
*15:74	B72								29										39				
*15:76	–							28											39				
*15:77	B62						27	28									37		39				
*15:80	B70								29					34			37		39				
*15:83	B62							28	29										39				
*15:85, 15:154, 15:165	B62, –							28									37		39				
*15:86	–							28											39				
*15:89	B77			w				28									37		39				
*15:91, 15:131	–								29			32	33						39				
*15:101	–							28															
*15:106	–				25			28	29								37		39				
*15:116	B62							28											39				
*15:124	B71								29					34					39				
*15:142	–		23					28											39				
*15:143	–							28	29										39				
*15:144	–							28									37		39				
*15:161	–								29			33	34						39				
*15:180	–		23						29										39				
*15:186	–								29					34					39		41		
*15:189	–						27	28						34					39				
*15:202	–							28										38	39				
*18:01:01-18:03, 18:05-18:08, 18:10- 18:11, 18:13-18:15, 18:17N-18:18, 18:20- 18:21, 18:23N-18:28, 18:30-18:36, 18:38- 18:43, 18:45-18:47, 18:50-18:52	B18, Null, –														35			38					
*18:04	B18														35			38					
*18:09	B18														35								
*18:12	B18														35			38					
*18:19	B18														35			38					
*18:22	B18														35		37	38					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

[illegible]

[illegible]

*39:01:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10-39:01:12, 39:03, 39:05:01-39:07, 39:09, 39:11, 39:14-39:15, 39:24, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51-39:57, 39:59-39:62	B 3901, B39, Low B39, B16, Null, –					26														41		43	
*39:01:09, 39:12, 39:45	B 3901, B39, –																			41		43	
*39:02:01, 39:08	B 3902, B39					26					32									41			
*39:02:02, 39:13:01-39:13:02, 39:23, 39:49	B 3902, B39, –										32									41			
*39:04	B39					27														41		43	
*39:10, 39:16, 39:20	B39																			41			
*39:17	B39																			41			
*39:18	B39					26														41		43	
*39:19:01	B39																38			41		43	
*39:19:02	B39					26											38			41		43	
*39:22	B39					26														41		43	
*39:36	–					26																43	
*39:39	B39					26					32									41			
*39:42	B39					26										37				41		43	
*39:58	–					26														41			
*40:01:01-40:01:06, 40:01:08-40:01:10, 40:01:12-40:01:17, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146-40:147, 40:150, 40:152-40:153, 40:155N	B60, B40, Null, –					25																	
*40:01:07, 40:53, 40:72:01-40:72:02, 40:77	B60, B61					25																	
*40:01:11	–					25																	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
*40:21	B40, B15																						
*40:25, 40:69, 40:106	B40, B60, –				25																		
*40:26	B21									30													
*40:28	–									30							37						
*40:30, 40:34	B60	22			25																		
*40:31, 40:45, 40:80	B60, –	22			25																		
*40:46	B60											32				36							
*40:47	B40				25																		
*40:48	B60				25						31												
*40:58	B40				25																		
*40:59	B40	22			25																		
*40:68	–																						
*40:71	B61			24																			
*40:74, 40:149	B60, –																						
*40:76	B60				25																		
*40:93	–											32				36							
*40:95	–																						
*40:96	–																						
*40:98	–	22																					
*40:109	–																						
*40:110	–				25																		
*40:117	–				25																		
*40:132	–				25																		
*40:135	–				25																		
*40:136	–				25				29					34									
*40:137	–	22																					
*40:138	–																						
*40:148	–																						
*40:157	–																						
*41:01, 41:03:01- 41:03:02, 41:06- 41:07, 41:12, 41:14- 41:16	B41, –				25				29														
*41:02:01-41:02:03, 41:11, 41:13	B41, –	22			25				29														
*41:04, 41:10	B41, –	22			25																		
*41:05	B41				25				29														
*41:08	B41				25				29					34									
*41:09, 45:02-45:03 ⁸	B45, –				25				29														
*42:01:01-42:02, 42:06, 42:10, 42:12	B42, –	22							29														
*42:04, 42:14	B42, –								29														
*42:05:01-42:05:02	B42	22							29					34									
*42:08	B42								29														
*42:09	B42	22																					
*42:11	–	22							29			32											
*42:13	–	22							29														
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

[illegible]

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
*48:04, 48:24	B48, -											32											
*48:05, 48:08, 48:10, 48:12, 48:14	B48, -	22										32											
*48:06	-	22																					
*48:13	-											32											
*48:16	B48	22										32											
*48:17	-	22										32											
*48:18	B48	22										32											
*48:21	-					26						32											
*49:01:01, 49:06, 49:08-49:10, 49:12-49:15	B49, -				25																		
*49:01:02, 49:11	B49, -				25																		
*49:02	B49				25																		
*49:03	B49																						
*49:04-49:05	B49				25																		
*49:07	-				25						31												
*50:01:01-50:01:02, 50:04-50:05, 50:08, 50:10-50:11	B50, -				25																		
*50:02	B45				25							32											
*50:06	-				25																		
*50:07	-																						
*50:09	-				25																		
*50:12	-				25																		
*51:01:01-51:01:03, 51:01:05-51:01:08, 51:01:10-51:02:03, 51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:26-51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:63, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:91, 51:94-51:96, 51:98N, 51:100, 51:102-51:103, 51:105	B51, B 5102, B 5103, Null, -									30													
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

[illegible]

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25- 57:30, 57:32-57:37	B57, - Null, -																37						
*57:01:05, 57:11, 57:13, 57:21, 57:31	B57, -																						
*57:04	B57						27										37						
*57:07, 57:24	B57, -																37						
*57:09	B57								29								37						
*57:12	B57																37						
*57:14	-																37						
*58:01:01-58:01:02, 58:01:04-58:01:08, 58:04-58:05, 58:10N- 58:15, 58:19, 58:21- 58:24, 58:28-58:31N	B58, - Null, -																37						
*58:01:03, 58:02, 58:06, 58:16, 58:25- 58:27	B58, -																						
*58:07	B58												33										
*58:08	B17, B5									30													
*58:09	B17																37						
*58:17N	Null																37						
*58:18	-			24																			
*58:20	-																			41			
*59:01, 59:05	B59, -																						
*59:02-59:03	B59, -																						
*59:04	-																						
*67:01:01, 67:03	B67, -					26														41			
*67:01:02-67:02	B67, -																			41			
*73:01-73:02	B73, -																						
*78:01, 78:02:02- 78:03, 78:07	B78, -									30								38					
*78:02:01, 78:04	B78, B35									30								38					
*78:05	-									30								38					
*78:06	-									30													
*81:01	B81	22																					
*81:02	B81	22																					
*81:03-81:05	Null, -	22																					
*82:01-82:03	B82, -								29			32						39					
*83:01	-																37						
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43
HLA-B allele ⁴	ser. ¹¹																						
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43

A*23:31, A*24:106, C*16:10						27																		
A*68:56, C*06:20, C*12:50																								
C*01:30				25																				
C*03:05, 03:25, 03:27			24																					
C*03:12, 03:19																								
C*03:102																37								
C*07:46								29																
C*15:02:04																								
C*15:25																								
C*15:39																								
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	

INTERPRETATION TABLE																											
HLA-B low resolution SSP																											
Amplification patterns of the B*07:02 to B*83:01 alleles																											
Well																											
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64							
160	180	105	325	115	130	90	90	145	120	430	145	300	160	180	180	90	95	115	360	350			Length of spec. PCR product(s)				
425				195	270	410	175	430	210				330	210		240		150									
				225																							
				260																							
1070	1070	1070	1070	1070	1070	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070	1070	1070			Length of int. pos. control ¹				
44	355	540	1 st I	355	209	41	209	48	357	49	206	1 st I	15	141	420	209	206	165	1 st I	1 st I			5'-primer(s) ²				
5' -ggC 3'	5' -TCA 3'	5' -gAC 3'	5' -CAG 3'	5' -TCA 3'	5' -ggC 3'	5' -CTg 3'	5' -ggC 3'	5' -gCC 3'	5' -Tgg 3'	5' -CAG 3'	5' -gAA 3'	5' -CAG 3'	5' -gCA 3'	5' -ATT 3'	5' -TTA 3'	5' -ggC 3'	5' -gAC 3'	5' -ACC 3'	5' -CAG 3'	5' -CAG 3'							
357	363			418	209	368	363	206					418	420		209		463									
5' -Tgg 3'	5' -Agg 3'			5' -Agg 3'	5' -ggg 3'	5' -gTT 3'	5' -AgC 3'	5' -gAA 3'					5' -Agg 3'	5' -TTA 3'		5' -ggA 3'		5' -TgA 3'									
				499	355											362											
				5' -TCT 3'	5' -TCA 3'											5' -ggT 3'											
302	499	603	282	538	299	282	259	309	435	309	311	259	175	311	559	256	259	272	317	311			3'-primer(s) ³				
5' -ggC 3'	5' -ggA 3'	5' -gTg 3'	5' -gCC 3'	5' -gTC 3'	5' -TCA 3'	5' -gCC 3'	5' -gTT 3'	5' -ATC 3'	5' -TCT 3'	5' -ATC 3'	5' -ggT 3'	5' -CTC 3'	5' -CCg 3'	5' -ggT 3'	5' -CAG 3'	5' -CCC 3'	5' -gTT 3'	5' -TgC 3'	5' -ggA 3'	5' -ggT 3'							
302		605		572	583	418	499		527			259	538	559		259		538									
5' -ggT 3'		5' -gCT 3'		5' -gCg 3'	5' -gTg 3'	5' -gTC 3'	5' -ggA 3'		5' -CCT 3'			5' -CTC 3'	5' -gTC 3'	5' -CgT 3'		5' -CTT 3'		5' -CCA 3'									
477												262				559											
5' -gCg 3'												5' -TgC 3'				5' -CgT 3'											
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			Well No.				

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																					ser. ¹¹	HLA-B allele ⁴
44			47																	64	B7, B42, Null, –	*07:02:01-07:02:20, 07:05:01-07:06, 07:09-07:15, 07:17- 07:18:02, 07:21- 07:23, 07:26, 07:29- 07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:48-07:49N, 07:51- 07:52, 07:55-07:59, 07:61-07:64, 07:66- 07:68:02, 07:70- 07:72, 07:74-07:76, 07:79-07:80, 07:82- 07:83, 07:86-07:99, 07:101-07:117
44																				64	B703, B7, –	*07:03, 07:08, 07:16, 07:32, 07:37
44		46	47																	64	B7	*07:04, 07:25
44			47																	64	B7, 42, –	*07:07, 07:19, 07:33, 07:40, 07:53, 07:77
44			47																	64	B7, blank	*07:20, 07:24
																			63	–	*07:27	
44			47																	64	B7	*07:28
			47																	63	B7, –	*07:36, 07:81
																				63	–	*07:38
44			w																	64	B7	*07:47
44																				64	–	*07:50
44			47									56								64	B7	*07:54
44			47																	64	–	*07:60
			47			50								58			w			64	–	*07:65
44																				64	–	*07:69, 07:85
44			47																	64	–	*07:73
44			47						53											64	–	*07:78
44			47						53											64	–	*07:84
44			47															62		64	–	*07:100
44		46																		64	B8, Null, –	*08:01:01-08:01:13, 08:05, 08:08N, 08:10 08:11, 08:15, 08:18- 08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33-08:35, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61, 08:63-08:64, 08:66
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
		46																	63		B8, –	*08:02-08:03, 08:52

44		46																64	B8, -	*08:04		
44		46																64	B8, -	*08:07, 08:14, 08:28, 08:37		
44		46																64	B8	*08:09, 08:12:01-08:12:03, 08:16		
44																		64	B8, -	*08:13, 08:20, 08:40		
44		46								56								64	-	*08:17, 08:54		
44		46																64	B8	*08:21		
44																		64	B8	*08:25		
44		46																64	B8, -	*08:26, 08:50, 08:62, 42:07 ⁵		
44		46													61			64	B8	*08:32		
		46																64	-	*08:36		
44		46								56								64	-	*08:38		
44		46		48														64	-	*08:49		
44																		64	-	*08:55		
44		46		48														64	-	*08:60		
44		46																64	-	*08:65		
										56								63	B13, -	*13:01:01-13:01:05, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29		
										56								63	B13, -	*13:02:01-13:02:09, 13:08Q, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40		
										56		59						63	B49, B15	*13:03		
										56								63	B49, B15, B13, -	*13:04, 13:10, 13:21, 13:35		
										56		58						63	B13	*13:06		
																		63	Null	*13:07N		
										56								64	-	*13:09		
										56								63	B13, -	*13:12-13:13, 13:26		
										56								63	B13	*13:14		
								53		56								63	B13	*13:16		
								w		56						62	63	-	*13:31			
					51					56								63	-	*13:36		
										56								64	-	*13:39		
44																		64	B64, -, Null	*14:01:01-14:01:02, 14:07N, 14:12, 14:14, 14:19		
44																		64	B65, B14, -	*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18		
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

44																	64	B65, –	*14:02:03, 14:02:06			
																	64	B65	*14:02:05			
44																	64	B14, –	*14:05-14:06:02, 14:13			
44																	64	B39, B16, Null, –	*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50 ⁶			
44				49													64	–	*14:10			
												56					62	64	B62, Null, – *15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:190N, 15:192-15:193, 15:201,			
											56						62	64	B62	*15:01:02		
											56						62	64	B62	*15:01:04		
												56						64	B62, –	*15:01:17, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:02, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199		
																	62	64	B62, Null, –	*15:01:19, 15:28, 15:120, 15:149N		
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

						51											64	B75, –	*15:02:01-15:02:05, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194			
											56						64	B72, B70, blank, –	*15:03:01-15:03:02, 15:47, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173			
											56						64	B72	*15:03:03			
							w				56					62	64	B62, –	*15:04, 15:137			
											56						62	64	B62, – 15:45, 15:126, 15:205, 15:207			
						51											62	64	B75, B62, Null, – 15:15, 15:148, 15:191, 15:209N			
																	64	B70	*15:09			
																	64	B71, B70, – 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133- 15:134, 15:153, 15:176, 15:197- 15:198, 15:200				
						51											64	B71	*15:10:02			
											56						62	64	B75 *15:11:04			
																	62	64	B76 *15:12, 15:14, 15:19			
						51											63	B77	*15:13			
							53										63	B63	*15:16, 15:67, 15:95			
																	63	B63, – 15:17:02, 15:162, 15:168, 15:177, 15:196, 15:208				
											56						64	B62, – 15:25:03, 15:204				
						51											64	B75 *15:21, 15:44				
																	63	B70, B5, – *15:23, 15:115				
											56						62	63	B62, B15, – *15:24, 15:87, 15:157			
																	64	B70 *15:29				
						51											64	B75 *15:31				
											56						63	– *15:36				
							53				56		58				64	B15 *15:42				
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
											56							62	63		B15	*15:43
											55	56						62		64	–	*15:46
											55	56						62		64	–	*15:53
												56						62		64	B72	*15:54
					51															64	B15	*15:55
			w									56						62		64	B62	*15:57
												56								64	–	*15:62
																				64	B71	*15:64
												56						62		64	B62	*15:66
												56						62		64	–	*15:68
												56								64	B72	*15:69
												56		59				62		64	B62	*15:73
												56								64	B72	*15:74
		47			w							56						62		64	–	*15:76
												56						62		64	B62	*15:77
																				64	B70	*15:80
												56								64	B62	*15:83
												56						62		64	B62, –	*15:85, 15:154, 15:165
												56		58						64	–	*15:86
					51														63		B77	*15:89
												56								64	–	*15:91, 15:131
		47										56						62		64	–	*15:101
											55	56								64	–	*15:106
45												56								64	B62	*15:116
45																				64	B71	*15:124
												56						62		64	–	*15:142
											55							62		64	–	*15:143
					51													62		64	–	*15:144
																				64	–	*15:161
																		62		64	–	*15:180
																				64	–	*15:186
					51													62		64	–	*15:189
												56						62		64	–	*15:202
																	61			64	B18, Null, –	*18:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:52
														58			61			64	B18	*18:04
																	61		63		B18	*18:09
												56								64	B18	*18:12
																	61	62		64	B18	*18:19
					49												61			64	B18	*18:22
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

44																		64	B 3901, B39, Low B39, B16, Null, –	*39:01:01:01- 39:01:01:02L, 39:01:03-39:01:08, 39:01:10-39:01:12, 39:03, 39:05:01- 39:07, 39:09, 39:11, 39:14-39:15, 39:24, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51- 39:57, 39:59-39:62		
44																		64	B 3901, B39, –	*39:01:09, 39:12, 39:45		
44										56								64	B 3902, B39	*39:02:01, 39:08		
44										56								64	B 3902, B39, –	*39:02:02, 39:13:01- 39:13:02, 39:23, 39:49		
44																		64	B39	*39:04		
44																		64	B39	*39:10, 39:16, 39:20		
44												58						64	B39	*39:17		
44																62		64	B39	*39:18		
44															61			64	B39	*39:19:01		
44															61			64	B39	*39:19:02		
44										56								64	B39	*39:22		
44																	62	64	–	*39:36		
44																		64	B39	*39:39		
44																		64	B39	*39:42		
44																		64	–	*39:58		
44	45										55	56						64	B60, B40, Null, –	*40:01:01-40:01:06, 40:01:08-40:01:10, 40:01:12-40:01:17, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146-40:147, 40:150, 40:152- 40:153, 40:155N		
44											55	56						64	B60, B61	*40:01:07, 40:53, 40:72:01-40:72:02, 40:77		
44	45										55	56		58				64	–	*40:01:11		
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

										55	56						64	B61, B41, Null, –	*40:02:01-40:04, 40:20, 40:39-40:40, 40:89, 40:94, 40:97, 40:122, 40:142N-40:144N		
										55	56						64	B4005, B50	*40:05		
44										55	56						64	B61, –	*40:06:01:01-40:06:03, 40:44, 40:70, 40:75, 40:83, 40:86, 40:103, 40:127, 40:131		
										55							64	B61, –	*40:08		
?										55	56						64	B61, –	*40:09, 40:11:01-40:11:02, 40:18, 40:24, 40:27, 40:29, 40:35, 40:37, 40:50, 40:56-40:57, 40:64, 40:78, 40:82, 40:85, 40:90-40:91, 40:99, 40:104-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:133Q, 40:145		
?	45									55	56						64	B60, Null, –	*40:10:01-40:10:02, 40:22N, 40:33, 40:36, 40:38, 40:42-40:43, 40:49, 40:52, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87-40:88, 40:92, 40:100-40:102, 40:108, 40:112-40:114, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:130, 40:134, 40:141, 40:151, 40:154, 40:156		
	45										56						64	B60, B48	*40:12		
						52		54		56							63	B61, –	*40:13, 40:19		
?										55	56						64	B60, –	*40:14:01-40:14:03, 40:121, 40:129, 40:158		
?										55	56						64	B40, –	*40:15, 40:32		
44										55	56						64	B61	*40:16		
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
?	45											56								64	B40, B15	*40:21
?	45										55									64	B40, B60, -	*40:25, 40:69, 40:106
?											55	56								64	B21	*40:26
?				49							55	56								64	-	*40:28
?	45										55	56								64	B60	*40:30, 40:34
?						51					55	56								64	B60, -	*40:31, 40:45, 40:80
?	45											56								64	B60	*40:46
	45											56							63		B40	*40:47
?	45										55	56								64	B60	*40:48
?											55	56		58						64	B40	*40:58
44	45										55	56								64	B40	*40:59
?											55									64	-	*40:68
?											55	56			59					64	B61	*40:71
?	45										55	56								64	B60, -	*40:74, 40:149
	45											56									B60	*40:76
44												56								64	-	*40:93
?								53			55	56								64	-	*40:95
44												56							63		-	*40:96
?											55	56								64	-	*40:98
44						52			54			56							63		-	*40:109
44												56							63		-	*40:110
	45					52			54			56							63		-	*40:117
?	45			48							55	56								64	-	*40:132
	45										55	56								64	-	*40:135
?		46									55	56								64	-	*40:136
?	45											56								64	-	*40:137
44	45										55	56								64	-	*40:138
44								53			55	56								64	-	*40:148
												56							63		-	*40:157
44		46									55	56								64	B41, -	*41:01, 41:03:01-41:03:02, 41:06-41:07, 41:12, 41:14-41:16
44		46									55	56								64	B41, -	*41:02:01-41:02:03, 41:11, 41:13
44		46									55	56								64	B41, -	*41:04, 41:10
44		46									?	56								64	B41	*41:05
44		46									55	56								64	B41	*41:08
44											55	56								64	B45, -	*41:09, 45:02-45:03 ⁸
44		46	47																	64	B42, -	*42:01:01-42:02, 42:06, 42:10, 42:12
44		46	47																	64	B42, -	*42:04, 42:14
44		46	47																	64	B42	*42:05:01-42:05:02
44			47																	64	B42	*42:08
44		46	47																	64	B42	*42:09
44		46										56								64	-	*42:11
44		46	47	48																64	-	*42:13
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
				48								56	57						63		B44, Null, –	*44:02:01:01-44:02:14, 44:11, 44:19N, 44:21-44:24, 44:27, 44:33-44:34, 44:48, 44:52N, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93, 44:101-44:102, 44:104, 44:112-44:113
				48								56							63		B44, Null, –	*44:03:01-44:03:06, 44:03:08-44:05:02, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35-44:36, 44:38-44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76-44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114-44:115
				48								56							63		B44, –	*44:03:07, 44:07, 44:81, 44:103, 44:111
				48						54			57				61		63		B44	*44:06
				48								56	57						63		B44	*44:08
				48							55	56	57						64		B45	*44:09
				48								56			59				63		B44	*44:10
				48									57						63		B44	*44:12
				48								56							63		B44	*44:14
												56			59				63		B12	*44:15
				48								56	57						63		B47, B44, –	*44:16, 44:91
				48								56							63		B44, –	*44:17, 44:43:01-44:43:02
							52					56			59				63		B44	*44:18
				48								56	57						63		B44, –	*44:20, 44:100
				48			52			54		56							63		B44	*44:25, 44:50
45				48								56							63		B44	*44:31
												56							63		B44, –	*44:37, 44:64:01-44:64:02
				48								56							63		B44	*44:40
				48								56	57						63		B44, –	*44:41, 44:53, 44:99

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.	
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.	
				48								56	57						63		B44	*44:44	
				48							55	56								64	B12, -	*44:46, 44:75	
				48								56							63		B44, -	*44:47, 44:65	
				48								56	57						63		B44	*44:49	
				48								56							63		B44	*44:54	
				48								56							63		-	*44:57	
				48								56							63		-	*44:60, 44:110	
				48								56							63		-	*44:62	
				48								56	57						63		-	*44:83	
			47	48		50					55	56	57							64	-	*44:90	
				48				52		54		56	57						63		-	*44:95	
			47	48		50						56	57						63		-	*44:97	
				48								56	57						63		-	*44:106	
44											55	56			59					64	B45, -	*45:01, 45:05, 45:07, 45:11-45:13	
44											55	56			59					64	B21	*45:04	
44			47								55	56			59					64	B45	*45:06	
44											55	56		58						64	-	*45:08	
44											55	56								64	-	*45:09	
44											55	56								64	-	*45:10	
			47		49														62		64	B46, Null, -	*46:01:01-46:01:05, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:22-46:24
			w		49														62		64	B46	*46:02
			47		49																64	B46, -	*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25-46:26
	45		47		49																64	-	*46:06
			47		49										59						64	B22	*46:11
			47		49														62		64	-	*46:12, 46:20
			47		49														62		64	-	*46:17
			47		49									58							64	-	*46:18
						50						56								63		B47, -	*47:01:01:01-47:01:01:02, 47:06-47:07
						50					55	56									64	B47	*47:02
						50					w	56								63	w	B47	*47:03
												56								63		B17, -	*47:04-47:05
44							51					56									64	B48, -	*48:01:01-48:01:03, 48:07, 48:09, 48:11, 48:15, 48:19-48:20, 48:22
44												56									64	B48, B70	*48:02:01-48:02:02
44	45											56									64	B48, -	*48:03:01-48:03:02, 48:23

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44							51					56								64	B48, -	*48:04, 48:24
44												56								64	B48, -	*48:05, 48:08, 48:10, 48:12, 48:14
44							51													64	-	*48:06
44												56								64	-	*48:13
44							51													64	B48	*48:16
44	45											56								64	-	*48:17
							51					56							63		B48	*48:18
44												56								64	-	*48:21
								52	53			56			59				63		B49, -	*49:01:01, 49:06, 49:08-49:10, 49:12-49:15
								52				56			59				63		B49, -	*49:01:02, 49:11
									53			56			59				63		B49	*49:02
									53			56			59				63		B49	*49:03
								52	53			56							63		B49	*49:04-49:05
								52	53			56							63		-	*49:07
44									53		55	56			59					64	B50, -	*50:01:01-50:01:02, 50:04-50:05, 50:08, 50:10-50:11
44									53		55	56			59					64	B45	*50:02
44									53			56			59					64	-	*50:06
44									53		55	56			59					64	-	*50:07
44									53		55	56								64	-	*50:09
44											55	56			59					64	-	*50:12
																					B51, B 5102, B 5103, Null, -	*51:01:01-51:01:03, 51:01:05-51:01:08, 51:01:10-51:02:03, 51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:26-51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:63, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:91, 51:94-51:96, 51:98N, 51:100, 51:102-51:103, 51:105
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

									53	54						61		63		B51, B5102, –	*51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92-51:93, 51:99
44									53	54						61		63		B51, –	*51:01:09, 51:24:03-51:24:04
				49						54						61		63		B51, –	*51:04, 51:46, 51:56
										54						61		63		B51, –	*51:06, 51:59
									53	54						61		63		B51, Null, –	*51:08, 51:20, 51:44N, 51:97
44										54						61		63		B51	*51:10
										54				59		61		63		–	*51:15, 51:62
44									53	54						61		63		B51	*51:24:02
									53	54						61		63		B51	*51:36
				48						54						61		63		–	*51:42
										54						61		63		–	*51:45, 51:81
									53							61		63		–	*51:54, 51:78
									w	54						61	62	63		–	*51:61
										54						61		63		–	*51:64
									53	54						61		63		–	*51:68
									53	54						61		63		–	*51:101
								51	53	54								63		–	*51:104
									53	54				59		61		63		–	*51:106
									53	54		56						63		B52, –	*52:01:01-52:01:04, 52:01:06-52:13, 52:15, 52:17-52:18, 52:22
									53	54		56						63		B52, –	*52:01:05, 52:14
									53	54		56						63		–	*52:16
									53	54		56						63		–	*52:19
									53			56						63		–	*52:20
									53	54		56					62	63		–	*52:21
										54						61		63		B53, –	*53:01:01-53:01:05, 53:04-53:05, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:24
				49						54						61		63		B53, B51	*53:02, 53:06
																61		63		B53, B35	*53:03, 53:09, 53:11-53:13
										54						61		63		–	*53:14
										54						61		63		–	*53:15
										54		56						63		–	*53:17
				48						54						61		63		–	*53:22
			47			50							57	58					64	B54, Null, –	*54:01-54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:23
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
			47			50							57		59					64	B22	*54:03
			47			50							57							64	B54, -	*54:04, 54:09, 54:14-54:15
			47			50							57					62		64	-	*54:06
44			47			50							57							64	B54	*54:11
			47			50		52					57	58					63		-	*54:12
			47			50			53				57	58						64	-	*54:20
			47			50			53					58						64	B55, -	*55:01:01-55:01:06, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45
			47			50								58						64	B55, B54, -	*55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 56:10 ⁹
			47			50														64	B55, B56, Null, -	*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 56:15, 56:18-56:19N, 56:22 ¹⁰
			47			50			53											64	B22, B55	*55:09, 55:24
			47			50										60				64	-	*55:14
			47			50			53											64	B55, -	*55:17, 55:28
												56		58						64	-	*55:18
	46	47				50														64	B22	*55:20
		47				50			w					58				62		64	B55	*55:21
		47				50														64	B22	*55:22
		47				50						56		58						64	-	*55:34
		47				50	w							58						64	-	*55:35
		47				50			53			56		58						64	-	*55:40
			47			50									59					64	B56, Null, -	*56:01:01-56:02, 56:04, 56:08, 56:14, 56:16, 56:20, 56:24, 56:26-56:30
			47			50												62		64	B22, 15	*56:03
			47			50			53					58						64	B56	*56:05:01
			47			50			53											64	B56	*56:05:02
			47			50			53					58			w			64	B78	*56:06
			47			50									59				63		B56	*56:07
			47			50														64	B56, B22	*56:09, 56:11-56:12
			47			50									59					64	B56	*56:13
															59					64	B56	*56:17
			47			50		52	53										63		-	*56:21
														58						64	-	*56:23
			47			50			53						59					64	-	*56:25

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
										54						60			63		B57, –	*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25- 57:30, 57:32-57:37
										54						60			63		B57, –	*57:01:05, 57:11, 57:13, 57:21, 57:31
										54						60			63		B57	*57:04
				48						54						60			63		B57, –	*57:07, 57:24
				48						54						60			63		B57	*57:09
																60				64	B57	*57:12
					49					54						60			63		–	*57:14
										54									63		B58, –	*58:01:01-58:01:02, 58:01:04-58:01:08, 58:04-58:05, 58:10N- 58:15, 58:19, 58:21- 58:24, 58:28-58:31N
										54									63		B58, –	*58:01:03, 58:02, 58:06, 58:16, 58:25- 58:27
										54									63		B58	*58:07
									53	54									63		B17, B5	*58:08
					49					54									63		B17	*58:09
										54									63		Null	*58:17N
										54									63		–	*58:18
										54									63		–	*58:20
								52						58					63		B59, –	*59:01, 59:05
								52											63		B59, –	*59:02-59:03
								52							59				63		–	*59:04
44			47																	64	B67, –	*67:01:01, 67:03
44			47																	64	B67, –	*67:01:02-67:02
44			47															62		64	B73, –	*73:01-73:02
									53					58			61			64	B78, –	*78:01, 78:02:02- 78:03, 78:07
									53								61			64	B78, B35	*78:02:01, 78:04
									53			56								64	–	*78:05
									53			56								64	–	*78:06
			47			50	51													64	B81	*81:01
44			47				51													64	B81	*81:02
			47				51													64	Null, –	*81:03-81:05
			47			50									59					64	B82, –	*82:01-82:03
			47	48		50							57							64	–	*83:01
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																					ser. ¹¹	HLA-B allele ⁴
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.

