

INTERPRETATION TABLE																					
HLA-B low resolution SSP																					
Amplification patterns of the B*07:02 to B*83:01 alleles																					
		Well																			
		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
Length of spec. PCR product(s)		285	215	140	130	185	190	290	165	165	180	105	195	105	115	80	150	135	60	210	170
				235	265	235			220	190		195		395	150				245		
									330	390				435					400		
Length of int. pos. control ¹		800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070	1070	1070	1070	800	1070	1070	800	800
5'-primer(s) ²		355	97	209	103	103	103	45	45	45	420	206	142	45	161	167	355	206	45	142	409
		5' -TCA 3'	5' -TCg 3'	5' -ggC 3'	5' -CCg 3'	5' -CCT 3'	5' -CCg 3'	5' -ggA 3'	5' -ggA 3'	5' -ggA 3'	5' -TTA 3'	5' -AgA 3'	5' -TCT 3'	5' -ggA 3'	5' -Cgg 3'	5' -gCT 3'	5' -TCA 3'	5' -gAC 3'	5' -ggA 3'	5' -TCT 3'	5' -ggC 3'
		363		363	103	103	418		357	435		420	419	540					368	368	420
		5' -AgC 3'		5' -AgC 3'	5' -CCg 3'	5' -CCT 3'	5' -Agg 3'		5' -Tgg 3'	5' -AAA 3'		5' -TTA 3'	5' -gTC 3'	5' -gAC 3'					5' -gTT 3'	5' -gTC 3'	5' -TTA 3'
				363	361	363			412										557		
				5' -AgC 3'	5' -AgT 3'	5' -AAT 3'			5' -ATA 3'										5' -ggA 3'		
3'-primer(s) ³		603	272	309	193	246	246	165	206	266	559	272	301	272	234	204	463	302	259	311	544
		5' -gTg 3'	5' -Tgg 3'	5' -gTg 3'	5' -CgT 3'	5' -TAT 3'	5' -TAT 3'	5' -Tgg 3'	5' -CCT 3'	5' -TCC 3'	5' -CTC 3'	5' -Tgg 3'	5' -gTC 3'	5' -TgC 3'	5' -TCT 3'	5' -TCT 3'	5' -gCT 3'	5' -ggC 3'	5' -gTT 3'	5' -ggg 3'	5' -ggT 3'
		605	272	559	193	246	572		538	559		572	301	309	272				259	311	
		5' -gCT 3'	5' -TgA 3'	5' -CAg 3'	5' -CgT 3'	5' -TAT 3'	5' -gCC 3'		5' -gTC 3'	5' -CAg 3'		5' -gCg 3'	5' -gTC 3'	5' -ATC 3'	5' -TgA 3'				5' -gTT 3'	5' -ggg 3'	
					583	559				583			570	605	272				266	538	
					5' -gTg 3'	5' -CTC 3'				5' -gTg 3'			5' -CCg 3'	5' -gCT 3'	5' -Tgg 3'				5' -TCC 3'	5' -gTC 3'	
													572						272		
													5' -gCg 3'						5' -Tgg 3'		
																			292		
																			5' -gTA 3'		
																			572		
																			5' -gCg 3'		
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41

[illegible]

[illegible]

[illegible]

[illegible]

*15:89	-			w				28								37		39			
*15:91, 15:131	-							29			32	33						39			
*15:101	-							28													
*15:106	-				25			28	29							37		39			
*15:116	-							28										39			
*15:124	-								29					34				39			
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*15:142	-		23					28										39			
*15:143	-							28	29									39			
*15:144	-							28									37	39			
*15:161	-								29			33	34					39			
*15:180	-		23																		
*15:186	-								29				34					39		41	
*15:189	-						27	28					34					39			
*15:202, 15:239	-							28										38	39		
*15:222	-									30				34				39			
*15:233	-						27	28										39			
*15:236	-							28													
*15:243	-								29	30				34				39			
*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71	B18, Null, -														35			38			
*18:04	B18														35			38			
*18:09	B18														35						
*18:12	-														35			38			
*18:19	-														35			38			
*18:22, 18:69	-														35		37	38			
*18:29, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151. 35:174 ⁸	B35, Null, -																	38			
*18:37	-										32				35						
*18:44	-				25										35			38			
*18:48	-														35						
*18:49	-														35						
*18:54	-											33			35				40		
*18:56	-											33			35			38			
*18:66	-														35						
*18:67	-														35						
*27:01	B27					26										36					
*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83	B27, Null, -					26										36					

*27:03-27:05:15, 27:05:17-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:41, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72-27:74, 27:78-27:80, 27:82-27:84, 27:88	B27, Null, -					26										36					
*27:05:16, 27:07:01-27:07:02, 27:11, 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70-27:71, 27:76	B27, -															36					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*27:08, 27:26, 27:40, 27:42, 27:44, 27:89	B2708, -					26										36					
*27:12	B27					26										36					
*27:14, 27:81	B27, -															36					
*27:16, 27:39, 27:59N	Null, -					26										36					
*27:18	-					26					32					36					
*27:19	B27															36	37				
*27:23	-					26										36					
*27:24	-															36					
*27:25	-															36					
*27:29	B27					26					32					36					
*27:30	B27															36	37				
*27:33	-															36					
*27:52, 27:67	-					26															
*27:62	-															36					
*27:75	-					26															
*27:77	-					26										36					
*35:01:01:01-35:01:09, 35:01:11-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06-35:08:04, 35:11:01-35:12:03, 35:15, 35:20:01-35:20:02, 35:23, 35:29:01-35:29:02, 35:33-35:34, 35:36, 35:39-35:42:01, 35:48, 35:50, 35:52, 35:54-35:55, 35:57, 35:59, 35:61, 35:65Q, 35:70, 35:74, 35:77-35:78, 35:82-35:85, 35:90-35:95, 35:98, 35:100-35:101:02, 35:103-35:108:02, 35:110-35:112, 35:116, 35:120-35:124, 35:126, 35:128-35:134N, 35:136-35:150, 35:152-35:153, 35:155-35:156, 35:158-35:162, 35:165N-35:168, 35:170-35:173N, 35:175-35:184 35:187 35:189	B35, Null, -																37	38			
*35:01:10, 35:04:02	B35									30							37	38			
*35:03:03, 35:27, 35:56	B35, -																37				

*35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199	B35, -	24																38							
*35:08:05	-													37											
*35:10, 35:13, 35:28, 35:69, 35:80	B35, -													37				38							
*35:14:01-35:14:02, 35:62, 35:163	B35, -													37				38							
*35:16	B35	24												37				38							
*35:17, 35:30, 35:113, 35:125, 35:164	B35, -	24												37				38							
*35:19	B35													37											
*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190	-													37				38							
*35:25, 35:49, 35:198	B35, -													37											
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41				
*35:26	-					27								37											
*35:35	B35													37				38				41			
*35:38, 35:115, 35:169	-													37				38							
*35:42:02	-													37				38							
*35:43-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135	B35, -																	38							
*35:45, 35:71	-									33				37				38							
*35:46	B35	25												37											
*35:47, 35:154	-	25												37											
*35:60	B35																	38							
*35:63	-	25												37											
*35:66	-	22																38							
*35:76	B35, 22													37				38							
*35:87	-	22				29												38							
*35:185	-																	38							
*35:186	-													37				38							
*35:191	-													37				38							
*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N	B37, Null, -									33				37				40							
*37:01:08, 37:09	B37, -									33								40							
*37:02	-					26				33								40							
*37:04:01-37:04:02	B37, -									33				37				40							
*37:05	-									33				37				40							
*37:07	-	22								33								40							
*37:08	B37									33				37				40							
*37:10, 37:22, 37:29	B37, -													37				40							
*37:11	-													37				38				40			
*37:12	-					29				33								40							
*37:14	-													37				40							
*37:28	-									32				33				37				40			
*38:01:01-38:02:02, 38:08-38:09, 38:11-38:16, 38:18-38:19, 38:21, 38:23-38:24, 38:27-38:29, 38:31-38:32, 38:34N. 38:38-38:40	B38(16), B16, Null, -					26																41			

[illegible]

[illegible]

[illegible]

*55:01:01-55:01:06, 55:01:08-55:01:09, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45. 55:52-55:55N	B55(22), - , B22, Null																				
*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 55:49, 56:15, 56:18-56:19N, 56:22, 56:31-56:32 ¹²	B55(22), B56(22), Null, -																				
*55:09, 55:24	B22, -									31											
*55:14	-																37				
*55:17, 55:28	-																				
*55:18	-										32										
*55:20	-								29												
*55:21	-																				
*55:22	B55(22)									31											
*55:34	-																				
*55:40	-																				
*55:51	-			24																	
*55:56	-							29													
*56:01:01-56:02, 56:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:35	B56(22), Null, -																				
*56:03	B22																				
*56:05:01	B56(22)									30											
*56:05:02	B56(22)									30											
*56:06	B78									30								38			
*56:07	B56(22)																				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41
*56:09, 56:11-56:12	B56(22), B55(22), -																37				
*56:13	B56(22)								29												
*56:17, 56:33	-																				
*56:20:02	-																				
*56:21	-									30											
*56:23	-																				
*56:25	-																				
*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50, 57:52-57:57. 58:36 ¹³	B57(17), Null, -																37				
*57:01:05, 57:11, 57:13, 57:21, 57:31, 57:47	B57(17), -																				
*57:04	B57(17)						27										37				
*57:07, 57:24	-																37				
*57:09	-								29								37				
*57:12	-																37				
*57:14	-																37				
*57:45, 57:51	-				25												37				

[illegible]

INTERPRETATION TABLE																											
HLA-B low resolution SSP																											
Amplification patterns of the B*07:02 to B*83:01 alleles																											
Well																											
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64					
110	395	160	180	105	325	115	130	90	90	145	120	430	145	300	160	180	180	90	95	115	360	350					
		425				195	270	410	175	430	210				330	210		240		150							
						225																					
						260																					
1070	1070	1070	1070	1070	800	1070	800	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070	1070	1070	1070				
246	44	44	355	540	1 st I	355	209	41	209	48	357	49	206	1 st I	15	141	420	209	206	165	1 st I	1 st I					
5'-gAA ^{3'}	5'-ggC ^{3'}	5'-ggC ^{3'}	5'-TCA ^{3'}	5'-gAC ^{3'}	5'-Cag ^{3'}	5'-TCA ^{3'}	5'-ggC ^{3'}	5'-CTg ^{3'}	5'-ggC ^{3'}	5'-gCC ^{3'}	5'-Tgg ^{3'}	5'-Cag ^{3'}	5'-gAA ^{3'}	5'-Cag ^{3'}	5'-gCA ^{3'}	5'-ATT ^{3'}	5'-TTA ^{3'}	5'-ggC ^{3'}	5'-gAC ^{3'}	5'-ACC ^{3'}	5'-Cag ^{3'}	5'-Cag ^{3'}					
		357	363			418	209	368	363	206					418	420		209		463							
		5'-Tgg ^{3'}	5'-Agg ^{3'}			5'-Agg ^{3'}	5'-ggg ^{3'}	5'-gTT ^{3'}	5'-AgC ^{3'}	5'-gAA ^{3'}					5'-Agg ^{3'}	5'-TTA ^{3'}		5'-ggA ^{3'}		5'-TgA ^{3'}							
						499	355											362									
						5'-TCT ^{3'}	5'-TCA ^{3'}											5'-ggT ^{3'}									
317	272	302	499	603	282	538	299	282	259	309	435	309	311	259	175	311	559	256	259	272	317	311					
5'-ggA ^{3'}	5'-TgC ^{3'}	5'-ggC ^{3'}	5'-ggA ^{3'}	5'-gTg ^{3'}	5'-gCC ^{3'}	5'-gTC ^{3'}	5'-TCA ^{3'}	5'-gCC ^{3'}	5'-gTT ^{3'}	5'-ATC ^{3'}	5'-TCT ^{3'}	5'-ATC ^{3'}	5'-ggT ^{3'}	5'-CTC ^{3'}	5'-CCg ^{3'}	5'-ggT ^{3'}	5'-Cag ^{3'}	5'-CCC ^{3'}	5'-gTT ^{3'}	5'-TgC ^{3'}	5'-ggA ^{3'}	5'-ggT ^{3'}					
		302		605		572	583	418	499		527			259	538	559		259		538							
		5'-ggT ^{3'}		5'-gCT ^{3'}		5'-gCg ^{3'}	5'-gTg ^{3'}	5'-gTC ^{3'}	5'-ggA ^{3'}		5'-CCT ^{3'}			5'-CTC ^{3'}	5'-gTC ^{3'}	5'-CgT ^{3'}		5'-CTT ^{3'}		5'-CCA ^{3'}							
		477												262				559									
		5'-gCg ^{3'}												5'-TgC ^{3'}				5'-CgT ^{3'}									
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64					

[illegible]

	43	44																				64	B14, -
	43	44																				64	B39(16), Null, -
	43	44					49															64	-
		44			47																	64	-
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
														56						62		64	B62(15), Null, -
														56						62		64	B62(15)
														56						62		64	B62(15)
														56								64	B62(15), B15, -
																				62		64	B62(15), Null, -
								51														64	B75(15), B15, -
														56								64	B72(70), B70, -
														56								64	B72(70)
											w			56						62		64	B62(15), -
														56						62		64	B62(15), -

									51									62		64	B75(15), B62(15), Null, -		
																			64		B70		
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
																					64	B71(70), B70, B15, Null, -	
																					64	B71(70)	
									51											62		64	B75(15)
														56						62		64	B76(15)
									51												63		B77(15), -
											53										63		B63(15), -
																					63		B63(15), -
														56								64	B62(15), -
									51													64	B75(15), -
																					63		-
														56						62	63		B62(15), -
									51													64	B15
														56								64	B75(15)
											53			56							63		-
														56		58						64	-
														56						62	63		-
													55	56						62		64	B72(70)
													55	56						62		64	-
														56						62		64	-
									51													64	B15
							w							56						62		64	-
														56								64	-
																						64	-
														56						62		64	-
														56						62		64	B35
														56								64	-
														56		59				62		64	B62(15)
														56								64	-
					47				w											62		64	-
														56						62		64	-
																						64	B70
														56								64	-
														56						62		64	-
														56		58						64	-

									51										63	-			
													56						64	-			
					47								56					62	64	-			
													55	56					64	-			
			45											56					64	-			
			45																64	-			
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
														56						62		64	-
													55							62		64	-
									51											62		64	-
																						64	-
																				62		64	-
																						64	-
									51											62		64	-
														56						62		64	-
																					63		-
														56						62		64	-
														56						62		64	-
														56								64	-
																			61			64	B18, Null, -
																58			61			64	B18
																			61		63		B18
														56								64	-
																			61	62		64	-
							49												61			64	-
																			61			64	B35, Null, -
																						64	-
																			61			64	-
														55								64	-
																			61			64	-
																			61		63		-
																			61			64	-
																						64	-
												54							61		63		-
					47			w						56						63			B27
					47			50				54		56						63			B27, Null, -

						47				50					56						63			B27, Null, -
						47				50					56						63			B27, -
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
					47			50						56								64	B2708, -	
														56								64	B27	
		44			47			50						56							63		B27, -	
														56							63		Null, -	
														56								64	-	
					47			50						56							63		B27	
																					63		-	
					47			50	51					56							63		-	
					47			50						56						62	63		-	
														56							63		B27	
					47			50				54		56							63		B27	
					47			50						56								64	-	
					47			50						56							63		-	
					47			50				54		56							63		-	
					47			50				54		56						62	63		-	
												54		56							63		-	
																			61			64	B35, Null, -	
																			61			64	B35	
																			61			64	B35, -	

																			61			64	B35, -	
																58						64	-	
																56						64	B35, -	
																			61	62		64	B35, -	
																56						64	B35	
																			61			64	B35, -	
																55						64	B35	
																					61		64	-
																							64	B35, -
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
																						64	-	
																			61			64	B35	
						48													61			64	-	
																58			61			64	-	
																			61	62		64	B35, -	
						48													61			64	-	
									51													64	B35	
														55								64	-	
																58			61			64	B35	
														55	56							64	-	
																			61			64	-	
						47			50										w			64	B35, 22	
				46															61			64	-	
																58			61	62		64	-	
						47			50													64	-	
																						64	-	
							48							56							63		B37, Null, -	
														56							63		B37, -	
														56							63		-	
							48	49						56							63		B37, -	
							48							56								64	-	
														56							63		-	
							48											61			63		B37	
							48							56							63		B37, -	
							48							56								64	-	
							48							56							63		-	
							48							56							63	w	-	
							48							56							63		-	
42	43																				63		B38(16), B16, Null, -	

	43																		63	-			
42													56						63	B16			
	43												56						63	-			
42	43																		63	B38(16), -			
																		61	63	-			
																		61	63	-			
42	43																		63	-			
42	43																		63	-			
42	43																		63	-			
42	43																		63	-			
42					47														63	-			
42	43																		63	-			
	43																		63	-			
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
	43	44																			64	B39(16), B3901, Low B39, Null, -	
	43	44																			64	B39(16), B3901, -	
		44											56								64	B39(16), B3902	
		44											56								64	B39(16), B3902, -	
	43	44																			64	B39(16)	
		44																			64	B39(16), -	
		44													58						64	-	
	43	44																	62		64	-	
	43	44																61			64	B39(16)	
	43	44																61			64	B39(16)	
	43	44											56								64	-	
	43																				64	-	
	43	44																	62		64	-	
		44																			64	-	
		44																			64	-	
		44																			64	-	
	43	44																			64	-	
		44																			64	-	
	43	44											55								64	-	
		44	45										55	56							64	B60(40), , Null	

						48							56	57						63		B44(12), Null, -	
						48								57						63		B44(12), -	
						48								56						63		B44(12), B21, Null, -	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
						48								56						63		B44(12), -	
						48						54			57				61	63		B44(12)	
						48								56	57					63		B44(12)	
						48							55	56	57						64	B45(12)	
						48								56			59			63		B44(12)	
						48								56						63		B44(12)	
														56			59			63		B12	
						48								56	57					63		B47, -	
						48								56						63		B44(12), -	
									52					56			59			63		-	
						48								56	57					63		-	
						48			52		54		56							63		-	
			45			48								56						63		B44(12)	
														56						63		-	
						48								56						63		B44(12), -	
						48								56	57					63		-, B44(12)	
						48								56	57					63		-	
						48							55	56							64	-	
						48								56						63		-	
						48								56	57					63		-	
						48								56						63		B44(12)	

		44												56							64	-
											52	53		56		59				63	B49(21), Null, -	
											52			56		59				63	B49(21), -	
												53		56		59				63	B49(21)	
												53		56		59				63	-	
											52	53		56						63	B49(21), -	
											52	53		56						63	-	
											52	53		56						63	-	
											52	53		56		59				63	-	
		44										53		55	56		59			64	B50(21), -	
		44										53		55	56		59			64	B45(12)	
		44										53			56		59			64	-	
		44										53		55	56		59			64	-	
		44										53		55	56					64	-	
		44												55	56		59			64	-	
		44										53		55	56					64	-	
		44										53		55			59			64	-	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64
												53	54						61	63	B51(5), B5102, B5103, Null, -	
											53	54						61	63	B51(5), B5102, -		
		44									53	54						61	63	B51(5), -		
							49					54						61	63	B51(5), -		
												54						61	63	B51(5), -		
											53	54						61	63	B51(5), Null, -		
		44										54						61	63	-		
												54				59		61	63	-		
		44									53	54						61	63	B51(5)		
											53	54						61	63	-		

							48							54						61		63		B44(12)
														54						61		63		-
													53							61		63		-
												w	54							61	62	63		-
													54							61		63		-
													53	54						61		63		-
													53	54						61		63		-
													53	54								63		B17, Null, -
													53	54				59		61		63		-
										52	53	54										63		-
												53	54		56							63		B52(5), -
42												53	54		56									B52(5), -
												53	54		56							63		-
												53	54		56							63		-
												53	54		56							63		-
												53			56							63		-
												53	54		56						62	63		-
												53	54		56		59					63		-
											53	54		56							63		-	
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
												54							61		63		B53, -	
							49					54							61		63		-	
																			61		63		B53, -	
												54							61		63		-	
												54							61		63		-	
												54		56							63		-	
							48					54							61		63		-	
					47			50							57	58						64	B54(22), Null, -	
					47			50								58						64	B55(22), B22, -	
					47			50							57	59						64	-	
					47			50							57							64	-	
					47			50							57				62			64	-	
		44			47			50							57							64	B54(22)	
					47			50		52					57	58					63		-	
					47			50			53				57	58						64	-	
					47			50							57							64	-	

						47				50			53														64	B55(22), , B22, Null
						47				50																	64	B55(22), B56(22), Null, -
						47				50			53														64	B22, -
						47				50								60									64	-
						47				50			53														64	-
															56		58										64	-
					46	47				50																	64	-
						47				50			w				58				62					64	-	
						47				50																	64	B55(22)
						47				50					56		58										64	-
						47				50			53		56		58										64	-
						47				50																	64	-
						47				50							58										64	-
																	58			w							64	-
																	59										64	-
						47				50							59										64	B56(22), Null, -
						47				50											62						64	B22
						47				50			53				58										64	B56(22)
						47				50			53														64	B56(22)
						47				50			53				58			w							64	B78
						47				50							59					63					64	B56(22)
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
						47				50																	64	B56(22), B55(22), -
						47				50							59										64	B56(22)
																	59										64	-
		44				47				50							59										64	-
						47				50		52	53												63			-
																	58										64	-
						47				50			53				59										64	-
												54						60								63		B57(17), Null, -
												54						60							63			B57(17), -
												54						60							63			B57(17)
							48					54						60							63			-
							48					54						60							63			-
																		60							64			-
								49				54						60							63			-
										52		54													63			-

[illegible]

[illegible]

Well No.
HLA-B allele ⁴
*07:02:01-07:02:31, 07:05:01-07:06, 07:09-07:15, 07:17-07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:48-07:49N, 07:51-07:52, 07:55-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:72, 07:74-07:76, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:138, 07:140-07:145, 07:147-07:148, 07:150, 07:154
*07:02:32
*07:03, 07:08, 07:16, 07:32, 07:37
*07:04, 07:25, 07:146
*07:07, 07:19, 07:33, 07:40, 07:53, 07:77, 07:122, 07:139
*07:20, 07:24, 07:131
*07:27
*07:28
*07:36, 07:81, 07:149
*07:38
*07:47
*07:50
*07:54, 07:123
*07:60
*07:65
*07:69, 07:85
*07:73
*07:78
*07:84
*07:100
*07:133
*07:134
*08:01:01-08:01:17, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33-08:35, 08:39, 08:41-08:48, 08:51, 08:53, 08:56-08:59, 08:61, 08:63-08:64, 08:66-08:69, 08:71-08:73, 08:75, 08:80-08:83, 08:86N
*08:02-08:03, 08:52, 08:78
*08:04
*08:07, 08:14, 08:28, 08:37, 08:74, 08:77
*08:09, 08:12:01-08:12:03, 08:16, 08:84
*08:13, 08:20, 08:40, 08:70, 08:79

*08:17, 08:54
Well No.
Well No.
*08:21
*08:25
*08:26, 08:50, 08:62, 08:85, 42:07 ⁵
*08:32
*08:36
*08:38
*08:49
*08:55
*08:60, 08:76
*08:65
*08:87
*08:88
*13:01:01-13:01:06, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29, 13:43, 13:50, 13:52
*13:02:01-13:02:12, 13:08Q, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40, 13:42, 13:44-13:45, 13:47, 13:49N, 13:54
*13:03
*13:04, 13:10, 13:21, 13:35, 44:135 ⁶
*13:06
*13:07N
*13:09
*13:12-13:13, 13:26, 13:51
*13:14
*13:16
*13:31
*13:36
*13:39
*13:41
*13:46
*13:48
*13:53
*14:01:01-14:01:02, 14:07N, 14:12, 14:14, 14:19, 14:26, 14:32
*14:02:01-14:02:02, 14:02:04, 14:03-14:04, 14:09, 14:11, 14:15-14:18, 14:20, 14:22- 14:23, 14:25, 14:27, 14:29, 14:31
*14:02:03, 14:02:06, 14:24, 14:28
*14:02:05, 14:30

*14:05-14:06:02, 14:13
*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50 ⁷
*14:10
*14:21
Well No.
Well No.
*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:28, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227, 15:228, 15:229
*15:01:02
*15:01:04
*15:01:17, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237
*15:01:19, 15:28, 15:120, 15:149N
*15:02:01-15:02:05, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213-15:214, 15:223
*15:03:01-15:03:02, 15:47:01-15:47:02, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242
*15:03:03
*15:04, 15:137
*15:07:01-15:07:02, 15:45, 15:126, 15:207

*15:08, 15:11:01-15:11:03, 15:11:05, 15:15, 15:148, 15:191, 15:209N, 15:215
*15:09
Well No.
Well No.
*15:10:01, 15:18:01-15:18:04, 15:37, 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133- 15:134, 15:153, 15:176, 15:197-15:198, 15:200, 15:221, 15:226N, 15:229, 15:238
*15:10:02
*15:11:04
*15:12, 15:14, 15:19
*15:13:01-15:13:02
*15:16:01-15:16:03, 15:67, 15:95
*15:17:01:01-15:17:02, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:230
*15:20, 15:25:01-15:25:03, 15:204, 15:240
*15:21, 15:44
*15:23, 15:115
*15:24, 15:87, 15:157
*15:29
*15:31
*15:36
*15:42
*15:43
*15:46
*15:53, 15:212
*15:54
*15:55
*15:57
*15:62
*15:64
*15:66
*15:68
*15:69
*15:73
*15:74
*15:76
*15:77
*15:80
*15:83
*15:85, 15:154, 15:165
*15:86, 15:224

*15:89
*15:91, 15:131
*15:101
*15:106
*15:116
*15:124
Well No.
Well No.
*15:142
*15:143
*15:144
*15:161
*15:180
*15:186
*15:189
*15:202, 15:239
*15:222
*15:233
*15:236
*15:243
*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71
*18:04
*18:09
*18:12
*18:19
*18:22, 18:69
*18:29, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151. 35:174 ⁸
*18:37
*18:44
*18:48
*18:49
*18:54
*18:56
*18:66
*18:67
*27:01
*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83

*27:03-27:05:15, 27:05:17-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:41, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72-27:74, 27:78-27:80, 27:82 27:84-27:88
*27:05:16, 27:07:01-27:07:02, 27:11, 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70-27:71, 27:76
Well No.
Well No.
*27:08, 27:26, 27:40, 27:42, 27:44, 27:89
*27:12
*27:14, 27:81
*27:16, 27:39, 27:59N
*27:18
*27:19
*27:23
*27:24
*27:25
*27:29
*27:30
*27:33
*27:52, 27:67
*27:62
*27:75
*27:77
*35:01:01:01-35:01:09, 35:01:11-35:03:02, 35:03:04-35:04:01, 35:04:03, 35:06-35:08:04, 35:11:01-35:12:03, 35:15, 35:20:01-35:20:02, 35:23, 35:29:01-35:29:02, 35:33-35:34, 35:36, 35:39-35:42:01, 35:48, 35:50, 35:52, 35:54-35:55, 35:57, 35:59, 35:61, 35:65Q, 35:70, 35:74, 35:77-35:78, 35:82-35:85, 35:90-35:95, 35:98, 35:100-35:101:02, 35:103-35:108:02, 35:110-35:112, 35:116, 35:120-35:124, 35:126, 35:128-35:134N, 35:136-35:150, 35:152-35:153, 35:155-35:156, 35:158-35:162, 35:165N-35:168, 35:170-35:173N, 35:175-35:184 35:187 35:189
*35:01:10, 35:04:02
*35:03:03, 35:27, 35:56

*35:05:01-35:05:02, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199
*35:08:05
*35:10, 35:13, 35:28, 35:69, 35:80
*35:14:01-35:14:02, 35:62, 35:163
*35:16
*35:17, 35:30, 35:113, 35:125, 35:164
*35:19
*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190
*35:25, 35:49, 35:198
Well No.
Well No.
*35:26
*35:35
*35:38, 35:115, 35:169
*35:42:02
*35:43-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135
*35:45, 35:71
*35:46
*35:47, 35:154
*35:60
*35:63
*35:66
*35:76
*35:87
*35:185
*35:186
*35:191
*37:01:01-37:01:07, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N
*37:01:08, 37:09
*37:02
*37:04:01-37:04:02
*37:05
*37:07
*37:08
*37:10, 37:22, 37:29
*37:11
*37:12
*37:14
*37:28
*38:01:01-38:02:02, 38:08- 38:09, 38:11-38:16, 38:18- 38:19, 38:21, 38:23-38:24, 38:27-38:29, 38:31-38:32, 38:34N. 38:38-38:40

*38:02:03
*38:03
*38:04, 38:25
*38:05, 38:33
*38:06
*38:07
*38:10, 38:36-38:37
*38:17
*38:20
*38:22
*38:26
*38:30
*38:35
Well No.
Well No.
*39:01:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10- 39:01:12, 39:03, 39:05:01- 39:07, 39:09, 39:11, 39:14- 39:15, 39:24:01, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40N-39:41, 39:44, 39:46, 39:48, 39:51-39:57, 39:59- 39:62, 39:64-39:68, 39:70- 39:71, 39:72
*39:01:09, 39:12, 39:45, 39:72
*39:02:01, 39:08
*39:02:02, 39:13:01-39:13:02, 39:23, 39:49
*39:04
*39:10:01-39:10:02, 39:16, 39:20
*39:17, 39:63
*39:18
*39:19:01
*39:19:02
*39:22
*39:24:02
*39:36
*39:39:01
*39:39:02
*39:42
*39:58
*39:69
*40:01:01-40:01:06, 40:01:08- 40:01:10, 40:01:12-40:01:17, 40:01:19-40:01:24, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146- 40:147, 40:150, 40:152- 40:153, 40:155N, 40:172,

*40:01:07, 40:01:18, 40:53, 40:72:01-40:72:02, 40:77
*40:01:11
*40:02:01-40:04, 40:06:05, 40:20, 40:39-40:40, 40:89, 40:94, 40:97, 40:122, 40:142N, 40:144N
*40:05
*40:06:01:01-40:06:04, 40:44, 40:70, 40:75, 40:83, 40:86, 40:103, 40:127, 40:131, 40:159, 40:162, 40:165, 40:167, 40:177, 40:190
*40:08
Well No.
Well No.
*40:09, 40:11:01-40:11:02, 40:18, 40:24, 40:27:01- 40:27:02, 40:29, 40:35, 40:37, 40:50, 40:56-40:57, 40:64, 40:78, 40:82, 40:85, 40:90- 40:91, 40:99, 40:104-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:133Q, 40:145, 40:169, 40:173, 40:176, 40:180-40:181.
*40:10:01-40:10:02, 40:135
*40:12
*40:13, 40:19
*40:14:01-40:14:03, 40:121, 40:129, 40:158, 40:164, 40:192
*40:15, 40:32
*40:16
*40:21
*40:22N, 40:33, 40:36, 40:38, 40:42-40:43, 40:49, 40:52, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100, 40:102, 40:108, 40:112- 40:114, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:130, 40:134, 40:141, 40:151, 40:154, 40:156, 40:163, 40:168, 40:171, 40:175, 40:179, 40:182-40:183, 40:187, 40:191, 40:193
*40:25, 40:69, 40:106
*40:26
*40:28
*40:30, 40:34, 40:160

*40:31, 40:45, 40:80
*40:46
*40:47, 40:188
*40:48
*40:58
*40:59
*40:68
*40:71
*40:74, 40:149, 40:170
*40:76
*40:93
*40:95
*40:96
*40:98
*40:109
*40:110
Well No.
Well No.
*40:117
*40:132, 40:178
*40:136
*40:137
*40:138
*40:148, 40:161
*40:157
*40:166
*40:174
*40:184
*40:185
*41:01, 41:03:01-41:03:02, 41:06-41:07, 41:12, 41:14- 41:17, 41:20
*41:02:01-41:02:04, 41:11, 41:13, 41:19
*41:04, 41:10, 41:18
*41:05
*41:08
*41:09, 45:02-45:03 ⁹
*42:01:01-42:02, 42:06, 42:10, 42:12, 42:15
*42:04, 42:14
*42:05:01-42:05:02
*42:08
*42:09, 42:16
*42:11
*42:13

*44:02:01-44:02:21, 44:02:23-44:02:24, 44:11, 44:19N, 44:21-44:24, 44:27:01- 44:27:02, 44:33-44:34, 44:48, 44:52N, 44:55, 44:59, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93, 44:101-44:102, 44:104, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126, 44:137, 44:142Q, 44:148-44:149N, 44:151-44:152
*44:02:22, 44:12
*44:03:01-44:03:06, 44:03:08- 44:05:03, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35- 44:36, 44:38-44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76- 44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114-44:115, 44:117, 44:120, 44:122-44:125, 44:128, 44:133, 44:136, 44:139, 44:141, 44:143
Well No.
Well No.
*44:03:07, 44:07, 44:81, 44:103, 44:111, 44:134
*44:06
*44:08
*44:09
*44:10
*44:14
*44:15
*44:16, 44:91, 44:132
*44:17, 44:43:01-44:43:02, 44:144
*44:18
*44:20, 44:100
*44:25, 44:50
*44:31
*44:37:01-44:37:02, 44:64:01- 44:64:02
*44:40, 44:130
*44:41:01-44:41:02, 44:53, 44:99, 44:138Q
*44:44
*44:46, 44:75, 44:129
*44:47, 44:65
*44:49
*44:54

*44:57
*44:60, 44:110
*44:62
*44:83, 44:145
*44:90
*44:95
*44:97
*44:106
*44:127
*44:131
*44:140
*44:146
*44:150
*45:01, 45:05, 45:07, 45:11-45:14
*45:04
*45:06
*45:08
*45:09
*45:10
*46:01:01-46:01:07, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:20, 46:22-46:24, 46:27-46:30
*46:02
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25-46:26
Well No.
Well No.
*46:06
*46:11
*46:12
*46:17
*46:18
*47:01:01:01-47:01:02, 47:06-47:08
*47:02
*47:03
*47:04-47:05
*48:01:01-48:01:03, 48:07, 48:09, 48:11, 48:15, 48:19-48:20, 48:22, 48:27
*48:02:01-48:02:02, 48:25
*48:03:01-48:03:02, 48:23
*48:04, 48:24
*48:05, 48:08, 48:10, 48:12, 48:14
*48:06
*48:13
*48:16
*48:17
*48:18

*48:21, 48:26
*49:01:01, 49:01:03, 49:06, 49:08-49:10, 49:12-49:17, 49:19N
*49:01:02, 49:11
*49:02
*49:03
*49:04-49:05
*49:07
*49:18
*49:20
*50:01:01-50:01:03, 50:04- 50:05, 50:08, 50:10-50:11, 50:13, 50:16
*50:02
*50:06
*50:07
*50:09
*50:12
*50:14
*50:15
Well No.
Well No.
*51:01:01-51:01:03, 51:01:05- 51:01:08, 51:01:10-51:01:20, 51:01:22-51:01:24, 51:01:26- 51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N- 51:12, 51:13:02-51:14, 51:16- 51:19, 51:21-51:23, 51:26- 51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57- 51:58, 51:60, 51:63, 51:65- 51:67, 51:69-51:77, 51:79- 51:80, 51:82-51:91, 51:94- 51:96, 51:98N, 51:100, 51:102- 51:103, 51:105, 51:107- 51:111, 51:113-51:117, 51:119, 51:121-51:123, 51:125- 51:127, 51:129-51:131
*51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92- 51:93, 51:99, 51:131
*51:01:09, 51:01:21, 51:01:25, 51:24:03-51:24:04
*51:04, 51:46, 51:56:01- 51:56:02
*51:06:01-51:06:02, 51:59
*51:08, 51:20, 51:44N, 51:97
*51:10
*51:15, 51:62
*51:24:02
*51:36

*51:42
*51:45, 51:81
*51:54, 51:78:01-51:78:02
*51:61
*51:64
*51:68
*51:101
*51:104, 51:118N, 58:08:01-58:08:02 ¹⁰
*51:106
*51:112
*52:01:01:01-52:01:04, 52:01:06-52:01:09, 52:01:11-52:13, 52:15, 52:17-52:18, 52:22-52:24, 52:27-52:29
*52:01:05, 52:14
*52:01:10
*52:16
*52:19
*52:20
*52:21
*52:25
*52:26
Well No.
Well No.
*53:01:01-53:01:07, 53:04-53:05, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:27
*53:02, 53:06
*53:03, 53:09, 53:11-53:13
*53:14
*53:15
*53:17:01-53:17:02
*53:22
*54:01:01, 54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:25
*54:01:02, 55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50, 56:10 ¹¹
*54:03
*54:04, 54:09, 54:14-54:15
*54:06
*54:11
*54:12
*54:20
*54:26

*55:01:01-55:01:06, 55:01:08-55:01:09, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45, 55:52-55:55N
*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 55:49, 56:15, 56:18-56:19N, 56:22, 56:31-56:32 ¹²
*55:09, 55:24
*55:14
*55:17, 55:28
*55:18
*55:20
*55:21
*55:22
*55:34
*55:40
*55:51
*55:56
*56:01:01-56:02, 56:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:35
*56:03
*56:05:01
*56:05:02
*56:06
*56:07
Well No.
Well No.
*56:09, 56:11-56:12
*56:13
*56:17, 56:33
*56:20:02
*56:21
*56:23
*56:25
*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50, 57:52-57:57, 58:36 ¹³
*57:01:05, 57:11, 57:13, 57:21, 57:31, 57:47
*57:04
*57:07, 57:24
*57:09
*57:12
*57:14
*57:45, 57:51

*58:01:01-58:01:02, 58:01:04-58:01:11, 58:04-58:05, 58:10N-58:15, 58:19, 58:21-58:24, 58:28-58:35
*58:01:03, 58:02, 58:06, 58:16, 58:25-58:27
*58:07
*58:09
*58:17N
*58:18
*58:20
*59:01:01:01-59:01:01:02, 59:05
*59:02-59:03
*59:04
*67:01:01, 67:03
*67:01:02-67:02
*73:01-73:02
*78:01:01-78:01:02, 78:02:02-78:03, 78:07
*78:02:01, 78:04
*78:05
*78:06
*81:01
*81:02
*81:03-81:05
*82:01-82:03
*83:01
A*23:31, A*24:106, C*07:231, C*16:10
A*24:174
Well No.
Well No.
A*26:68, A*68:56, C*02:56, C*06:20, C*12:50
C*01:30
C*03:05, C*03:25, C*03:27, C*03:143
C*03:12, C*03:19
C*03:102
C*03:129
C*06:72
C*07:02:30, C*08:16:02
C*07:46
C*15:02:04
C*15:25
C*15:39
C*15:51
Well No.