

INTERPRETATION TABLE																									
HLA-B low resolution SSP																									
Amplification patterns of the B*07:02 to B*83:01 alleles																									
		Well																							
		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		
Length of spec.		285	215	140	130	185	180	290	165	165	180	105	195	105	115	80	150	140	60	210	170	110	395		
PCR product(s)				235	265	235			220	390		195		435	150				400						
									330																
Length of int. pos. control ¹		800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070	1070	1070	1070	800	1070	1070	800	800	1070	1070		
5'-primer(s) ²		355	97	209	103	103	103	45	45	45	420	206	142	45	161	167	355	206	45	142	409	246	44		
		5'-TCA ³	5'-TCG ³	5'-ggC ³	5'-CCg ³	5'-CCT ³	5'-CCg ³	5'-ggA ³	5'-ggA ³	5'-ggA ³	5'-TTA ³	5'-AgA ³	5'-TCT ³	5'-ggA ³	5'-Cgg ³	5'-gCT ³	5'-TCA ³	5'-gAC ³	5'-ggA ³	5'-TCT ³	5'-ggC ³	5'-gAA ³	5'-ggC ³		
		363		363	103	103	418		357	435		420	419	540					368	368	420				
		5'-AgC ³		5'-AgC ³	5'-CCg ³	5'-CCT ³	5'-Agg ³		5'-Tgg ³	5'-AAA ³		5'-TTA ³	5'-gTC ³	5'-gAC ³					5'-gTT ³	5'-gTC ³	5'-TTA ³				
				363	361	363			412										557						
				5'-AgC ³	5'-AgT ³	5'-AAT ³			5'-ATA ³									5'-ggA ³							
3'-primer(s) ³		603	272	309	193	246	246	165	206	266	559	272	301	272	234	204	463	302	259	311	544	317	272		
		5'-gTg ³	5'-Tgg ³	5'-gTg ³	5'-CgT ³	5'-TAT ³	5'-TAT ³	5'-Tgg ³	5'-CCT ³	5'-TCC ³	5'-CTC ³	5'-Tgg ³	5'-gTC ³	5'-TgC ³	5'-TCT ³	5'-TCT ³	5'-gCT ³	5'-ggC ³	5'-gTT ³	5'-ggg ³	5'-ggT ³	5'-ggA ³	5'-TgC ³		
		605	272	312	193	246	572		538	559		572	301	309	272			312	259	311					
		5'-gCT ³	5'-TgA ³	5'-gCC ³	5'-CgT ³	5'-TAT ³	5'-gCC ³		5'-gTC ³	5'-CAG ³		5'-gCg ³	5'-gTC ³	5'-ATC ³	5'-TgA ³			5'-AgT ³	5'-gTT ³	5'-ggg ³					
				559	583	559			603	583			570	605	272				266	538					
				5'-CAG ³	5'-gTg ³	5'-CTC ³			5'-gTg ³	5'-gTg ³			5'-CCg ³	5'-gCT ³	5'-Tgg ³				5'-TCC ³	5'-gTC ³					
													572						272						
													5'-gCg ³						5'-Tgg ³						
																			292						
																			5'-gTA ³						
																		572							
																		5'-gCg ³							
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		

*14:02:03, 14:02:06, 14:24, 14:28	B65, -				25													43
*14:02:05, 14:02:08, 14:30	B65, -				25	27												43
*14:02:09, 14:05-14:06:02, 14:13	B14, -					27												43
*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50, 39:74, 39:82 ⁸	B39, B16, Null, -					26												43
*14:10	-					26							37					43
*14:21	-				25													
*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:30, 15:01:32, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227-15:228, 15:231-15:232, 15:234, 15:241, 15:244-15:245Q, 15:247, 15:251, 15:257-15:258N, 15:260-15:262N, 15:264, 15:267, 15:269, 15:272N	B62, Null, -					28								39				
*15:01:02	B62					28	30								39			
*15:01:04	B62														39			
*15:01:17, 15:01:31, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237, 15:248-15:249	B62, blank, -					28									39			
*15:01:19, 15:28, 15:120, 15:149N, 15:246N	B62, Null, -					28									39			
*15:02:01-15:02:06, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213-15:214, 15:223, 15:265	B75, -					28						37		39				

*15:83	B62					28	29								39			
*15:85, 15:154, 15:165	B62, -					28							37		39			
*15:86, 15:224	-					28									39			
*15:89	B77			w		28								37	39			
*15:91, 15:131	-						29			32	33				39			
*15:101	-					28												
*15:106	-			25		28	29							37	39			
*15:116	B62					28									39			
*15:124	B71						29					34			39			
*15:142	-		23			28									39			
*15:143	-					28	29								39			
*15:144	-					28								37	39			
*15:161	-						29				33	34			39			
*15:180	-		23															
*15:186	-						29					34			39	41		
*15:189	-					27	28					34			39			
*15:202, 15:239	-					28									38	39		
*15:222	-							30				34			39			
*15:233	-					27	28								39			
*15:236	-					28												
*15:243	-						29	30				34			39			
*15:255	-						29											
*15:259	-					28	29			32					39			
*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71, 18:73-18:78, 18:80-18:89	B18, blank, Null, -											35			38			
*18:04	B18											35			38			
*18:09	B18											35						
*18:12	B18											35			38			
*18:19	B18											35			38			
*18:22, 18:69	B18, -											35	37		38			
*18:29, 18:72, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151, 35:174, 35:205 ⁹	B35, Null, -														38			
*18:37	-									32		35						
*18:44:01-18:44:02	-			25								35			38			
*18:48	-											35						
*18:49, 18:79	-											35						
*18:54	-										33	35				40		
*18:56	-										33	35			38			
*18:66	-											35						
*18:67	-											35						
*27:01	B27				26								36					
*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83, 27:95, 27:102	B27, Null, -				26								36					

*35:01:10, 35:04:02	B35							30						37	38				
*35:03:03, 35:27, 35:56, 35:203	B35, -													37					
*35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199	B35, -			24											38				
*35:08:05	-													37	38				
*35:10, 35:13, 35:28, 35:69, 35:80, 35:226	B35, -													37	38				
*35:14:01-35:14:02, 35:62, 35:163, 35:213	B35, -													37	38				
*35:16	B35			24										37	38				
*35:17, 35:30, 35:113, 35:125, 35:164	B35, -			24										37	38				
*35:19	B35, blank													37					
*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190	B35, B78, -													37	38				
*35:25, 35:49, 35:198, 35:210	B35, -													37					
*35:26	B35					27								37					
*35:35	B35													37	38		41		
*35:38, 35:115, 35:169	B35, -													37	38				
*35:42:02	-													37	38				
*35:43:01-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135	B35, -														38				
*35:45, 35:71	B35									33				37	38				
*35:46	B35			25										37					
*35:47, 35:154, 35:217	B35, -			25										37					
*35:60	B35														38				
*35:63	-			25										37					
*35:66	B35	22													38				
*35:76	B35, 22													37	38				
*35:87	-	22					29								38				
*35:185	-														38				
*35:186	-													37	38				
*35:191	-													37	38				
*35:207	-			25										37					
*35:218	-		23											37	38				
*35:222	-													37	38				
*35:229	-													37	38				
*37:01:01-37:01:07, 37:01:09, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N, 37:35-37:36	B37, Null, -									33				37			40		
*37:01:08, 37:09	B37, -									33							40		
*37:02	B37, B27			26						33							40		
*37:04:01-37:04:02	B37, -									33				37			40		
*37:05	Blank									33				37			40		
*37:07	-	22								33							40		
*37:08	B37									33				37			40		
*37:10, 37:22, 37:29	B37, -													37			40		
*37:11, 37:37	-													37	38		40		
*37:12	B37					29				33							40		
*37:14	-													37			40		
*37:28	-									32	33			37			40		
*37:34	-													37			40		

[illegible]

[illegible]

[illegible]

*48:16	B48	22							32							
*48:17	-	22							32							
*48:18	B48	22							32							
*48:21, 48:26	-				26				32							
*48:28	-	22							32		35					
*49:01:01, 49:01:03-49:01:04, 49:06, 49:08-49:10, 49:12-49:17, 49:19N, 49:22-49:24	B49, Null, -			25												
*49:01:02, 49:11	B49, -			25												
*49:02	B49			25												
*49:03	B49															
*49:04-49:05	B49			25												
*49:07	-			25				31								
*49:18	-			25			30									
*49:20	-			25				32								
*49:21	-			25												
*49:25	-	23														
*49:26	-									36					42	43
*50:01:01-50:01:03, 50:04-50:05, 50:08, 50:10-50:11, 50:13, 50:16, 50:18-50:19, 50:32	B50, -			25												
*50:02	B45			25				32								
*50:06	-			25												
*50:07	-															
*50:09, 50:20	-			25												
*50:12	-			25												
*50:14	-			25			30									
*50:15	-			25												
*50:17	-			25							37					
*50:31	-			25				31								
*51:01:01-51:01:03, 51:01:05-51:01:08, 51:01:10-51:01:20, 51:01:22-51:01:24, 51:01:26-51:01:39, 51:01:41-51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:26-51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:63, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:91, 51:94-51:96, 51:98N, 51:100, 51:102-51:103, 51:105, 51:107-51:111, 51:113-51:117, 51:119-51:130, 51:132-51:134, 51:136-51:138, 51:140, 51:142, 51:145-51:146, 51:149N-	B51, B5102, B5103, Null, blank, -							30								
*51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92-51:93, 51:99, 51:131, 51:143-51:144	B51, B5102, -															
*51:01:09, 51:01:21, 51:01:25, 51:24:03-51:24:04	B51, -							30								
*51:01:40, 51:54, 51:78:01-51:78:02	-							30								
*51:04, 51:46, 51:56:01-51:56:02, 51:139	B51, -							30			37					

[illegible]

INTERPRETATION TABLE																											
HLA-B low resolution SSP																											
Amplification patterns of the B*07:02 to B*83:01 alleles																											
Well																											
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64							
160	180	105	325	80	130	90	90	145	120	430	145	300	160	180	180	90	95	115	360	350					Length of spec.		
425				115	270	410	175	430	210				330	210		240		150							PCR product(s)		
				195																							
				260																							
1070	1070	1070	800	1070	800	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070	1070	1070					Length of int.		
44	355	540	1 st	355	209	41	209	48	357	49	206	1 st	15	141	420	209	206	165	1 st	1 st					pos. control ¹		
																									5'-primer(s) ²		
5'-ggC ^{3'}	5'-TCA ^{3'}	5'-gAC ^{3'}	5'-CAG ^{3'}	5'-TCA ^{3'}	5'-ggC ^{3'}	5'-CTg ^{3'}	5'-ggC ^{3'}	5'-gCC ^{3'}	5'-Tgg ^{3'}	5'-CAG ^{3'}	5'-gAA ^{3'}	5'-CAG ^{3'}	5'-gCA ^{3'}	5'-ATT ^{3'}	5'-TTA ^{3'}	5'-ggC ^{3'}	5'-gAC ^{3'}	5'-ACC ^{3'}	5'-CAG ^{3'}	5'-CAG ^{3'}							
357	363			418	209	368	363	206					418	420		209		463									
5'-Tgg ^{3'}	5'-Agg ^{3'}			5'-Agg ^{3'}	5'-ggg ^{3'}	5'-gTT ^{3'}	5'-AgC ^{3'}	5'-gAA ^{3'}	5'-Tgg ^{3'}	5'-CAG ^{3'}	5'-gAA ^{3'}	5'-CAG ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}		5'-ggA ^{3'}		5'-TgA ^{3'}									
				499	355											362											
				5'-TCT ^{3'}	5'-TCA ^{3'}											5'-ggT ^{3'}											
302	499	603	282	538	299	282	259	309	435	309	311	259	175	311	559	256	259	272	317	311					3'-primer(s) ³		
5'-ggC ^{3'}	5'-ggA ^{3'}	5'-gTg ^{3'}	5'-gCC ^{3'}	5'-gTC ^{3'}	5'-TCA ^{3'}	5'-gCC ^{3'}	5'-gTT ^{3'}	5'-ATC ^{3'}	5'-TCT ^{3'}	5'-ATC ^{3'}	5'-ggT ^{3'}	5'-CTC ^{3'}	5'-CCg ^{3'}	5'-ggT ^{3'}	5'-CAG ^{3'}	5'-CCC ^{3'}	5'-gTT ^{3'}	5'-TgC ^{3'}	5'-ggA ^{3'}	5'-ggT ^{3'}							
302		605		572	583	418	499		527			259	538	559		259		538									
5'-ggT ^{3'}		5'-gCT ^{3'}		5'-gCg ^{3'}	5'-gTg ^{3'}	5'-gTC ^{3'}	5'-ggA ^{3'}	5'-CCT ^{3'}				5'-CTC ^{3'}	5'-gTC ^{3'}	5'-CgT ^{3'}		5'-CTT ^{3'}		5'-CCA ^{3'}									
477												262				263											
5'-gCg ^{3'}												5'-TgC ^{3'}				5'-gTT ^{3'}											
																559											
																5'-CgT ^{3'}											
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64					Well No.		

44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																					ser. ¹⁷	HLA-B allele ⁴
44			47																	64	B7, B42, Null, blank, -	*07:02:01-07:02:31, 07:02:33-07:02:38, 07:05:01-07:06, 07:09-07:15, 07:17-07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:48-07:49N, 07:51-07:52, 07:55-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:72, 07:74-07:76, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:138, 07:140-07:145, 07:147-07:148, 07:150-07:164, 07:166-07:173, 07:175-07:178, 07:181N-07:190
44			47																	64	-	*07:02:32
44																				64	B703, B7, blank	*07:03, 07:08, 07:16, 07:32, 07:37
44		46	47																	64	B7, -	*07:04, 07:25, 07:146
44			47																	64	B7, 42, -	*07:07, 07:19, 07:33:01, 07:40, 07:53, 07:77, 07:122, 07:139, 07:165
44			47																	64	B7, blank, -	*07:20, 07:24, 07:131
																			63	-	*07:27	
44			47																	64	B7	*07:28
			47																	63	B7, -	*07:36, 07:81, 07:149
																				63	-	*07:38
44			w																	64	B7	*07:47
44																				64	-	*07:50
44			47									56								64	B7, -	*07:54, 07:123
44			47																	64	-	*07:60
			47			50								58			w			64	-	*07:65
44																				64	-	*07:69, 07:85
44			47																	64	-	*07:73
44			47						53											64	-	*07:78
44			47						53											64	-	*07:84
44			47															62		64	-	*07:100
44			47								55									64	-	*07:133
44			47														w			64	-	*07:134
			47			51														64	B81, -	*07:174, 81:02 ⁵
																				64	-	*07:179
																			63	-	*07:180	
44		46																		64	B8, Null, -	*08:01:01-08:01:20, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:22-08:24, 08:27, 08:29-08:31, 08:33-08:35, 08:39, 08:41-08:48, 08:51, 08:53:01-08:53:02, 08:56:01-08:59:02, 08:61, 08:63-08:64, 08:66-08:69, 08:71-08:73, 08:75, 08:80-08:83, 08:86N, 08:90-08:93, 08:05-08:100, 08:104, 08:106
		46																		63	B8, -	*08:02-08:03, 08:52, 08:78
44		46																		64	B8, blank	*08:04

44	46									64	B8, -	*08:07, 08:14, 08:28, 08:37, 08:74, 08:77, 08:108
44	46									64	B8, -	*08:09, 08:12:01-08:12:03, 08:16, 08:84
44										64	B8, -	*08:13, 08:20, 08:40, 08:70
44	46					56				64	-	*08:17, 08:54
44	46									64	B8	*08:21
44										64	B8	*08:25
44	46									64	B8, -	*08:26, 08:50, 08:62, 08:85, 08:94, 42:07 ⁶
44	46							61		64	B8	*08:32
	46									64	-	*08:36
44	46					56				64	-	*08:38, 08:101
44	46	48								64	-	*08:49
44										64	-	*08:55
44	46	48								64	-	*08:60, 08:76
44	46									64	-	*08:65
44										64	-	*08:79
44	46	48				57				64	-	*08:87
44	46									64	-	*08:88, 08:102
44	46									64	-	*08:89, 08:107
44										64	-	*08:103
						56				63	B13, Null, -	*13:01:01-13:01:06, 13:17, 13:20, 13:22:01-13:23, 13:25, 13:28-13:29, 13:43, 13:50, 13:52, 13:57, 13:60-13:61, 13:63N
						56				63	B13, Null, -	*13:02:01-13:02:15, 13:08, 13:11, 13:15, 13:18-13:19, 13:27, 13:30, 13:32-13:34, 13:37-13:38, 13:40, 13:42, 13:44-13:45, 13:47, 13:49N, 13:54-13:56N. 13:64-13:66
						56	59			63	B49, B15	*13:03
						56				63	B49, B15, B13, -	*13:04, 13:10, 13:21, 13:35, 13:59, 44:135, 44:158 ⁷
						56	58			63	B13	*13:06
										63	Null	*13:07N
						56				64	-	*13:09, 13:67
						56				63	B13, -	*13:12:01-13:13, 13:26, 13:51
						56				63	B13	*13:14
				53		56				63	B13	*13:16
				w		56			62	63	-	*13:31
		51								63	-	*13:36
						56				64	-	*13:39
						56			62	63	-	*13:41
						56				63	-	*13:46
			53			56	59			63	-	*13:48
						56	58			63	-	*13:53
						56				63	-	*13:58
				53		56				63	-	*13:62
44										64	B64, Null, -	*14:01:01-14:01:03, 14:07N, 14:12, 14:14, 14:19, 14:26, 14:32
44										64	B65, B14, -	*14:02:01-14:02:02, 14:02:04, 14:02:07, 14:03-14:04, 14:09, 14:11, 14:15-14:18, 14:20, 14:22-14:23, 14:25, 14:27, 14:29, 14:31, 14:33

44																64	B65, -	*14:02:03, 14:02:06, 14:24, 14:28	
																64	B65, -	*14:02:05, 14:02:08, 14:30	
44																64	B14, -	*14:02:09, 14:05-14:06:02, 14:13	
44																64	B39, B16, Null, -	*14:08, 39:25N, 39:30, 39:32-39:34, 39:43, 39:47, 39:50, 39:74. 39:82 ⁸	
44				49												64	-	*14:10	
44			47													64	-	*14:21	
																	B62, Null, -	*15:01:01:01-15:01:01:02N, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:30, 15:01:32, 15:26N-15:27:03, 15:30, 15:32, 15:34-15:35, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:152, 15:159-15:160, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227-15:228, 15:231-15:232, 15:234, 15:241, 15:244-15:245Q, 15:247, 15:251, 15:257-15:258N, 15:260-15:262N, 15:264, 15:267, 15:269. 15:272N	
								56								62	64	B62	*15:01:02
								56								62	64	B62	*15:01:04
																	64	B62, blank, -	*15:01:17, 15:01:31, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:48, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237, 15:248-15:249
																62	64	B62, Null, -	*15:01:19, 15:28, 15:120, 15:149N, 15:246N
						51											64	B75, -	*15:02:01-15:02:06, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213-15:214. 15:223. 15:265

[illegible]

					56					64	B62	*15:83
					56				62	64	B62, -	*15:85, 15:154, 15:165
					56	58				64	-	*15:86, 15:224
			51							63	B77	*15:89
					56					64	-	*15:91, 15:131
	47				56				62	64	-	*15:101
					55 56					64	-	*15:106
45					56					64	B62	*15:116
45										64	B71	*15:124
					56				62	64	-	*15:142
					55				62	64	-	*15:143
			51						62	64	-	*15:144
										64	-	*15:161
									62	64	-	*15:180
										64	-	*15:186
			51						62	64	-	*15:189
					56				62	64	-	*15:202, 15:239
										63	-	*15:222
					56				62	64	-	*15:233
					56				62	64	-	*15:236
					56					64	-	*15:243
	47				56					64	-	*15:255
					56				62	64	-	*15:259
								61		64	B18, blank, Null, -	*18:01:01:01-18:03, 18:05-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:18, 18:20-18:21, 18:23N-18:28, 18:30-18:36, 18:38-18:43, 18:45-18:47, 18:50-18:53, 18:55, 18:57-18:65, 18:68, 18:70-18:71, 18:73-18:78, 18:80-18:89
						58		61		64	B18	*18:04
								61		63	B18	*18:09
					56					64	B18	*18:12
								61	62	64	B18	*18:19
		49						61		64	B18, -	*18:22, 18:69
								61		64	B35, Null, -	*18:29, 18:72, 35:09:01-35:09:03, 35:18, 35:31-35:32:02, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:75, 35:88, 35:99, 35:118-35:119, 35:127, 35:151, 35:174, 35:205 ⁹
										64	-	*18:37
								61		64	-	*18:44:01-18:44:02
					55					64	-	*18:48
								61		64	-	*18:49, 18:79
								61	63		-	*18:54
								61		64	-	*18:56
										64	-	*18:66
				54				61	63		-	*18:67
	47		w		56				63		B27	*27:01
	47		50		54 56					63	B27, Null, -	*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83, 27:95, 27:102

																			B27, Null, -	*27:03-27:05:15, 27:05:17- 27:05:19, 27:05:21-27:06: 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35- 27:37, 27:41, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72- 27:74, 27:78-27:80, 27:82, 27:84-27:88, 27:90:01-27:91, 27:93-27:94N, 27:96-27:101, 27:103-27:105
																			B27, -	*27:05:16, 27:05:20, 27:07:01- 27:07:03, 27:11, 27:15, 27:21, 27:28, 27:32, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70- 27:71. 27:76
																		64	B2708, B27, B7, -	*27:08, 27:26, 27:40, 27:42, 27:44, 27:89
																		64	-	*27:12
44																		63	B27, -	*27:14, 27:81
																		63	B27, Null, -	*27:16, 27:39, 27:59N
																		64	-	*27:18
																		63	B27	*27:19
																		63	-	*27:23, 27:92
																		63	B27	*27:24
																		62	B27	*27:25
																		63	B27	*27:29
																		63	B27	*27:30
																		64	-	*27:33
																		63	-	*27:52, 27:67
																		63	-	*27:62
																		62	63	- *27:75
																		63	-	*27:77
																			B35, Null, -	*35:01:01:01-35:01:09, 35:01:11-35:03:02, 35:03:04- 35:04:01, 35:04:03, 35:06- 35:08:04, 35:08:06, 35:11:01- 35:12:03, 35:15, 35:20:01- 35:20:02, 35:23, 35:29:01- 35:29:02, 35:33-35:34, 35:36, 35:39-35:42:01, 35:48, 35:50, 35:52, 35:54-35:55, 35:57, 35:59, 35:61:01-35:61:02, 35:65Q, 35:70, 35:74, 35:77- 35:78, 35:82-35:85, 35:90- 35:95, 35:98, 35:100- 35:101:02, 35:103-35:108:02, 35:110-35:112, 35:116, 35:120- 35:124, 35:126, 35:128- 35:134N, 35:136-35:150, 35:152-35:153, 35:155-35:156, 35:158-35:162, 35:165N- 35:168, 35:170-35:173N, 35:175-35:184, 35:187, 35:189, 35:192-35:197, 35:200- 35:202, 35:204, 35:206, 35:208-35:209, 35:211-35:212, 35:214-35:216N, 35:219- 35:221, 35:223-35:225, 35:227- 35:228

									61			64	B35	*35:01:10, 35:04:02
									61			64	B35, -	*35:03:03, 35:27, 35:56, 35:203
									61			64	B35, -	*35:05:01-35:05:03, 35:22, 35:51, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199
								58				64	-	*35:08:05
							56					64	B35, -	*35:10, 35:13, 35:28, 35:69, 35:80, 35:226
									61	62		64	B35, -	*35:14:01-35:14:02, 35:62, 35:163, 35:213
							56					64	B35	*35:16
									61			64	B35, -	*35:17, 35:30, 35:113, 35:125, 35:164
							55					64	B35, blank	*35:19
			49						61			64	B35, B78, -	*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190
												64	B35, -	*35:25, 35:49, 35:198, 35:210
												64	B35	*35:26
									61			64	B35	*35:35
			48						61			64	B35, -	*35:38, 35:115, 35:169
								58	61			64	-	*35:42:02
									61	62		64	B35, -	*35:43:01-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135
			48						61			64	B35	*35:45, 35:71
				51								64	B35	*35:46
							55					64	B35, -	*35:47, 35:154, 35:217
								58	61			64	B35	*35:60
							55	56				64	-	*35:63
									61			64	B35	*35:66
			47		50				w			64	B35, 22	*35:76
		46							61			64	-	*35:87
								58	61	62		64	-	*35:185
			47		50							64	-	*35:186
												64	-	*35:191
							56					64	-	*35:207
									61			64	-	*35:218
									61			-	*35:222	
			47		50				61			64	-	*35:229
			48				56				63	B37, Null, -	*37:01:01-37:01:07, 37:01:09, 37:03N, 37:06, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N, 37:35-37:36	
							56				63	B37, -	*37:01:08, 37:09	
							56				63	B37, B27	*37:02	
			48	49			56				63	B37, -	*37:04:01-37:04:02	
			48				56				64	Blank	*37:05	
							56				63	-	*37:07	
			48						61		63	B37	*37:08	
			48				56				63	B37, -	*37:10, 37:22, 37:29	
			48				56				64	-	*37:11, 37:37	
			48				56				63	B37	*37:12	
			48				56				63	w	-	*37:14
			48				56				63	-	*37:28	
			48			54	56				63	-	*37:34	

										63		B38, Null, -	*38:01:01-38:02:02, 38:02:04-38:02:05, 38:08-38:09, 38:11-38:16, 38:18-38:19, 38:21, 38:23-38:24, 38:27-38:29, 38:31-38:32, 38:34N, 38:38-38:43
										63		-	*38:02:03
						56				63		B38	*38:03
						56				63		B38, -	*38:04, 38:25
										63		B38, -	*38:05, 38:33
									61	63		B38	*38:06
									61	63		B38	*38:07
										63		B38, -	*38:10, 38:36-38:37
										63		-	*38:17
										63		-	*38:20
										63		-	*38:22
		47								63		-	*38:26
										63		-	*38:30
										63		-	*38:35
44											64	B3901, B39, Low B39, B16, Null, -	*39:01:01:01-39:01:01:03, 39:01:03-39:01:08, 39:01:10-39:01:14, 39:01:16, 39:03, 39:05:01-39:07, 39:09, 39:11, 39:14-39:15, 39:24:01, 39:26-39:29, 39:31, 39:35, 39:37-39:38Q, 39:40:01N-39:41, 39:44, 39:46, 39:48, 39:51-39:57, 39:59-39:62, 39:64-39:68, 39:70-39:71, 39:73, 39:75, 39:77, 39:79, 39:81
44											64	B3901, B39, -	*39:01:09, 39:12, 39:45, 39:72, 39:76, 39:78
44												-	*39:01:15
44						56					64	B3902, B39	*39:02:01, 39:08
44						56					64	B3902, B39, -	*39:02:02, 39:13:01-39:13:02, 39:23, 39:49
44											64	B39	*39:04
44											64	B39, -	*39:10:01-39:10:02, 39:16, 39:20
44								58			64	B39, -	*39:17, 39:63
44										62	64	B39	*39:18
44									61		64	B39	*39:19:01
44									61		64	B39	*39:19:02
44						56					64	B39	*39:22
											64	-	*39:24:02
44										62	64	-	*39:36
44											64	B39	*39:39:01
44											64	-	*39:39:02
44											64	B39	*39:42
44											64	-	*39:58
44						55					64	-	*39:69
44	45					55	56				64	B60, B40, Null, -	*40:01:01-40:01:06, 40:01:08-40:01:10, 40:01:12-40:01:17, 40:01:19-40:01:29, 40:07, 40:23, 40:51, 40:60, 40:66, 40:73, 40:79, 40:124:01, 40:139-40:140, 40:146-40:147, 40:150, 40:152-40:153, 40:155N, 40:172, 40:186
44						55	56				64	B60, B61, -	*40:01:07, 40:01:18, 40:53, 40:72:01-40:72:02, 40:77

44	45									55	56		58					64	-	*40:01:11
										55	56							64	B61, B41, Null, -	*40:02:01-40:04, 40:06:05, 40:20, 40:39-40:40, 40:89, 40:94, 40:97, 40:122, 40:142N-40:144N, 40:214, 40:229
										55	56							64	B4005, B50	*40:05
44										55	56							64	B61, -	*40:06:01:01-40:06:04, 40:06:06-40:06:07, 40:44, 40:70, 40:75, 40:83, 40:86, 40:103, 40:127, 40:131, 40:159, 40:162, 40:165, 40:167, 40:177, 40:190, 40:230
										55								64	B61, blank	*40:08
?										55	56							64	B61, -	*40:09, 40:11:01-40:11:02, 40:18, 40:24, 40:27:01-40:27:02, 40:29, 40:35, 40:37, 40:50, 40:56-40:57, 40:64, 40:78, 40:82, 40:85, 40:90-40:91, 40:99, 40:104-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:133Q, 40:145, 40:169, 40:173, 40:176, 40:180-40:181, 40:189, 40:200, 40:202-40:203, 40:205-40:206, 40:211, 40:220, 40:224-40:226
	45									55	56							64	B60, -	*40:10:01-40:10:02, 40:135
	45										56							64	B60, B48	*40:12
						52		54		56								63	B61, -	*40:13, 40:19
?										55	56							64	B60, -	*40:14:01-40:14:03, 40:121, 40:129, 40:158, 40:164, 40:192, 40:197, 40:219, 40:233, 40:237
?										55	56							64	B40, -	*40:15, 40:32
44										55	56							64	B61	*40:16
?	45										56							64	B40, B15	*40:21
?	45									55	56							64	Null, B60, -	*40:22N, 40:33, 40:36, 40:38, 40:42-40:43, 40:49, 40:52, 40:54-40:55, 40:61-40:63, 40:65, 40:67, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:108, 40:112-40:114, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:130, 40:134, 40:141, 40:151, 40:154, 40:156, 40:163, 40:168, 40:171, 40:175, 40:179, 40:182-40:183, 40:187, 40:191, 40:193-40:196, 40:198-40:199, 40:204, 40:207, 40:212-40:213, 40:215-40:218, 40:221-40:223, 40:227-40:228, 40:234-40:236, 40:238

?	45								55						64	B40, B60, -	*40:25, 40:69, 40:106
?									55	56					64	B21	*40:26
?				49					55	56					64	-	*40:28
?	45								55	56					64	B60, -	*40:30, 40:34, 40:160:01-40:160:02
?					51				55	56					64	B60, -	*40:31, 40:45, 40:80
?	45									56					64	B60	*40:46
	45									56				63		B40, -	*40:47, 40:188
?	45								55	56					64	B60	*40:48
?									55	56		58			64	B40	*40:58
44	45								55	56					64	B40	*40:59
?									55						64	-	*40:68, 40:232
?									55	56		59			64	B61	*40:71
?	45								55	56					64	B60, -	*40:74, 40:149, 40:170, 40:210
	45									56						B60	*40:76
44										56					64	-	*40:93
?						53			55	56					64	-	*40:95
44										56				63		-	*40:96
?									55	56					64	-	*40:98
44					52		54			56				63		-	*40:109
44										56				63		-	*40:110
	45				52		54			56				63		-	*40:117
?	45		48						55	56					64	-	*40:132, 40:178
?		46							55	56					64	-	*40:136
?	45									56					64	-	*40:137
44	45								55	56					64	-	*40:138
44						53			55	56					64	-	*40:148, 40:161
										56				63		-	*40:157, 40:201
?	45		47						55						64	-	*40:166
?									55	56					64	-	*40:174
?									55	56					64	-	*40:184
?	45									56					64	-	*40:185
?	45									56					64	-	*40:208
?										56					64	-	*40:209
?	45	46		48					55	56					64	-	*40:231
44		46							55	56					64	B41, -	*41:01, 41:03:01-41:03:02, 41:06-41:07, 41:12, 41:14-41:17, 41:20-41:22, 41:25
44		46							55	56					64	B41, -	*41:02:01-41:02:05, 41:11, 41:13, 41:19, 41:23-41:24, 41:27
44		46							55	56					64	B41, -	*41:04, 41:10, 41:18
44		46							?	56					64	B41	*41:05
44		46							55	56					64	B41	*41:08
44									55	56					64	B45, -	*41:09, 45:02-45:03 ¹⁰
44		46								56						-	*41:26
44		46	47												64	B42, -	*42:01:01-42:02, 42:06, 42:10, 42:12, 42:15, 42:17-42:18
44		46	47												64	B42, -	*42:04, 42:14
44		46	47												64	B42	*42:05:01-42:05:02
44			47												64	B42	*42:08
44		46	47												64	B42, -	*42:09, 42:16, 42:20
44		46								56					64	-	*42:11
44		46	47	48											64	-	*42:13
44			47												64	-	*42:19

				48								56	57					63	B44, Null, -	*44:02:01:01-44:02:21, 44:02:23-44:02:29, 44:11, 44:19N, 44:21-44:24, 44:27:01-44:27:02, 44:33-44:34:02, 44:48, 44:52N, 44:55, 44:59:01-44:59:02, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:84, 44:86-44:87, 44:89, 44:93, 44:101-44:102, 44:104, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126:01-44:126:02, 44:137, 44:142, 44:148-44:149N, 44:151-44:152, 44:162, 44:168, 44:170, 44:171N	
				48									57					63	B44, Null, -	*44:02:22, 44:12, 44:171N	
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.
																					ser. ¹⁷ HLA-B allele ⁴
				48								56						63	B44, Null, -	*44:03:01-44:03:06, 44:03:08-44:03:15, 44:03:18, 44:04-44:05:04, 44:13, 44:26, 44:28:01-44:30, 44:32, 44:35-44:36, 44:38-44:39, 44:42, 44:45, 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76-44:79, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114-44:115, 44:117, 44:120, 44:122-44:125, 44:128, 44:133, 44:136, 44:139, 44:141, 44:143, 44:147, 44:153-44:155, 44:157, 44:159-44:161, 44:163-44:165, 44:167, 44:174-44:175	
				48								56						63	B44, -	*44:03:07, 44:03:16, 44:03:19, 44:07, 44:81, 44:103, 44:111, 44:134	
												56						63	-	*44:03:17	
				48						54			57			61		63	B44	*44:06	
				48								56	57					63	B44	*44:08	
				48							55	56	57					64	B45	*44:09	
				48								56			59			63	B44	*44:10	
				48								56			59			63	B44	*44:14	
												56			59			63	B12	*44:15	
				48								56	57					63	B47, B44, -	*44:16, 44:91, 44:132	
				48								56						63	B44, -	*44:17, 44:43:01-44:43:02, 44:144	
							52					56			59			63	B44	*44:18	
				48								56	57					63	B44, -	*44:20, 44:100	
				48			52			54		56						63	B44	*44:25, 44:50	
45				48								56						63	B44	*44:31	
												56						63	B44, -	*44:37:01-44:37:02, 44:64:01-44:64:02	
				48								56						63	B44, -	*44:40, 44:130	
				48								56	57					63	B44, -	*44:41:01-44:41:02, 44:53, 44:99, 44:138Q	
				48								56	57					63	B44	*44:44	
				48							55	56						64	B12, -	*44:46, 44:75, 44:129	

				48								56						63	B44, -	*44:47, 44:65		
				48								56	57					63	B44	*44:49		
				48								56						63	B44	*44:54		
				48								56						63	-	*44:57		
				48								56						63	-	*44:60, 44:110		
				48								56						63	-	*44:62		
				48								56	57					63	-	*44:83, 44:145		
			47	48	50					55	56	57						64	-	*44:90		
				48			52		54		56	57						63	-	*44:95		
			47	48	50						56	57						63	-	*44:97		
				48							56	57						63	-	*44:106		
				48							56	57							-	*44:127		
				48						55	56	57					64	-	*44:131			
				48							56			59				63	-	*44:140		
				48							56						62	63	-	*44:146		
											56							63	-	*44:150		
				48							56							63	-	*44:156		
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.
																					ser. ¹⁷	HLA-B allele ⁴
		46										56						63	-	*44:166		
	45			48								56	57					63	-	*44:169		
44											55	56			59				64	B45, -	*45:01, 45:05, 45:07, 45:11-45:14	
44											55	56			59				64	B21	*45:04	
44			47								55	56			59				64	B45	*45:06	
44											55	56		58					64	-	*45:08	
44											55	56							64	-	*45:09	
44											55	56							64	-	*45:10	
			47		49													62		64	B46, Null, -	*46:01:01-46:01:09, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:20, 46:22-46:24, 46:26-46:31, 46:34-46:38
			w		49													62		64	B46	*46:02
			47		49															64	B46, -	*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25
45		47		49																64	-	*46:06
		47		49											59					64	B22	*46:11
		47		49														62		64	-	*46:12
		47		49														62		64	-	*46:17
		47		49										58						64	-	*46:18
		47		49					w									62		64	-	*46:32
		47		49																64	-	*46:33
					50							56						63	B47, -	*47:01:01-47:01:02, 47:06-47:08		
					50						55	56								64	B47	*47:02
					50						w	56						63	w	B47	*47:03	
												56						63		B17, -	*47:04-47:05	
44						51						56								64	B48, -	*48:01:01-48:01:04, 48:07, 48:09, 48:11, 48:15, 48:19-48:20, 48:22, 48:27, 48:29-48:30
44												56								64	B48, B70, -	*48:02:01-48:02:03, 48:25
44	45											56								64	B48, -	*48:03:01-48:03:02, 48:23
44						51						56								64	B48, -	*48:04, 48:24
44												56								64	B48, -	*48:05, 48:08, 48:10, 48:12, 48:14
44						51														64	-	*48:06
44												56								64	-	*48:13

44					51								64	B48	*48:16		
44	45								56				64	-	*48:17		
					51				56				63	B48	*48:18		
44									56				64	-	*48:21, 48:26		
44					51				56				64	-	*48:28		
						52	53			56		59		63	B49, Null, -	*49:01:01, 49:01:03-49:01:04, 49:06, 49:08-49:10, 49:12-49:17, 49:19N, 49:22-49:24	
						52				56		59		63	B49, -	*49:01:02, 49:11	
							53			56		59		63	B49	*49:02	
							53			56		59		63	B49	*49:03	
						52	53			56				63	B49	*49:04-49:05	
						52	53			56				63	-	*49:07	
						52	53			56				63	-	*49:18	
						52	53			56		59		63	-	*49:20	
						52	53			56		59			-	*49:21	
						52	53			56		59		63	-	*49:25	
							53					59		63	-	*49:26	
44							53		55	56		59			64	B50, -	*50:01:01-50:01:03, 50:04-50:05, 50:08, 50:10-50:11, 50:13, 50:16, 50:18-50:19, 50:32
44							53		55	56		59			64	B45	*50:02
44							53			56		59			64	-	*50:06
44							53		55	56		59			64	-	*50:07
44							53		55	56					64	-	*50:09, 50:20
44									55	56		59			64	-	*50:12
44							53		55	56					64	-	*50:14
44							53			55		59			64	-	*50:15
44									55	56		59			64	-	*50:17
44							53		55	56					64	-	*50:31
								53	54				61		63	B51, B5102, B5103, Null, blank, -	*51:01:01-51:01:03, 51:01:05-51:01:08, 51:01:10-51:01:20, 51:01:22-51:01:24, 51:01:26-51:01:39, 51:01:41-51:02:03, 51:02:05-51:03, 51:05, 51:07:01-51:07:02, 51:09:01-51:09:02, 51:11N-51:12, 51:13:02-51:14, 51:16-51:19, 51:21-51:23, 51:26-51:33, 51:37-51:41N, 51:43, 51:48-51:53, 51:55, 51:57-51:58, 51:60, 51:63, 51:65-51:67, 51:69-51:77, 51:79-51:80, 51:82-51:91, 51:94-51:96, 51:98N, 51:100, 51:102-51:103, 51:105, 51:107-51:111, 51:113-51:117, 51:119-51:130, 51:132-51:134, 51:136-51:138, 51:140, 51:142, 51:145-51:146, 51:149N-
								53	54				61		63	B51, B5102, -	*51:01:04, 51:02:04, 51:13:01, 51:24:01, 51:34-51:35, 51:92-51:93, 51:99, 51:131, 51:143-51:144
44								53	54				61		63	B51, -	*51:01:09, 51:01:21, 51:01:25, 51:24:03-51:24:04
								53					61		63	-	*51:01:40, 51:54, 51:78:01-51:78:02
				49					54				61		63	B51, -	*51:04, 51:46, 51:56:01-51:56:02, 51:139

						54				61	63	B51, -	*51:06:01-51:06:02, 51:59, 51:135
					53	54				61	63	B51, Null, -	*51:08, 51:20, 51:44N, 51:97, 51:141
44						54				61	63	B51	*51:10
						54			59	61	63	-	*51:15, 51:62
44					53	54				61	63	B51	*51:24:02
					53	54				61	63	B51	*51:36
		48				54				61	63	-	*51:42
						54				61	63	-	*51:45, 51:81
					53	54				61	62 63	-	*51:61:01-51:61:02
						54				61	63	-	*51:64
					53	54				61	63	-	*51:68
					53	54				61	63	-	*51:101
					53	54					63	B17, B5, Null, -	*51:104, 51:118N, 51:147, 58:08:01-58:08:02 ¹¹
					53	54			59	61	63	-	*51:106
				52	53	54					63	-	*51:112
						54				61	63	-	*51:148, 53:14 ¹²
					53	54	56				63	B52, -	*52:01:01:01-52:01:04, 52:01:06-52:01:09, 52:01:11-52:01:14, 52:01:16-52:13, 52:15, 52:17-52:18, 52:22-52:24. 52:27-52:31
					53	54	56				63	B52, -	*52:01:05, 52:01:15, 52:14
					53	54	56				63	-	*52:01:10
					53	54	56				63	-	*52:16
					53	54	56				63	-	*52:19
					53		56				63	-	*52:20
					53	54	56				62 63	-	*52:21
					53	54	56		59		63	-	*52:25
					53	54	56				63	-	*52:26
						54				61	63	B53, -	*53:01:01-53:01:08, 53:04-53:05, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:27. 53:29
			49			54				61	63	B53, B53, B51	*53:02, 53:06
										61	63	B53, B35	*53:03, 53:09, 53:11-53:13
						54				61	63	-	*53:15
						54	56				63	-	*53:17:01-53:17:02
		48				54				61	63	-	*53:22
			49			54	56				63	-	*53:28
				52		54					63	-	*53:30, 57:45, 57:51 ¹³
		47		50				57	58		64	B54, Null, -	*54:01:01, 54:01:03-54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:25. 54:27-54:28
		47		50					58		64	B55, B54, -	*54:01:02, 55:01:07, 55:02:01-55:02:06, 55:07, 55:10, 55:12, 55:16, 55:19, 55:26, 55:30, 55:35, 55:37, 55:39, 55:41-55:43, 55:47-55:48, 55:50, 55:57. 56:10 ¹⁴
		47		50				57	59		64	B22	*54:03
		47		50				57			64	B54, -	*54:04, 54:09, 54:14-54:15
		47		50				57		62	64	-	*54:06
44		47		50				57			64	B54	*54:11
		47		50	52			57	58		63	-	*54:12
		47		50		53		57	58		64	-	*54:20
		47		50				57			64	-	*54:26

		47		50		53			58			64	B55, Null, -	*55:01:01-55:01:06, 55:01:08-55:01:10, 55:03, 55:05, 55:11, 55:15, 55:25, 55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45. 55:52-55:55N
		47		50								64	B55, B56, Null, -	*55:04, 55:08, 55:13, 55:23, 55:27, 55:32, 55:46, 55:49, 56:01:05, 56:15, 56:18-56:19N. 56:22. 56:31-56:32 ¹⁵
		47		50		53						64	B22, B55	*55:09, 55:24
		47		50					60			64	-	*55:14
		47		50		53						64	B55, -	*55:17, 55:28, 55:58
								56	58			64	-	*55:18
	46	47		50								64	B22	*55:20
		47		50		w			58		62	64	B55	*55:21
		47		50								64	B22	*55:22
		47		50				56	58			64	-	*55:34
		47		50		53		56	58			64	-	*55:40
		47		50								64	-	*55:51
		47		50					58			64	-	*55:56
		47		50					59			64	B56, Null, -	*56:01:01-56:01:04, 56:02, 56:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:36. 56:38N
		47		50							62	64	B22, 15	*56:03
		47		50		53			58			64	B56	*56:05:01
		47		50		53						64	B56	*56:05:02
		47		50		53			58		w	64	B78	*56:06
		47		50					59			63	B56	*56:07
		47		50								64	B56, B22	*56:09, 56:11-56:12
		47		50					59			64	B56	*56:13
									59			64	B56, -	*56:17, 56:33
44		47		50					59			64	-	*56:20:02
		47		50		52 53						63	-	*56:21
									58			64	-	*56:23
		47		50		53			59			64	-	*56:25, 56:37
		47		50	w				59			64	-	*56:39
								54		60		63	B57, Null, -	*57:01:01-57:01:04, 57:01:06-57:03:02, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25-57:30, 57:32-57:44, 57:46, 57:48-57:50, 57:52-57:58, 57:60-57:64, 58-36 ¹⁶
								54		60		63	B57, -	*57:01:05, 57:11, 57:13, 57:21, 57:31, 57:47
								54		60		63	B57	*57:04
			48					54		60		63	B57, -	*57:07, 57:24
			48					54		60		63	B57	*57:09
										60		64	B57	*57:12
				49				54		60		63	-	*57:14
45								54		60		63	-	*57:59
								54				63	B58, Null, -	*58:01:01-58:01:02, 58:01:04-58:01:12, 58:04-58:05, 58:10N-58:15, 58:19, 58:21-58:24, 58:28-58:29, 58:31N-58:35, 58:37, 58:39N-58:40, 58:42
								54				63	B58, -	*58:01:03, 58:02, 58:06, 58:16:01, 58:25-58:27, 58:38
								54				63	B58	*58:07

				49					54							63	B17	*58:09				
									54							63	Null	*58:17N				
									54							63	-	*58:18				
									54							63	-	*58:20				
							52					58				63	B59, -	*59:01:01:01-59:01:01:02, 59:05				
							52									63	B59, -	*59:02-59:03				
							52					59				63	-	*59:04				
44			47													64	B67, -	*67:01:01, 67:03				
44			47													64	B67, -	*67:01:02-67:02				
44			47												62	64	B73, -	*73:01-73:02				
								53				58		61		64	B78, -	*78:01:01-78:01:02, 78:02:02-78:03, 78:07				
								53						61		64	B78, B35	*78:02:01, 78:04				
								53			56					64	-	*78:05				
								53			56					64	-	*78:06				
			47			50	51									64	B81	*81:01				
			47				51									64	Null, -	*81:03-81:06				
			47			50							59			64	B82, -	*82:01-82:03				
			47	48		50						57				64	-	*83:01				
																		A*23:31, A*24:106, C*07:231, C*16:10				
																		A*24:174				
															62			A*26:68, A*68:56, C*02:56, C*06:20, C*12:50				
																		C*01:30, C*12:87				
					49													C*01:73				
																		C*03:05, C*03:25, C*03:27, C*03:143, C*03:167				
													59					C*03:19				
													59					C*03:102				
45																		C*03:129, C*17:16				
																		C*03:179				
														60				C*06:72				
																		C*07:02:30, C*08:16:02				
																		C*07:46				
																		C*07:302				
													58					C*15:02:04				
				48														C*15:25				
				49														C*15:39				
																		C*15:51				
44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		Well No.