

INTERPRETATION TABLE			
HLA-B*27 SSP typing			
Amplification pattern of the B*2701 to 2759N alleles ¹			
	Well		
	1	2	
Length of spec.	145	95	Length of spec.
PCR product			PCR product
Length of int.	430	515	Length of int.
pos. control ²			pos. control
5'-primer ³	167	363	5'-primer ³
	5'-gCT ^{3'}	5'-AAT ^{3'}	
3'-primer ⁴	272	418	3'-primer ⁴
	5'-TgC ^{3'}	5'-gTC ^{3'}	
Well No.	1	2	Well No.
HLA-B allele ⁵			HLA-B allele ⁵
*2701-270508, 270510-270512, 2708, 2710, 2713, 2715, 2717, 2725, 2728, 2736-2740, 2742, 2744, 2745, 2747-2759N	1	2	*2701-270508, 270510-270512, 2708, 2710, 2713, 2715, 2717, 2725, 2728, 2736-2740, 2742, 2744, 2745, 2747-2759N
*270509, 2712, 2716, 2718, 2723, 2726, 2729, 2731, 3702, 4704, 4705		2	*270509, 2712, 2716, 2718, 2723, 2726, 2729, 2731, 3702, 4704, 4705
*2706, 2707, 2709, 2711, 2714, 2719-2721, 2724, 2727, 2730, 2732-2735, 2741, 2743, 2746	1		*2706, 2707, 2709, 2711, 2714, 2719-2721, 2724, 2727, 2730, 2732-2735, 2741, 2743, 2746
HLA-B allele ⁵			HLA-B allele ⁵
Well No.	1	2	Well No.