

INTERPRETATION TABLE		
HLA-B*27 SSP typing		
Amplification pattern of the B*27:01 to 27: 86 alleles ¹		
	Well	
	1	
Length of spec.	145	Length of spec.
PCR products	205	PCR products
Length of int.	430	Length of int.
pos. control ²		pos. control
5'-primer ³	167	5'-primer ³
	5' -gCT 3'	
	363	
	5' -AAT 3'	
3'-primer ⁴	272	3'-primer ⁴
	5' -TgC 3'	
	272	
	5' -TgC 3'	
	527	
	5' -CCT 3'	
Well No.	1	Well No.
HLA-B allele ⁵		HLA-B allele ⁵
*27:01-27:21, 27:24-27:28, 27:30-27:74, 27:76, 27:78-27:84, 27:86	1	*27:01-27:21, 27:24-27:28, 27:30-27:74, 27:76, 27:78-27:84, 27:86
*44:97	1	*44:97
HLA-B allele ⁵		HLA-B allele ⁵
Well No.	1	Well No.