

INTERPRETATION TABLE	
HLA-B*27 SSP typing	
Amplification pattern of the B*27:01 to 27:100 alleles ¹	
	Well
	1
Length of spec.	145
PCR products	205
Length of int.	430
pos. control ²	
5'-primer ³	167
	5' -gCT 3'
	363
	5' -AAT 3'
3'-primer ⁴	272
	5' -TgC 3'
	272
	5' -TgC 3'
	527
	5' -CCT 3'
Well No.	1
HLA-B allele ⁵	
*27:01-27:21, 27:24-27:28, 27:30-27:74, 27:76, 27:78-27:84, 27:86-27:100	1
*44:97	1
HLA-B allele ⁵	
Well No.	1