

INTERPRETATION TABLE								
DPA1 SSP typing								
Amplification patterns of the DPA1 alleles								
	Well ⁵							
	1	2	3	4	5	6	7	8
Length of spec.	85	250	200	115	105	155	100	100
PCR product						195		
Length of int.	515	515	430	430	430	515	430	430
pos. control ¹								
5'-primer(s) ²	15	11	28	4	84	31	11	11
	5'-ACg ^{3'}	5'-C gC ^{3'}	5'-gAA ^{3'}	5'-Cg g ^{3'}	5'-AAT ^{3'}	5'-g CA ^{3'}	5'-C gC ^{3'}	5'-C AT ^{3'}
3'-primer(s) ³	31	83	83	28	intr ⁶	69	31	31
	5'-CAT ^{3'}	5'-ggT ^{3'}	5'-ggT ^{3'}	5'-TC g ^{3'}	5'-ggC ^{3'}	5'-gTC ^{3'}	5'-CTg ^{3'}	5'-CTg ^{3'}
						83		
						5'-ggT ^{3'}		
Well No.	1	2	3	4	5	6	7	8
DPA1 allele ⁴								
*010301-010302	1	2	3					
*010303		2	3					
*0104	1	2		4				
*0105	1				5			
*0106		2	3			6	7	
*0107	1	2	3					
*0108	1	2		4				
*0109	1	2	3					
*020101-020106					5		7	
*020201-020203					5			8
*0203					5			
*0204					5	6		8
*0301			3					
*0302			3					
*0303				4				
*0401	1				5			
DPA1 allele ⁴								
Well No.	1	2	3	4	5	6	7	8

INTERPRETATION TABLE								
DPA1 SSP typing								
Amplification patterns of the DPA1 alleles								
Well ⁵								
9	10	11	12	13	14	15	16	
205	85	90	205	135	140	245	220	Length of spec.
								PCR product
430	515	515	430	430	430	430	430	Length of int.
								pos. control ¹
11	15	66	18	51	50	15	23	5'-primer(s) ²
5'-C AT ^{3'}	5'-ACC ^{3'}	5'-A TC ^{3'}	5'-gA A ^{3'}	5'-CA T ^{3'}	5'-C Cg ^{3'}	5'-ACC ^{3'}	5'-T AC ^{3'}	
	15					15		
	5'-ACC ^{3'}					5'-ACC ^{3'}		
66	31	83	73	83	83	83	83	3'-primer(s) ³
5'-T CA ^{3'}	5'-CAT ^{3'}	5'-ggT ^{3'}	5'-AgC ^{3'}	5'-ggT ^{3'}	5'-ggT ^{3'}	5'-ggT ^{3'}	5'-ggT ^{3'}	
9	10	11	12	13	14	15	16	Well No.
								DPA1 allele ⁴
								*010301-010302
	10					15		*010303
								*0104
								*0105
								*0106
				13				*0107
					14			*0108
							16	*0109
								*020101-020106
9								*020201-020203
	10							*0203
9								*0204
	10	11				15		*0301
9	10					15		*0302
	10	11				15		*0303
			12					*0401
								DPA1 allele ⁴
9	10	11	12	13	14	15	16	Well No.