

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
DPA1 allele ^{1,2,3}																				
*01:03:01:01-01:03:02, 01:03:04-01:03:05	1	2	3																	
*01:03:03		2	3							10					15					
*01:04	1	2		4																
*01:05	1				5															
*01:06:01-01:06:02		2	3			6	7													
*01:07	1	2	3								11			14						
*01:08	1	2		4										14						
*01:09	1	2	3										13							
*01:10	1	2	3			6														
*01:11	1	2	3											14						
*01:12	1	2	3								11									
*01:13	1		3			6														
*01:14	1	2	3																19	
*02:01:01:01-02:01:02:02, 02:01:08:01-02:01:08:02					5		7										17			
*02:01:03-02:01:07					5		7									?	?			
*02:02:02:01-02:02:02:04					5			8	9								17			
*02:02:03-02:02:06					5			8	9							?	?			
*02:03					5				10							?	?			
*02:04					5	6		8	9							?	?			
*02:05					5			8	9					14		?	?			
*02:06					5			8	9				13				17			
*02:07:01:01-02:07:01:03					5			8	9							16	17			
*02:08					5		7										17		19	
*02:09					5		7										17	18		
*03:01:01			3							10	11				15		17			
*03:01:02			3							10	11				15					
*03:02			3						9	10					15		?			
*03:03				4						10	11				15		?			
*04:01	1				5							12								
*04:02	1				5							12					17			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20

Negative Control