

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
DPA1 allele ^{1,2,3}																						
*01:03:01:01-01:03:02, 01:03:04-01:03:05	1	2	3																			
*01:03:03		2	3						10						15							
*01:04	1	2		4																		
*01:05	1				5																	
*01:06:01-01:06:02		2	3			6	7															
*01:07	1	2	3							11				14								
*01:08	1	2		4										14								
*01:09	1	2	3										13									
*01:10	1	2	3			6																
*01:11	1	2	3											14								
*01:12	1	2	3							11												
*01:13	1		3			6																
*01:14	1	2	3																19			
*01:15	1	2	3															18				
*01:16	1	2	3				7															
*01:17	1	2	3																	20		
*02:01:01:01-02:01:02:02, 02:01:08:01-02:01:08:02					5		7										17					
*02:01:03-02:01:07					5		7										?					
*02:02:02:01-02:02:02:04					5			8	9								17					
*02:02:03-02:02:06					5			8	9								?					
*02:03					5					10							?					
*02:04					5	6		8	9								?					
*02:05					5			8	9					14			?					
*02:06					5			8	9				13				17					
*02:07:01:01-02:07:01:03, 02:12					5			8	9						16		17					
*02:08					5		7										17		19			
*02:09					5		7										17	18				
*02:10					5			8	9								17					21
*02:11			3		5		7										?					
*03:01:01			3							10	11				15		17					
*03:01:02			3							10	11				15							
*03:02			3							9	10				15		?					
*03:03				4						10	11				15		?					
*04:01	1				5							12										
*04:02	1				5							12					17					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22

Negative Control