

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
DQB1 allele ^{1,2}	ser. ³															
*02:01:01-02:02:12, 02:04-02:76, 02:78-02:167N	DQ2, Null, -			3												
*02:03:01, 02:77	DQ2, -			3			6									
*02:03:02	-		2	3												
*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:01:22, 03:01:03-03:01:06, 03:01:09-03:01:48, 03:04:01-03:04:04, 03:09, 03:10:03, 03:16, 03:19:01:01-03:19:04, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:57, 03:59-03:60, 03:69, 03:71, 03:76-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:119-03:122, 03:127-03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:148, 03:150, 03:152, 03:157-03:160, 03:163-03:167, 03:169-03:173, 03:180, 03:182, 03:186-03:188, 03:191-03:193, 03:196-03:198:02, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231, 03:235-03:236, 03:241-03:243, 03:246, 03:253-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275-03:276N, 03:281, 03:284-03:286, 03:288, 03:290, 03:292-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311-03:312, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353, 03:358N, 03:360-03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:381, 03:385N, 03:387, 03:389-03:390, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408	DQ7, Null, -				4			7				12				
*03:01:02, 03:01:08, 03:73, 03:75, 03:109, 03:118N, 03:154, 03:162, 03:252, 03:291, 03:335, 03:354N, 03:391	DQ7, Null, -				4			7								
*03:01:07, 03:13, 03:58, 03:194, 03:232	DQ7, -				4			w					12			
*03:02:01:01-03:02:02, 03:02:04-03:02:05, 03:02:07-03:02:12, 03:02:14-03:02:17, 03:02:19-03:02:29, 03:02:31, 03:05:04, 03:07, 03:11, 03:18, 03:32, 03:45:01-03:45:02, 03:63-03:64, 03:66N-03:68, 03:85, 03:106, 03:146, 03:153, 03:161, 03:174-03:175, 03:178, 03:184-03:185, 03:190, 03:199, 03:203, 03:205, 03:210-03:211, 03:213N-03:215, 03:220-03:221, 03:224, 03:229, 03:240, 03:245, 03:247, 03:251, 03:261, 03:263, 03:265, 03:269N, 03:273-03:274, 03:277-03:279, 03:287, 03:289, 03:295-03:296, 03:298-03:301, 03:310N, 03:315, 03:320-03:324, 03:333-03:334N, 03:339N, 03:343-03:345, 03:348-03:349, 03:362, 03:364, 03:367-03:369, 03:379, 03:383, 03:386, 03:388, 03:392, 03:403N, 03:409-03:410, 03:412	DQ8, Null, -					5		7					12			
*03:02:03, 03:02:13, 03:05:02, 03:37, 03:189, 03:204, 03:223, 03:226, 03:233, 03:237N, 03:352, 03:371	DQ8, Null, -					5		7								
*03:02:06, 03:02:18, 03:02:30, 03:125, 03:228	-					5							12			
*03:03:02:01-03:03:02:05, 03:03:04-03:03:10, 03:03:12-03:03:15, 03:03:17-03:03:25, 03:15, 03:23:03, 03:26, 03:30-03:31, 03:33-03:34, 03:38:01, 03:39-03:41, 03:43, 03:74, 03:79, 03:86-03:91Q, 03:95N-03:98, 03:111, 03:113, 03:117, 03:123-03:124, 03:126, 03:136, 03:141, 03:145, 03:149, 03:155, 03:168, 03:176-03:177, 03:200, 03:209, 03:212, 03:222, 03:227, 03:230, 03:234, 03:239, 03:248-03:249, 03:270, 03:280, 03:282N, 03:304, 03:313, 03:316, 03:319, 03:332, 03:337, 03:356N-03:357N, 03:359, 03:363, 03:365, 03:374-03:375N, 03:382, 03:384, 03:393, 03:395, 03:397-03:398, 03:402, 03:405-03:406	DQ9, Null, -						6	7					12			
*03:03:03, 03:03:16, 03:20, 03:38:02, 03:99Q, 03:105, 03:112, 03:156, 03:238, 03:258, 03:336, 03:351	DQ9, -						6	7								
*03:03:11	-						6	7				11	12			
*03:03:26, 03:62, 03:78, 03:81, 03:110, 03:151, 03:225, 03:244, 03:272, 03:308, 03:314	-							7					12			
*03:05:01, 03:05:05, 03:61, 03:181, 03:250, 03:346	DQ8, -					5		7				11				
*03:05:03	DQ8					5		7				11	12			
*03:06, 03:25:01-03:25:02, 03:411N	DQ3, Null, -						6						12			
*03:08, 03:107	DQ8, -					5		w					12			

*03:10:01-03:10:02:02, 03:14:01-03:14:02, 03:195	DQ3, -				4				7				12			15
*03:12	DQ9							6	7				12			15
*03:17:01	DQ9								7			11	12			
*03:17:02	-								7							
*03:23:01	-			2					7							
*03:23:02, 03:217, 03:259	-			2					7				12			
*03:65, 03:137, 03:283	-							6	w				12			
*03:70, 03:179	-						5		7				12			15
*03:72, 03:100	-								7			11				
*03:104	-						5	6	7							
*03:132	-						5		7	8						
*03:147	-				4								12			
*03:183	-				4				7							15
*03:262	-						5		w							
*03:355	-			2		4			7				12			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
*04:01:01:01-04:01:02, 04:01:04-04:02:01:01, 04:02:01:04-04:02:15, 04:02:17-04:02:18, 04:02:20, 04:06-04:09, 04:11-04:21, 04:23, 04:25N-04:32, 04:34-04:37, 04:39-04:48, 04:50-04:74, 04:76-04:77	DQ4, Null, -								8	9						
*04:01:03, 04:02:16	-								8	9			12			
*04:02:19, 04:04-04:05, 04:33, 04:49	-								8							
*04:03:01-04:03:02	-						6		8	9						
*04:03:03	-						6			9			12			
*04:10	-		2						8	9				13		
*04:22, 04:24, 04:38	-								8	w						
*04:75	-					5			8	9						
*05:01:01:01-05:01:15, 05:01:17-05:03:23, 05:05:01-05:43:02, 05:45-05:51, 05:53, 05:55-05:59, 05:61-05:71, 05:73-05:76, 05:78-05:81, 05:84-05:97, 05:99-05:104, 05:106-05:113, 05:115, 05:117-05:127, 05:129-05:131, 05:133-05:145, 05:147-05:174, 05:177-05:206N, 05:208N-05:217, 05:219-05:243, 06:325 ⁴	DQ5, Null, -	1									10					
*05:01:16, 05:60, 05:82-05:83, 05:98, 05:116, 05:175, 05:218	-										10					
*05:04, 05:44, 05:52, 05:54, 05:72, 05:77, 05:105, 05:114, 05:128N, 05:132Q, 05:146, 05:207	DQ5, Null, -	1														
*05:176	-	1	2													
*06:01:01:01-06:01:06, 06:01:08-06:01:26, 06:43, 06:54N-06:58, 06:98, 06:100-06:105, 06:108, 06:120, 06:132, 06:140, 06:153:01-06:153:02, 06:157, 06:167, 06:177, 06:181, 06:194, 06:205, 06:214, 06:229, 06:239, 06:243, 06:245-06:247, 06:251, 06:257-06:258, 06:260, 06:263, 06:268, 06:274, 06:277, 06:285, 06:302, 06:305, 06:307, 06:309-06:310, 06:312, 06:321, 06:323, 06:330N, 06:342, 06:359	DQ6, Null, -		2													15
*06:01:07, 06:03:11, 06:04:01:01-06:05:01, 06:07:01-06:07:02, 06:09:01:01-06:09:12, 06:12, 06:17, 06:21-06:22:03, 06:25, 06:28, 06:34-06:36, 06:38-06:39, 06:42, 06:45, 06:52-06:53:02, 06:69:01-06:69:02, 06:85-06:																

[illegible]