

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DRB1 allele^{1,2,3}																								
*01:01:01-01:01:31, 01:38, 01:41, 01:45, 01:47-01:49, 01:53, 01:56, 01:59-01:60, 01:64-01:66, 01:69-01:73, 01:75, 01:78-01:80, 01:82, 01:84	1		3																					
*01:02:01-01:02:05, 01:02:07-01:02:13, 01:43, 01:46, 01:57, 01:76, 01:81, 01:83, 01:86-01:88		2	3																				23	
*01:02:06, 01:58			3																				23	
*01:03:01-01:03:02	1			4			7																	
*01:04			3		5				9															
*01:05	1		3			6																		
*01:06			3						9	10														
*01:07	1		3					8																
*01:08	1		3						9															
*01:09	1		3							10														
*01:10	1		3	4							11													
*01:11:01-01:11:02	1		3		5																			
*01:12, 01:37	1		3									12												
*01:13	1		3										13											
*01:14	1		3											14										
*01:15	1						7			10														
*01:16	1		3		5										15									
*01:17	1		3				7																	
*01:18:01-01:18:02	1		3								11													
*01:19	1		3													16								
*01:20:01, 01:85			3						9															
*01:20:02		2	3						9															
*01:21	1		3												15	16								
*01:22	w		w														17							
*01:23		2											13	14									23	
*01:24:01-01:24:02	1													14		16								
*01:25	1		3															18						
*01:26		2	3									12											23	
*01:27	1		3					8									17							
*01:28	1		3										13		15									
*01:29:01	1		3		5	6																		
*01:29:02	1					6																		
*01:30	1		3														17							
*01:31	1		3			6						12												
*01:32	1		3																		21			
*01:33N																			20					
*01:34		2	3																		21		23	
*01:35			3		5																		23	
*01:36	1		3																			22		
*01:39N	1			4			7												19					
*01:40N, 01:67 ⁴	1		3															19						
*01:42				4			7																23	
*01:44:01-01:44:02		2	3								11												23	
*01:50	1		3					8										18						
*01:51	1												13	14										
*01:52N	1		3							10												22		
*01:54		2												14		16							23	
*01:55	1				5										15									
*01:61				4			7		9															
*01:62N	1		3	4																	21			
*01:63			3																					
*01:68N	1		3								11			14										
*01:74		2	3		5	6																	23	
*01:77	1		3															19	20					
*01:89	1			4	5		7																	
*13:165, DRB3*02:63																	17							
*14:112								8																
*14:196																		19	20					
DRB3*01:47								8									17							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control