

INTERPRETATION TABLE			
DRB1*10 SSP subtyping			
Amplification patterns of the DRB1*10 alleles			
	Tube		
	1	2	3
Length of spec.	210	205	205
PCR product			
Length of int.	515	515	430
pos. control¹			
5'-primer²	31	31	31
	5'-gC g ^{3'}	5'-gC g ^{3'}	5'-gC g ^{3'}
3'-primer³	87	86	86
	5'-CTC ^{3'}	5'-C AC ^{3'}	5'-C CA ^{3'}
Tube No.	1	2	3
DRB1 allele			
*100101-100102	+	+	
*1002	+		+
DRB1 allele			
Tube No.	1	2	3