

INTERPRETATION TABLE				
DRB1*10 SSP subtyping				
Amplification patterns of the DRB1*10 alleles				
	Well			
	1	2	3	4
Length of spec.	210	205	205	70
PCR product				
Length of int.	515	515	430	430
pos. control¹				
5'-primer²	31(178)	31(178)	31(178)	31(178)
	5' -gCg 3'	5' -gCg 3'	5' -gCg 3'	5' -gCg 3'
3'-primer³	87(346)	86(344)	86(344)	40(206)
	5' -CTC 3'	5' -CAC 3'	5' -CCA 3'	5' -CgA 3'
Well No.	1	2	3	4
DRB1 allele				
*10:01:01-10:01:03	1	2		
*10:02	1		3	
*10:03	1	2		4
DRB1 allele				
Well No.	1	2	3	4