

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
DRB1 allele^{1,2,3}																						
*07:01:01:01-07:01:23, 07:31, 07:33, 07:35- 07:38, 07:41, 07:43, 07:45, 07:47, 07:49, 07:52-07:55, 07:57, 07:59-07:60, 07:62- 07:63, 07:65-07:67, 07:69-07:111	1	2																				
*07:03	1		3																			
*07:04	1			4																		
*07:05	1	2			5																	
*07:06, 07:32⁴	1	2				6																
*07:07	1	2					7															
*07:08	1	2						8														
*07:09	1	2	3																			
*07:10N	1	2		4																		
*07:11	1				5																	
*07:12	1	2							9													
*07:13	1	2								10												
*07:14, 07:24⁴	1	2									11											
*07:15	1	2										12										
*07:16	1	2											13									
*07:17, 07:27⁴	1	2												14								
*07:18	1	2							10						15							
*07:19	1	2													15							
*07:20	1	2												14	15							
*07:21	1	2										12	13									
*07:22		2							9													
*07:23	1	2							9			12										
*07:25				4	5																	
*07:26N	1	2	3					8														
*07:28	1	2					7								15							
*07:29	1	2				6	7															
*07:30, 07:42⁴	1	2																	19			
*07:34	1	2								10	11											
*07:39, 07:68N⁴	1	2																		20		
*07:40, 07:56⁴	1	2															17					
*07:44, 07:48⁴	1	2																			21	
*07:46	1	2																18				
*07:50	1																					
*07:51	1	2						8					13									
*07:58N	1	2													16							
*07:61																		18				
*07:64		2																				
*04:90, 09:08									9													
*09:35										10	11											

Negative Control