

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DRB1 allele^{1,2,3}																								
*08:01:01-08:01:02, 08:01:04-08:01:07, 08:50, 08:77, 08:79, 08:86, 08:88, 08:90, 08:93, 08:96	1		3		5		7																	
*08:02:01:01-08:02:04, 08:42, 08:72-08:73, 08:82, 08:84-08:85, 08:89N		2	3		5		7																	
*08:03:02:01-08:03:11, 08:46, 08:51, 08:53, 08:58, 08:61-08:62, 08:66, 08:68:01-08:69, 08:71, 08:74, 08:76, 08:91, 08:95, 08:97	1			4	5		7									16								
*08:04:01, 08:04:04-08:04:07, 08:59, 08:83, 08:94		2	3		5			8									17						23	
*08:04:02-08:04:03		2	3		5			8									w							
*08:05	1		3			6	7																	
*08:06, 08:54, 08:70	1		3		5			8									17			20			23	
*08:07			3		5		7											19						
*08:08			3		5		7				11													
*08:09		2	3		5		7											18						
*08:10	1			4	5			8								16	17			20			23	
*08:11			3		5		7						13											
*08:12	1			4	5			8		10						16	17			20				
*08:13		2			5		7					12												
*08:14	1			4			7		9							16	17							
*08:15				4	5		7				11													
*08:16	1		3		5		7												19					
*08:17	1		3		5		7							14							21			
*08:18, 08:47	1			4		6	7									16								
*08:19				4	5		7								15							22		
*08:20, 14:12:01-14:12:02, 14:84, 14:156, 14:173⁴					5												17	18						
*08:21		2	3				7											18						
*08:22	1		3		5			8		10							17			20				
*08:23	1			4	5		7									16						22		
*08:24:01-08:24:02		2	3			6	7																	
*08:25				4		6	7									15								
*08:26	1		3		5		7		9															
*08:27	1			4	5		7					12				16								
*08:28		2	3		5			8						14			17				21		23	
*08:29	1			4			7									16						22		
*08:30:01-08:30:03		2		4	5		7																	
*08:31			3			6		8			11												23	
*08:32	1			4	5		7									16		18						
*08:33	1			4	5		7						13			16								
*08:34				4	5		7								15									
*08:35	1			4	5		7		9							16		18						
*08:36:01-08:36:02	1			4	5		7		9							16								
*08:37	1			4	5		7							14		16					21			
*08:38	1			4	5		7									16			19					
*08:39	1		3		5		7						13											
*08:40, 08:92	1			4		6	7			10						16								
*08:41			3			6	7				11													
*08:43	1		3				7					12											23	
*08:44		2	3		5		7					12												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*08:45:01-08:45:02		2		4	5		7							14							21			
*08:48	1				5		7					12												
*08:49	1			4	5		7								15	16			19					
*08:52			3		5		7								15									
*08:55, 08:64	1		3		5																			
*08:56	1			4			7									16								
*08:57					5																			
*08:60N	1			4	5		7							14		16								
*08:63				4	5		7				11			14							21			
*08:65	1			4	5											16								
*08:67		2	3					8									17						23	
*08:75			3		5			8			11						17						23	
*08:78N	1		3		5		7														21			
*08:80					5												17							
*08:81		2			5		7																	
*08:87	1		3		5		7											18						
*03:12, 03:158, 14:65, 14:169-14:170									9											20				
*03:92		2						8						14									23	

Negative Control

[illegible]

*14:78					5				9							16	17	18		20				
*14:85					5				9									18						
*14:126:01-14:126:02		2						8														22	23	
*14:138								8			11			14								22	23	
*14:175		2						8						14								22	23	
*15:21																	17				21			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

[illegible]