

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
DRB1 allele ^{1,2,3}												
*09:01:02-09:01:11, 09:26	1	2										
*09:02:01-09:02:02	1		3									
*09:03, 09:09, 09:13, 09:15, 09:17-09:18 ⁴	1	2		4								
*09:04	1	2			5							
*09:05		2				6						
*09:06	1	2					7					
*09:07, 09:20		2						8				
*09:08	1				5							
*09:10-09:11, 09:14	1	2				6						
*09:12, 09:19	1	2						8				
*09:16		2		4								
*09:21	1	2								10		
*09:22	1	2							9			
*09:23	1										11	
*09:24	1											
*01:65					5	6						
*04:90, 04:155, 04:167, 07:50, DRB5*01:12, DRB5*01:15											11	
*07:52								8				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Negative Control