

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB1 allele^{1,2,3}																
*09:01:02:01 -09:01:11, 09:26, 09:28-09:30, 09:33-09:41	1	2														
*09:02:01-09:02:02	1		3													
*09:03, 09:15	1	2		4												
*09:04	1	2			5											
*09:05		2				6										
*09:06	1	2					7									
*09:07, 09:20		2						8								
*09:08	1				5											
*09:09, 09:13, 09:17-09:18 ⁴	1	2												14		
*09:10-09:11, 09:14	1	2				6										
*09:12, 09:19	1	2						8								
*09:16		2		4												
*09:21	1	2								10						
*09:22	1	2							9							
*09:23	1										11					
*09:24, 09:32	1															
*09:25	1	2										12				
*09:27	1	2											13			
*09:31	1	2													15	
*01:65:01					5	6										
*04:90, 04:155, 04:167, 07:50, 13:220, 16:42, <i>DRB5*01:12</i> , <i>DRB5*01:15</i>											11					
*07:52								8								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control