

INTERPRETATION TABLE			
DRA SSP typing			
Amplification patterns of the DRA alleles			
	well		
	1	2	
Length of spec.	65	100	Length of spec.
PCR product			PCR product
Length of int.	515	430	Length of int.
pos. control			pos. control
5'-primer ²	217	197	5'-primer ²
	5'-gA g ^{3'}	5'-CC C ^{3'}	
3'-primer ³	224	217	3'-primer ³
	5'-g TT ^{3'}	5'-CAA ^{3'}	
Well No.	1	2	Well No.
DRA allele			DRA allele
*0101	1		*0101
*010201-010202		2	*010201-010202
DRA allele			DRA allele
Well No.	1	2	Well No.