

INTERPRETATION TABLE			
DRB1*04 Add-on SSP typing			
	Well		
	1	2	3
Length of spec.	170	75	200
PCR product			
Length of int.	515	430	430
pos. control ¹			
5'-primer(s) ²	13(125)	196(674)	13(125)
	5' -ACA 3'	5' -ACA 3'	5' -ACA 3'
3'-primer(s) ³	56(256)	206(706)	66(286)
	5' -ATC 3'	5' -CAT 3'	5' -gAg 3'
	56(256)		
	5' -ATC 3'		
	58(261)		
	5' -TCA 3'		
Well No.	1	2	3
DRB1 allele			
*04:01:01-04:04:06, 04:04:08, 04:06:01-04:08:02, 04:11:03, 04:13-04:14, 04:16, 04:18-04:23, 04:25, 04:27, 04:31-04:33, 04:35-04:36, 04:38-04:44, 04:46-04:47, 04:49-04:56, 04:59-04:60, 04:63, 04:65, 04:68, 04:70-04:76, 04:78-04:79, 04:85, 04:88, 04:93-04:98:02, 04:100-04:102, 04:105:01-04:105:02, 04:109-04:115, 04:117-04:124, 04:127-04:130, 04:133-04:135, 04:139-04:144, 04:148-04:151, 04:153, 04:155-04:158N	1		
*04:26	W		
*04:37, 04:58, 04:132	1		3
*04:92	1	2	
DRB1 allele			
Well No.	1	2	3