

INTERPRETATION TABLE		
DRB1*12 Add-on SSP typing		
	Well	
	1	2
Length of spec.	215	80
PCR product		
Length of int.	515	430
pos. control ¹		
5'-primer(s) ²	26(165)	-16(40)
	5' -TTA 3'	5' -CAA 3'
3'-primer(s) ³	85(341)	-3(79)
	5' -CAg 3'	5' -AgC 3'
Well No.	1	2
DRB1 allele		
*12:01:01-12:02:05, 12:04-12:07, 12:09, 12:11-12:12, 12:14-12:15, 12:17-12:18, 12:20-12:21, 12:24N-12:26, 12:28-12:35	1	
*12:13	w	
*12:10	1	2
DRB1 allele		
Well No.	1	2