

INTERPRETATION TABLE								
DRB4 SSP subtyping								
Amplification patterns of the DRB4 alleles								
	Well							
	1	2	3	4	5	6	7 ⁸	8
Length of spec.	185	140	130	245	155	190	155	290
PCR product								
Length of int.	515	430	430	515	430	430	430	515
pos. control ¹								
5'-primer ²	28(170)	42(213)	105(401) ⁴	1 st I ⁵	1 st I ⁷	28(170)	96(375)	2 nd I ⁹
	5'-gAT 3'	5'-AgT 3'	5'-AAA 3'	5'-ggg 3'	5'-CAA 3'	5'-gAT 3'	5'-CAA 3'	5'-TgA 3'
3'-primer ³	76(314)	76(314)	135(490) ⁴	5 ⁶	42(213)	77(317)	135(490) ⁴	2 nd I ⁹
	5'-TgT 3'	5'-TgC 3'	5'-gCT 3'	5'-TgC 3'	5'-TCA 3'	5'-AgT 3'	5'-gCC 3'	5'-TTC 3'
Well No.	1	2	3	4	5	6	7 ⁸	8
DRB4 allele								
*01010101	1		3	4				8
*0102		2		4			7	
*01030101, 010302-010304	1			4			7	
*01030102N	1				5		7	
*0104			?	4		6	?	?
*0105	1		?	?			?	?
*0106	1		3	4				?
*0107	1		?	4			?	?
*0201N			3	4				?
*0301N			3					8
DRB4 allele								
Well No.	1	2	3	4	5	6	7 ⁸	8

INTERPRETATION TABLE				
DRB4 SSP subtyping				
Amplification patterns of the DRB4 alleles				
Well				
9	10	11	12	
155	85	110	210	Length of spec. PCR product
515	515	430	430	Length of int. pos. control ¹
42(213)	28(170)	112(421) ⁴	28(170)	5'-primer ²
5' -AgT 3'	5' -gAT 3'	5' -ACT 3'	5' -gAT 3'	
81(328)	42(213)	135(490) ⁴	84(337)	3'-primer ³
5' -gTg 3'	5' -TCA 3'	5' -gCT 3'	5' -CCg 3'	
9	10	11	12	Well No.
				DRB4 allele
				*01010101
				*0102
				*01030101, 010302-010304
				*01030102N
				*0104
9				*0105
		11		*0106
			12	*0107
	10			*0201N
				*0301N
				DRB4 allele
9	10	11	12	Well No.