

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
DRB4 allele ^{1,2,3}																																		
*01:01:01:01, 01:01:01:03-01:01:01:04, 01:01:03-01:01:04, 01:01:06-01:01:10, 01:75, 01:78-01:79, 01:95, 01:125, 01:127, 01:136, 01:141, 01:144	1		3	4																														
*01:01:02	1		?	4																														
*01:01:05, 01:70			3	4																														
*01:02		2		4			7								?		?																	
*01:03:01:01, 01:03:01:03-01:03:01:12, 01:03:01:14-01:03:14, 01:03:16-01:03:18, 01:03:20-01:03:22, 01:63, 01:72, 01:74, 01:76-01:77, 01:80N-01:81, 01:84N-01:89, 01:94, 01:96-01:97, 01:103-01:105, 01:107, 01:112, 01:117-01:124, 01:126, 01:128N-01:132, 01:134-01:135, 01:137-01:138, 01:142-01:143, 01:145N-01:149N	1			4			7																											
*01:03:01:02N, 01:03:01:13N	1				5		7																											
*01:03:15, 01:90	1			4																														
*01:03:19				4			7																											
*01:04			?	4		6	?																											
*01:05	1		?	?			?		9												?													
*01:06	1		3	4							11																							
*01:07:01	1		?	4			?					12									?													
*01:07:02	1		?	4			?					12																						
*01:08. 01:39 ⁴	1		?	4			?						13																					
*01:09	1			4			7							14																				
*01:10	1			4			7								15																			
*01:11	1		?	4			?																							29				
*01:12	1		?	4		6	?						13																					
*01:13	1		?	4			?									16																		
*01:14N	1				5	6	7					12																						
*01:15	1		?	4			?																			25								
*01:16N	1		?	4			?																						28					
*01:17	1		?	4			?									16	17																	
*01:18	1		?	4			?																											
*01:19	1		?	4			?																											
*01:20. 01:36 ⁴	1		?	4			?														21													
*01:21	1		?	4			?								15		17																	
*01:22	1		?	4			?										17			20														
*01:23	1		?	4			?													20														
*01:24. 01:37 ⁴	1		?	4			?																23											
*01:25	1		?	4			?																											
*01:26	1		?	4			?										17										26		28					
*01:27			?	4			?													19	20													
*01:28	1		?	4			?													19														
*01:29	1		?	4			?												18	19														
*01:30	1		?	4			?																											
*01:31	1		?	4			?												18															
*01:32. 01:43. 01:54:01N-01:54:02N ⁴	1		?	4			?																	24										
*01:33. 01:38N. 01:40 ⁴	1		?	4			?														22													
*01:34	1		?	4			?																					26	27					
*01:35	1		?	4			?																	24	25									

Negative Control

[illegible]