

INTERPRETATION TABLE		
DRB4*01030102N SSP typing		
	Well	
	1	2
Length of spec.	155	245
PCR product		
Length of int.	515	430
pos. control ¹		
5'-primer(s) ²	1 st I ⁴	1 st I ⁵
	5'-CAA ^{3'}	5'-ggg ^{3'}
3'-primer(s) ³	42(213)	5 ⁶
	5'-TCA ^{3'}	5'-TgC ^{3'}
Well No.	1	2
DRB4 allele		
*01030102N	1	
*01010101-01030101, 010302-0104, 0106, 0107, 0201N		2
*0105	?	
DRB4 allele		
Well No.	1	2