

INTERPRETATION TABLE		
DRB4*01:03:01:02N SSP typing		
	Well	
	1	2
Length of spec.	155	245
PCR product		
Length of int.	515	430
pos. control ¹		
5'-primer(s) ²	1st I ⁴	1st I ⁵
	5' -CAA 3'	5' -ggg 3'
3'-primer(s) ³	42(213)	5 ⁶
	5' -TCA 3'	5' -TgC 3'
Well No.	1	2
DRB4 allele		
*01:03:01:02N	1	
*01:01:01:01-01:03:01:01, 01:03:02-01:04, 01:06-01:08, 02:01N		2
*01:05		?
DRB4 allele		
Well No.	1	2