

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB5 allele ^{1,2,3}																
*01:01:01	1	2	3	4	5											
*01:01:02	1	2	3	4												
*01:02	1	2	3			6	7									
*01:03	1	2				6	7	8								
*01:04	1	2	3	4					9							
*01:05	1	2	3		5	6										
*01:06				4				8			11					
*01:07	1	2	3	4	5					10						
*01:08N	1	2	3			6	7						13			
*01:09	1	2	3	4	5									14		
*01:10N	1	2	3			6	7								15	
*01:11	1			4				8			11					
*01:12	1	2	3							10					15	
*01:13	1	2	3		5				9							
*01:14	1	2	3											14		
*01:15	1	2	3	4						10					15	
*02:02								8			11	12				
*02:03	1							8			11					
*02:04		2						8				12				
*02:05			3				7					12				
*02:06				4				8			11	12				
DRB1*09:07	1	2	3												15	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control