

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB5 allele ^{1,2,3}																
*01:01:01	1	2	3	4	5											
*01:01:02	1	2	3	4									?			
*01:01:03	1	2	3	4												
*01:02	1	2	3			6	7									
*01:03	1	2				6	7	8					?			
*01:04	1	2	3	4					9							
*01:05	1	2	3		5	6										
*01:06				4				8			11					
*01:07	1	2	3	4	5					10			?			
*01:08N	1	2	3			6	7						13			
*01:09	1	2	3	4	5			8					?			
*01:10N	1	2	3			6	7								15	
*01:11	1			4				8			11					
*01:12	1	2	3							10					15	
*01:13	1	2	3		5				9							
*01:14	1	2	3								11					
*01:15	1	2	3	4						10					15	
*01:16	1	2	3	4	5				9							
*01:17	1	2	3	4	5		7									
*01:18	1	2	3	4	5								13			
*02:02								8			11	12				
*02:03	1							8			11					
*02:04		2						8				12	?			
*02:05			3				7					12				
*02:06				4				8			11	12				
*02:07								8			11	12		14		
DRB1*09:07	1	2	3												15	
DRB1*15:02:03, DRB1*16:01:02, DRB1*16:02:02, DRB1*16:05:01, DRB4*01:05, DRB4*01:07, DRB7*01:01:02, DRB8*01:01:01:01													?			
DRB1*15:86													13			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control