

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
DRB5 allele^{1,2,3}																			
*01:01:01 , 01:19, 01:22-01:24, 01:26, 01:29-01:30	1	2	3	4	5														
*01:01:02	1	2	3	4									?						
*01:01:03	1	2	3	4															
*01:02 , 01:25, 01:28	1	2	3			6	7												
*01:03	1	2				6	7	8					?						
*01:04	1	2	3	4					9										
*01:05	1	2	3		5	6													
*01:06				4				8			11								
*01:07	1	2	3	4	5				10				?						
*01:08N	1	2	3			6	7						13						
*01:09	1	2	3	4	5			8					?						
*01:10N	1	2	3			6	7								15				
*01:11	1			4				8			11								
*01:12	1	2	3							10					15				
*01:13	1	2	3		5				9										
*01:14	1	2	3								11								
*01:15	1	2	3	4					10						15				
*01:16	1	2	3	4	5				9										
*01:17	1	2	3	4	5		7												
*01:18	1	2	3	4	5								13						
*01:20	1	2	3			6	7									16			
*01:21				4				8			11					16			
*01:27N	1	2	3			6	7										17		
*01:31		2	3	4	5														
*02:02 , 02:09, 02:11, 02:14-02:16								8			11	12							
*02:03	1							8			11								
*02:04		2						8				12	?						
*02:05			3				7					12							
*02:06				4				8			11	12							
*02:07								8			11	12		14					
*02:08		2	3			6	7					12				16			
*02:10								8			11	12				16			
*02:12		2	3			6	7					12					17		
*02:13								8			11							18	
<i>DRB1*03:112, DRB1*09:06</i>																			18
<i>DRB1*08:59, DRB1*11:210</i>																16			
<i>DRB1*09:07</i>	1	2	3												15				
<i>DRB1*15:02:03, DRB1*16:01:02, DRB1*16:02:02, DRB1*16:05:01, DRB4*01:05, DRB4*01:07:01</i>													?						
<i>DRB1*15:86</i>													13						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

Negative Control