

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
DRB5 allele ^{1,2,3}																					
*01:01:01:01:01:01:02, 01:19, 01:22-01:24, 01:26, 01:29-01:30, 01:33-01:34, 01:36-01:38, 01:40, 01:42, 01:45, 01:50	1	2	3	4	5																
*01:01:02	1	2	3	4									?								
*01:01:03	1	2	3	4																	
*01:02, 01:25, 01:28, 01:32, 01:35	1	2	3			6	7														
*01:03	1	2				6	7	8					?								
*01:04	1	2	3	4					9												
*01:05	1	2	3		5	6															
*01:06				4				8			11										
*01:07	1	2	3	4	5					10			?								
*01:08N	1	2	3			6	7						13								
*01:09	1	2	3	4	5			8					?								
*01:10N	1	2	3			6	7								15						
*01:11	1			4				8			11										
*01:12	1	2	3							10					15						
*01:13, 01:41	1	2	3		5				9												
*01:14	1	2	3								11										
*01:15	1	2	3	4						10					15						
*01:16	1	2	3	4	5				9												
*01:17	1	2	3	4	5		7														
*01:18	1	2	3	4	5								13								
*01:20	1	2	3			6	7								16						
*01:21				4				8			11				16						
*01:27N, 01:46	1	2	3			6	7										17				
*01:31, 01:44		2	3	4	5																
*01:39						6	7														
*01:43		2		4				8								16					
*01:47	1	2	3	4	5										16						
*01:48N	1	2	3	4	5					10										20	
*01:49N	1	2	3	4	5														19		
*01:51	1	2	3		5																
*02:02:01:02:02:03, 02:09, 02:11, 02:14-02:16, 02:18, 02:20-02:23								8			11	12									
*02:03	1							8			11										
*02:04		2						8				12	?								
*02:05			3				7					12									
*02:06				4				8			11	12									
*02:07								8			11	12		14							
*02:08		2	3			6	7					12			16						
*02:10								8			11	12			16						
*02:12		2	3			6	7					12					17				
*02:13								8			11							18			
*02:17								8			11										
*02:19N								8			11	12					17	18			
*02:24								8				12									
DRB1*03:112, DRB1*09:06																			18		
DRB1*08:59, DRB1*11:210															16						
DRB1*09:07	1	2	3												15						
DRB1*15:02:03, DRB1*16:01:02, DRB1*16:02:02, DRB1*16:05:01, DRB4*01:05, DRB4*01:07:01													?								
DRB1*15:86													13								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

Negative Control