

INTERPRETATION TABLE		
DRB5*01:08N SSP typing		
	Well	
	1	2
Length of spec.	195	175
PCR product		
Length of int.	515	430
pos. control¹		
5'-primer(s)²	107(409)	13(125)
	5' -AgA 3'	5' -gTA 3'
3'-primer(s)³	159(565)	57(258)
	5' -CAT 3'	5' -gCg 3'
		58(260)
		5' -CCT 3'
Well No.	1	2
DRB5 allele		
*01:08N	1	2
*01:01:01-01:07, 01:09-01:14, 02:02-02:06		2
DRB5 allele		
Well No.	1	2