

INTERPRETATION TABLE				
HLA-A*80 SSP typing				
	Well			
	1	2	3	
Length of spec.	135	155	165	
PCR product				
Length of int.	800	1070	1070	
pos. control ¹				
5'-primer ²	176	176	176	
	5' -gCA 3'	5' -gCA 3'	5' -gCA 3'	
3'-primer ³	270	292	299	
	5' -ACA 3'	5' -gTT 3'	5' -CCA 3'	
Well No.	1	2	3	
HLA-A allele				
*80:01	1	2		
*80:02	1		3	
HLA-A allele				
Well No.	1	2	3	