

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
HLA-B allele ^{1,2,3}																																
*13:01:01-13:01:10, 13:51, 13:60, 13:73, 13:77-13:80	1	2	3																													
*13:02:01-13:02:18, 13:54, 13:65-13:66, 13:68, 13:74-13:75, 13:81-13:82	1	2		4																												
*13:03	1				5																											
*13:04	1					6					11																					
*13:06, 13:53	1						7																									
*13:07N	1	2	3																												30	
*13:08	1	2		4					9																							
*13:09, 13:67		2		4						10																						
*13:10	1	2									11																					
*13:11	1	2										12																				
*13:12:01, 13:52 ⁴	1	2	3										13																			
*13:12:02		2	3										13																			
*13:13:01-13:13:02	1	2	3											14																		
*13:14	1	2		4											15																	
*13:15, 13:33 ⁴	1	2		4												16																
*13:16	1	2		4	5																			24								
*13:17	1	2	3													16																
*13:18	1	2		4				8																								
*13:19	1	2		4									13																			
*13:20	1	2	3		5																			24								
*13:21	1	2	3							10				14																		
*13:22:01	1		3		5									15																		
*13:22:02	1		3		5																											
*13:23	1	2	3																											29		
*13:25	1	2	3																			22										
*13:26:01-13:26:02	1	2	3								11																					
*13:27, 13:49N ⁴	1	2		4													17															
*13:28	1	2	3															18														
*13:29	1	2	3																			21										
*13:30	1	2		4																		21										
*13:31	1	2		4	5			8																24								
*13:32	1	2		4														18							24							
*13:34, 13:38 ⁴	1	2		4															19													
*13:35, 13:59	1	2		4										14																		
*13:36	1	2	3																			20										
*13:37	1	2		4																				23								
*13:39		2	3							10																						
*13:40	1	2		4																		21	22									
*13:41	1	2		4	5			8																								
*13:42	1	2		4					9															23								
*13:43	1	2	3															17														
*13:44	1	2		4														18						23								
*13:45	1			4	5																											
*13:46	1					6																										
*13:47	1	2		4																						26						
*13:48	1				5																			24								
*13:50	1	2	3																							26						
*13:55	1	2		4																										29		
*13:56:01N-13:56:02N	1	2		4		6																										
*13:57	1	2	3										13			16																
*13:58		2		4																												
*13:61	1	2	3																									27				
*13:62	1																							24								
*13:63N	1	2	3																										28			
*13:64	1	2																								25						
*13:69	1	2		4																								27				
*13:70	1	2		4																						25						
*13:71	1	2								10				14																		
*13:72	1	2		4		6					11																					
*13:76N	1	2	3																19													
*13:83	1	2	3												15																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	

Negative Control

[illegible]

[illegible]

[illegible]

[illegible]